Package 'espadon'

September 26, 2024

```
Type Package
```

Title Easy Study of Patient DICOM Data in Oncology

Version 1.8.0

Description Exploitation, processing and 2D-3D visualization of DICOM-

RT files (structures, dosimetry, imagery) for medical physics and clinical research, in a patient-oriented perspective.

License GPL-3

URL https://espadon.cnrs.fr

Encoding UTF-8

RoxygenNote 7.3.2

Imports colorspace, DT, graphics, grDevices, igraph, js, mathjaxr, Matrix, methods, misc3d, openxlsx, progress, qs, Rcpp (>= 1.0.10), Rdpack, rgl (>= 1.2.1), Rvcg (>= 0.22.1), shiny, shinyWidgets, sodium, stats

LinkingTo Rcpp

RdMacros mathjaxr, Rdpack

Depends R (>= 4.3.0)

Suggests knitr, rmarkdown, sf

VignetteBuilder knitr

NeedsCompilation yes

Author Cathy Fontbonne [aut, cre] (LPC-CAEN, France), Jean-Marc Fontbonne [aut] (LPC-CAEN, France),

Nathan Azemar [ctb] (LPC-CAEN, France)

Maintainer Cathy Fontbonne <contact.espadon@lpccaen.in2p3.fr>

Repository CRAN

Date/Publication 2024-09-26 13:00:02 UTC

2 Contents

Contents

add.margin
add.shape
bin.closing
bin.clustering
bin.cuboid
bin.dilation
bin.erosion
bin.from.roi
bin.from.vol
bin.intersection
bin.inversion
bin.opening
bin.subtraction
bin.sum
castlow.str
castup.str
dicom.browser
dicom.parser
dicom.patient.anonymiser
dicom.raw.data.anonymizer
dicom.raw.data.loader
dicom.set.tag.value
dicom.tag.dictionary
dicom.tag.parser
dicom.to.Rdcm.converter
dicom.viewer
display.2D.histo
display.3D.contour
display.3D.mesh
display.3D.sections
display.3D.stack
display.DVH
display.DVH.pc
display.dV_dx
display.histo
display.kplane
display.legend
display.obj.links
display.palette
display.plane
espadon.class
export
fan.beam
fan.planar
fan.sphere
fan to voxel

Contents 3

get.extreme.pt	
get.ijk.from.index	74
get.ijk.from.xyz	75
get.line	
get.obj.connection	77
get.plane	78
get.rigid.M	80
get.roi.connection	81
get.value.from.ijk	
get.value.from.mesh	
get.value.from.xyz	
get.volume.from.bin	
get.volume.from.roi	
get.xyz.from.index	
grid.equal	
histo.2D	
histo.DVH	
histo.from.bin	
histo.from.roi	
histo.vol	
load.obj.data	
load.obj.from.dicom	
load.obj.from.Rdcm	
load.patient.from.dicom	
load.patient.from.Rdcm	
load.Rdcm.raw.data	
load.T.MAT	
mesh.from.bin	
mesh.in.new.ref	
mesh.repair	
1	
mesh.spheric.proj	
nesting.bin	
nesting.cube	
nesting.roi	
obj.create	
orientation.create	
pal.rainbow	
pal.RVV	
plot	
Rdcm.inventory	
Rdcm.upgrade	
ref.add	
ref.cutplane.add	
ref.remove	
ref.srctodest.add	
rt.chi.index	
rt.gamma.index	
rt.indices.from.bin	133

4 add.margin

	rt.indices.from.roi	6
	save.T.MAT	3
	save.to.Rdcm	4
	select.names	5
	set.reference.obj	6
	sp.similarity.from.bin	7
	sp.similarity.from.mesh	8
	struct.clustering	0
	struct.from.bin	1
	struct.from.mesh	3
	struct.in.new.ref	4
	struct.merge	5
	study.deployment	6
	toy.dicom.raw	8
	toy.load.patient	9
	vector.product	0
	vol.copy	0
	vol.create	1
	vol.from.bin	2
	vol.gradient	3
	vol.in.new.ref	4
	vol.median	5
	vol.oversampling	6
	vol.regrid	7
	vol.repair	8
	vol.subsampling	9
	vol.sum	1
	xlsx.from.dcm	2
	xlsx.from.Rdcm	3
Index	17:	5

add.margin

Adding or removing a margin to a volume

Description

The add.margin function adds or subtracts a margin of the rectangular parallelepiped circumscribed by a volume.

```
add.margin(vol, xyz.margin, alias = "", description = NULL)
```

add.shape 5

Arguments

vol	"volume" class object.
xyz.margin	Vector of the 3 positive or negative x, y and z margins in mm, in the frame of reference of volume cutting planes.
alias	Character string, \$alias of the created object
description	Character string, describing the created object. If description = NULL (default value), it will be set to vol\$description

Value

Returns a "volume" class object (see espadon.class for class definitions), in which 3D volume is restricted or increased by the requested margins. If the created volume exceeds the initial volume, new voxels are set to NA.

See Also

nesting.cube, nesting.roi and nesting.bin.

Examples

add.shape

Adding shape to a 3D volume.

Description

The add. shape function adds the shape defined by espadon volume object of the modality "binary" or "weight" to a 3D volume.

```
add.shape(obj, shape.bin, shape.mean, shape.sd)
```

6 bin.closing

Arguments

obj Object of the "volume" class and a modality other than "binary" or "weight".

Shape.bin Object of the "volume" class and "binary" or "weight" modality, sharing the same voxels as obj..

Shape.mean, shape.sd

Positive numbers, representing the mean value and the standard deviation of voxels identified by shape.bin.

Value

Returns the 'volume' class object obj, in which the voxels identified by shape.bin have been replaced by a linear interpolation between the initial value and new values. These new values follow a normal distribution with mean shape.mean and standard deviation shape.sd. The weights of the linear interpolation are defined by the voxels values of shape.bin.

See Also

bin.cuboid, bin.cylinder, bin.ellipsoid, bin.from.roi.

Examples

bin.closing

Binary volume closing

Description

The bin.closing function performs a morphological operation of closing, using a sphere, on a "volume" class object of "binary" modality. Closing is useful for:

- filling holes that are smaller than the radius,
- merging two shapes close to each other.

```
bin.closing(vol, radius = 10, alias = "", description = NULL)
```

bin.clustering 7

Arguments

vol	"volume" class object, of "binary" modality
radius	Positive number, or xyz-vector of 3 positive numbers.By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object.alias, "closing r =", radius).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data has been transformed by the closing operation.

See Also

bin.dilation, bin.erosion, bin.opening, add.margin, nesting.cube.

Examples

bin.clustering

Binary volume clustering

Description

The bin.clustering function groups and labels TRUE voxels that have a 6-connectivity (i.e. sharing a common side).

```
bin.clustering(vol, alias = "", description = NULL)
```

8 bin.cuboid

Arguments

vol	"volume" class object, of "binary" modality
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object.alias,"clustering")

Value

Returns a "volume" class object (see espadon.class for class definitions), of "cluster" modality. This object contains the \$cluster.info field, detailing the label and volumes in cm^3 of the different clusters. Note that the label "0" is used for the background.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = "ct",</pre>
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]
# generation of a binary volume
b \leftarrow bin.from.vol(CT, min = -80, max = 20)
\# Display of the n = 3 largest volumes
n <- 3
cluster.b <- bin.clustering (b)</pre>
col <- c ("#00000000", rainbow (n))</pre>
breaks <- seq (-0.5, n + 0.5, length.out = n+2)
display.plane (CT, top = b, main = "Before clustering",
               view.coord = 50, top.col = col, top.breaks = breaks,
               interpolate = FALSE)
display.plane (CT, top = cluster.b, main = "After clustering",
               view.coord = 50, top.col = col, top.breaks = breaks,
                interpolate = FALSE)
```

bin.cuboid

Create a volume containing predefined shapes

Description

These functions create espadon objects of class "volume", and of modality "binary" or "weight", by selecting the voxels defining a rectangular cuboid, an elliptical cylinder or an ellipsoid.

bin.cuboid 9

Usage

```
bin.cuboid(
 back.vol,
 pt.min,
 pt.max,
 ref.shape = back.vol$ref.pseudo,
 T.MAT = NULL,
 modality = "weight",
 alias = "",
 description = NULL
)
bin.cylinder(
 back.vol,
 center,
 radius,
 height,
 ref.shape = back.vol$ref.pseudo,
 T.MAT = NULL,
 modality = "weight",
 alias = "",
 description = NULL
bin.ellipsoid(
 back.vol,
 center,
 radius,
 ref.shape = back.vol$ref.pseudo,
 T.MAT = NULL
 modality = "weight",
 alias = "",
 description = NULL,
)
```

back.vol	"volume" class object.
pt.min, pt.max	Minimum and maximum x , y , z coordinates of the vertices of the rectangular cuboid.
ref.shape	Character string. Pseudonym of the frame of reference in which the requested shape is defined.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. IfT.MAT = NULL, ref.shape must be equal to back.vol\$ref.pseudo.
modality	modality ("binary" or "weight") of the generated object.
alias	Character string, \$alias of the created object.

10 bin.cuboid

description	Character string, describing the created object.
center	Numeric vector of length 3, representing the xyz-center of the cylinder or the ellipsoid.
radius	Positive number, or xy-vector or xyz-vector of 2 or 3 positive numbers, representing the radius of the cylinder or the ellipsoid.
height	Positive number representing the height of the cylinder.
	Additional arguments, such as error, representing the relative error on the volume of the ellipsoid, if back vol $\$$ orientation = $c(1, 0, 0, 0, 1, 0)$.

Value

Returns a "volume" class object of "binary" or "weight" modality (see espadon.class for class definitions), with the same grid as back.vol.

- In the "binary" modality, voxels with 50 percent of their volume within the requested shape are set to TRUE.
- In the "weight" modality, the value of each voxel is its volume fraction included in the requested shape.

See Also

add.shape

```
# Creation of back.vol
CT <- vol.create (c(80, 80,40), c(1.2, 1.2, 2),
                  pt000 = c(-50.4, -50.4, -75), modality = "ct",
                  default.value = as.integer(-997), value.sd = 1)
# Creation of a cuboid
cuboid <- bin.cuboid(CT, pt.min = c(-45.3, -40.2, -60.4),
                     pt.max = c(-15.6, -20.2, -40.2))
# Definition of the shape frame of reference
angle.yz =30
Myz <- cbind(c(1, 0, 0, 0),
             c(0, cos(pi * angle.yz / 180), sin(pi * angle.yz / 180), 0),
             c(0, -\sin(pi * angle.yz / 180), \cos(pi * angle.yz / 180), 0),
             c(0, 0, 0, 1)
angle.xy <- 59
Mxy <- cbind(c(cos(pi * angle.xy / 180), sin(pi * angle.xy / 180), 0, 0) ,</pre>
             c(-sin(pi * angle.xy / 180), cos(pi * angle.xy / 180), 0, 0),
             c(0,0,1,0), c(0,0,0,1)
t.mat <- ref.srctodest.add("ref1", "shaperef", TM = Myz %*% Mxy)</pre>
# Creation of a cylinder
cylinder \leftarrow bin.cylinder(CT, center = c(10.3, 30.6, -30.7),
                         radius = c(10, 20), height = 50,
```

bin.dilation 11

```
ref.shape ="shaperef", T.MAT = t.mat)
# Creation of an ellipsoid
ellipsoid \leftarrow bin.ellipsoid(CT, center = c(-20.1, 0.1, -15),
                       radius = c(19.3, 20.2, 15.3)
# Display
k.idx <- unique(which(cuboid$vol3D.data>0, arr.ind = TRUE)[,3]) - 1
display.3D.stack(cuboid, k.idx, border = FALSE,
               col = c("#FFFFFF00", "#EBDFDFFF", "#D8BFBFFF", "#C59F9FFF",
                      "#B27F7FFF", "#9F5F5FFF", "#8C3F3FFF", "#791F1FFF",
                      "#660000FF"))
k.idx <- unique(which(cylinder$vol3D.data>0, arr.ind = TRUE)[,3]) - 1
display.3D.stack(cylinder, k.idx, border = FALSE,
               "#006600FF"))
k.idx <- unique(which(ellipsoid$vol3D.data>0, arr.ind = TRUE)[,3]) - 1
display.3D.stack(ellipsoid, k.idx, border = FALSE,
               "#000066FF"))
```

bin.dilation

Binary volume dilation

Description

The bin.dilation function enlarges a "volume" class object, of "binary" modality, by means of convolution with a sphere. Dilation is useful for :

- filling holes that are smaller than the radius,
- enlarging capes,
- filling narrow channels,
- merging two shapes close to each other.

Usage

```
bin.dilation(vol, radius = 10, alias = "", description = NULL)
```

vol	"volume" class object, of "binary" modality
radius	Positive number, or xyz-vector of 3 positive numbers.By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = $NULL$ (default value), it will be set to paste ($vol\$object.alias$, "dilataion $r=$ ", radius).

bin.erosion

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which the selected volume has been enlarged by the radius.

See Also

bin.erosion, bin.opening, bin.closing, add.margin, nesting.cube.

Examples

bin.erosion

Binary volume erosion

Description

The bin.erosion function decreases a "volume" class object, of "binary" modality, by means of convolution with a sphere. Erosion is useful for :

- removing volumes that are smaller than the radius,
- eliminating narrow capes,
- enlarging channels,
- turning peninsulas into islands.

Usage

```
bin.erosion(vol, radius = 10, alias = "", description = NULL)
```

vol	"volume" class object, of "binary" modality
radius	Positive number, or xyz-vector of 3 positive numbers. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object.alias, "erosion r =", radius).

bin.from.roi 13

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which the selected volume has been reduced by the radius.

See Also

bin.dilation, bin.opening, bin.closing, add.margin, nesting.cube.

Examples

bin.from.roi

Creation of a binary volume according to RoI

Description

The bin.from.roi function creates a "volume" class object, of modality "binary" or "weight", by selecting the voxels defined by the RoI.

```
bin.from.roi(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  within = TRUE,
  alias = "",
  description = NULL,
  modality = c("binary", "weight"),
  ...
)
```

14 bin.from.roi

Arguments

vol	"volume" class object.
struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default ${\tt roi.name}$ = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo or set to NULL.
within	Boolean, defaults to TRUE. If within = TRUE, the contours included in a RoI are managed, depending on their \$level field. If within = FALSE, only the \$level = 0 fields of the RoI are used (i.e. only the external outlines).
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to struct\$roi.info\$roi.pseudo[roi.idx].
modality	modality ("binary" or "weight") of the generated object.
	additional argument.

Details

roi.name, roi.sname, and roi.idx must select only one RoI.

Value

Returns a "volume" class object of "binary" or "weight" modality (see espadon.class for class definitions), with the same grid as vol.

- In the "binary" modality, voxels with 50 percent of their volume within the ROI are set to TRUE.
- In the "weight" modality, the value of each voxel is its volume fraction included in the ROI.

See Also

bin.from.vol.

bin.from.vol

```
# "optical nerve" binary, with modality "binary" and "weight"
binl <- bin.from.roi (CT, struct = S, roi.sname = "left optical",</pre>
                      alias = "left_optical_nerve", verbose = FALSE,
                      modality = "binary")
binr <- bin.from.roi (CT, struct = S, roi.sname = "right optical",</pre>
                      alias = "right_optical_nerve", verbose = FALSE,
                      modality = "weight")
view.coord <- S$roi.info[S$roi.info$roi.pseudo == "leftopticalnerve",]$Gz</pre>
palette <- grey.colors(100, start = 0, end = 1,</pre>
                       alpha = c(0, rep(1,99))
plot(S, view.coord = view.coord, main ="optical nerves")
plot(bin1, view.coord = view.coord, col = palette,
    cut.interpolate = FALSE, add = TRUE)
plot(binr, view.coord = view.coord, col =palette,
     cut.interpolate = FALSE, add = TRUE)
plot(S, view.coord = view.coord, lwd = 2, add= TRUE)
## Not run:
# with a smaller step
step <- 1
patient <- toy.load.patient (modality = c("ct", "rtstruct"),</pre>
                              roi.name = c("eye", "optical nerve", "brain"),
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]
S <- patient$rtstruct[[1]]</pre>
#' view.coord <- S$roi.info[S$roi.info$roi.pseudo == "leftopticalnerve",]$Gz</pre>
# "optical nerve" binary without inclusions management
bin <- bin.from.roi (CT, struct = S, roi.sname = "left optical",</pre>
                     alias = "left_optical_nerve", within = FALSE,
                     verbose = FALSE)
display.plane (CT, top = bin, struct = S, view.coord = view.coord,
               legend.shift = -80, interpolate = FALSE, main = "Left nerve")
# "optical nerve" binary with inclusions management
bin <- bin.from.roi (CT, struct = S, roi.sname = "left optical",</pre>
                      alias = "left_optical_nerve", within = TRUE,
                      verbose = FALSE)
display.plane (CT, top = bin, struct = S, view.coord = view.coord,
               legend.shift = -80, interpolate = FALSE, main = "Left nerve")
## End(Not run)
```

16 bin.from.vol

Description

The bin.from.vol function creates a "volume" class object, of "binary" modality, in which the voxels fulfilling a condition on their value are selected.

Usage

```
bin.from.vol(
  vol,
  min = -Inf,
  max = Inf,
  in.selection = TRUE,
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object.
min	Minimum value of the selected voxel. Default to -Inf.
max	Maximum value of the selected voxel. Default to +Inf.
in.selection	Boolean, defaults to TRUE. If in.selection = FALSE, the selected pixels are those whose value is not between min and max.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (min, vol\$object.alias, max, sep = "<=") or if in.selection = FALSE, paste ("!(", description, ")").

Value

Returns a "volume" class object of "binary" modality, with the same grid as vol, in which the selected voxels (i.e. set to TRUE) are those satisfying the following conditions:

- If in.selection = TRUE, then min <= vol\$vol3D.data <= max.
- If in.selection = FALSE, then vol\$vol3D.data < min or max < vol\$vol3D.data

bin.intersection 17

bin.intersection

Intersection of two binaries

Description

The bin.intersection function creates a "volume" class object, of "binary" modality, representing the intersection (logical AND) of two binary objects.

Usage

```
bin.intersection(vol1, vol2, alias = "", description = NULL)
```

Arguments

vol1, vol2 "volume" class objects, of "binary" modality.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1\$object.alias, "&", vol2\$object.alias).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, intersection of vol1 and vol2.

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("mr", "rtstruct"),</pre>
                              roi.name = c("brain", "labyrinth processing unit"),
                              dxyz = rep (step, 3))
MR <- patient$mr[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
z.brain <- S$roi.info$Gz[S$roi.info$roi.pseudo == "brain"]</pre>
# Try to discriminate the processing unit with binary selections
bin.brain <- bin.from.roi (MR, struct = S, roi.name = "brain",
                            alias = "brain", T.MAT = patient$T.MAT,
                            verbose = FALSE)
bin.pu.density <- bin.from.vol (MR, min = 160)</pre>
display.plane (MR, top = bin.pu.density, display.ref = S$ref.pseudo,
               view.coord = z.brain, T.MAT = patient$T.MAT,
               interpolate = FALSE, main = "before brain intersection")
bin.pu <- bin.intersection (vol1 = bin.pu.density, vol2 = bin.brain,</pre>
                             alias = "processing unit")
display.plane (MR, top = bin.pu, display.ref = S$ref.pseudo,
               view.coord = z.brain, T.MAT = patient$T.MAT,
               interpolate = FALSE, main = "after brain intersection")
```

18 bin.inversion

bin.inversion	Inversion of a binary
0111.111761.31011	miversion of a binary

Description

The bin.inversion function creates a "volume" class object, of "binary" modality, representing the inverse (logical NOT) of another binary object.

Usage

```
bin.inversion(vol, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste ("!", vol\$object.alias, sep = "").

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, inverse of vol.

bin.opening 19

ne opening

Description

The bin opening function performs a morphological operation of opening, using a sphere, on a "volume" class object of "binary" modality. Opening is useful for:

- removing volumes that are smaller than the radius,
- · smoothing shapes.

Usage

```
bin.opening(vol, radius = 10, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality.
radius	Positive number, or xyz-vector of 3 positive numbers. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object, alias, "opening r =", radius).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data has been transformed by the opening operation.

Note

Opening can be time consuming, try to reduce the binary volume to the strict minimum, before any operations.

See Also

bin.dilation, bin.erosion, bin.closing, add.margin, nesting.cube.

20 bin.subtraction

bin.subtraction

Subtraction of two binaries

Description

The bin.subtraction function creates a "volume" class object of "binary" modality, representing the subtraction of two binary objects.

Usage

```
bin.subtraction(vol1, vol2, alias = "", description = NULL)
```

Arguments

```
vol1, vol2 "volume" class objects of "binary" modality.

alias Character string, $alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1$object.alias, "-", vol2$object.alias).
```

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, in which vol2 is subtracted from vol1.

bin.sum 21

bin.sum

Sum of two binaries

Description

The bin.sum function creates a "volume" class object of "binary" modality, representing the sum (logical OR) of two binary objects.

Usage

```
bin.sum(vol1, vol2, alias = "", description = NULL)
```

Arguments

```
vol1, vol2 "volume" class objects of "binary" modality.

alias Character string, $alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1$object.alias, "+", vol2$object.alias).
```

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, sum of vol1 and vol2.

22 castup.str

castlow.str

Cast of a character string

Description

The castlow.str function converts a word to lowercase, without accents and spaces.

Usage

```
castlow.str(st)
```

Arguments

st

character string

Value

Returns the ASCII//TRANSLIT transcription of the word st, without accents, spaces and in lower-case letters.

See Also

castup.str.

Examples

```
castlow.str (st = c("Right eye", "Left_Lung", "Right-Lung"))
```

castup.str

Cast of a character string

Description

The castup.str function converts a word to upper case, without accents and spaces.

Usage

```
castup.str(st)
```

Arguments

st

character string

dicom.browser 23

Value

Returns the ASCII//TRANSLIT transcription of the word st, without accents, spaces and in capitals.

See Also

castlow.str.

Examples

```
castup.str (st = c("Right eye", "Left_Lung", "Right-Lung"))
```

dicom.browser

DICOM raw data browser

Description

the dicom. browser function creates a dataframe describing the tags contained in the raw data of a DICOM file, as well as the information to access them.

Usage

```
dicom.browser(
  dicom.raw.data,
  nbTAG = 0,
  stop.tag = "",
  stop.level = 0,
  full.info = FALSE,
  tag.dictionary = dicom.tag.dictionary()
)
```

dicom.raw.data	Raw vector, representing the binary extraction of the DICOM file.
nbTAG	Integer. If $nbTAG = 0$ (default), and stop. tag = "", all the DICOM raw data is browsed. Otherwise, the function only browses the first $nbTAG$ tags.
stop.tag	Character string, representing the tag that stops the browse of the dicom.raw.data.
stop.level	$Positive\ integer, specifying\ the\ encapsulation\ level\ of\ the\ {\tt stop.tag}\ in\ {\tt dicom.raw.data}.$
full.info	Boolean. If TRUE, more information about the DICOM data is returned.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

24 dicom.parser

Value

Returns a dataframe if dicom. raw. data is DICOM raw data, NULL otherwise.

If full.info = FALSE, dataframe columns are

- tag: the tags contained in dicom.raw.data,
- VR : value representation of the content of the tag,
- endian: the endianness of the tag content,
- start: the start address in dicom. raw. data of the tag content.
- stop: the stop address in dicom. raw. data of the tag content.

If full.info = TRUE, the following columns are added:

- encaps.load: If the tag contains nested data, this column gives the number of bytes remaining until the end of the nesting. If there are several levels of nesting, these numbers are collapsed and separated by a space.
- load.start: the start address in dicom.raw.data of the tag load size.
- load.stop: the stop address in dicom.raw.dataof the tag load size.
- tag.start: the start address in dicom.raw.data of the tag.

See Also

dicom.raw.data.loader, dicom.tag.parser

Examples

```
# DICOM information dataframe of the dummy raw data toy.dicom.raw ()
df <- dicom.browser (toy.dicom.raw (), full.info = TRUE)
str (df)</pre>
```

dicom.parser

Conversion of DICOM raw data into a dataframe or a list of DICOM TAG information

Description

The dicom.parser function creates a dataframe or a list from DICOM raw data. The created dataframe or list provides information about the content of the DICOM TAGs included in the raw data.

dicom.parser 25

Usage

```
dicom.parser(
  dcm,
  as.txt = TRUE,
  nested.list = FALSE,
  try.parse = FALSE,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary(),
  ...
)
```

Arguments

dcm	espadon object of class "volume", "rtplan", "struct" provided by DICOM files, or DICOM filename, or Rdcm filename, or raw vector representing the binary extraction of the DICOM file.
as.txt	Boolean. If as.txt = TRUE, the function returns a dataframe, a list otherwise.
nested.list	Boolean. Only used if as.txt = FALSE. If nested.list = FALSE, the returned list consists of nested lists.
try.parse	Boolean. If TRUE, the tag with unknown DICOM VR (value representation) is converted into string if possible.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
	Additional argument dicom. browser when previously calculated by dicom.browser. Argument dicom.raw.data (deprecated) replaced by dcm argument. Argument nb or dicom.nb representing the number of DICOM file, when dcm contains multiple DICOM files.

Value

Returns a list of elements or a dataframe, depending on as.list.

If it returns a dataframe, the columns are names TAG, VR (value representation), VM (value multiplicity), loadsize and Value. The field \$Value is a string representation of the true value.

If it returns a list, each of its elements, named by a TAG, is either a vector or a string, depending of the TAG included in dicom.raw.data.

See Also

dicom.raw.data.loader, dicom.tag.parser, dicom.viewer,xlsx.from.dcm,xlsx.from.Rdcm

Examples

```
# content of the dummy raw data toy.dicom.raw (), as a list.
L <- dicom.parser (toy.dicom.raw (), as.txt = FALSE)
str(L[40:57])

L <- dicom.parser (toy.dicom.raw (), as.txt = FALSE, nested.list = TRUE)
str(L[40:45])

# content of the dummy raw data toy.dicom.raw (), as a dataframe.
L <- dicom.parser (toy.dicom.raw (), as.txt = TRUE)
str (L)</pre>
```

dicom.patient.anonymiser

Anonymisation of a patient's DICOM files

Description

the dicom.patient.anonymiser function anonymises all DICOM files in a patient's directory.

Usage

```
dicom.patient.anonymiser(
  dcm.files,
  pat.dest.dir,
  offset = 0,
  new.PIN = "Anonymous patient",
  reset.private.tag = FALSE,
  new.UID = FALSE,
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

	dcm.files	String vector, representing the list of the full names of the DICOM files of the same patient, or its directories.
	pat.dest.dir	Character string,representing the full name of the patient's directory, which will contain the patient's anonymized files.
	offset	Integer, default to 0. Each date of the DICOM will be shifted by this offset expressed in days.
	new.PIN	Character string, representing the PIN remplacing the old one.
reset.private.tag		
		Boolean, if TRUE, the value of tags that are not in the tag. $\operatorname{dictionary}$ is removed.
	new.UID	Boolean. If TRUE, new UID are generated and replace the old ones.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

verbose Boolean. If TRUE, a progress bar indicates loading progress.

Value

Creation of the pat.dest.dir directory, with anonymous DICOM files

Note

The files are anonymized as follows:

- Each date of the DICOM file will be shifted by offset expressed in days.
- Each patient's name, and patient'ID are remplaced by new.PIN
- · All other patient data are deleted, except age, weight, height, gender and shifted birthday.
- All address, phone, physician, operator, author, reviewer, service.
- If reset.private.tag = TRUE, the values of the tags not contained in the tag.dictionary are deleted.

File names are composed of their modality and a part of the SOP UID.

```
# First, save toy.dicom.raw () raw data to a temporary file pat.dir for testing.
temp <- tempdir()</pre>
pat.dir <- file.path (temp, "toy_dcm")</pre>
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "toyrtplan", tmpdir = pat.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
# Files anonymisation
anonymous.pat.dir <- file.path (temp, "anonymous")</pre>
dicom.patient.anonymiser (dcm.files = pat.dir, pat.dest.dir = anonymous.pat.dir,
                           offset = 0, new.PIN = "Anonymous patient",
                           reset.private.tag = TRUE)
lf <- list.files(anonymous.pat.dir, full.names = TRUE)</pre>
dp <- dicom.parser(lf[1])</pre>
dp[grep("^[(]0008|^[(]0010", dp$TAG),]
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

```
dicom.raw.data.anonymizer
```

DICOM anonymizer

Description

the dicom.raw.data.anonymizer function anonymizes dicom.raw.data.

Usage

```
dicom.raw.data.anonymizer(
   dicom.raw.data,
   offset = 0,
   new.PIN = "Anonymous ",
   reset.private.tag = FALSE,
   new.UID = FALSE,
   tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

offset Integer, default to 0. Each date of the DICOM will be shifted by this offset

expressed in days.

new.PIN Character string, representing the PIN remplacing the old one.

reset.private.tag

Boolean, if TRUE, the value of tags that are not in the tag.dictionary is re-

moved.

new.UID Boolean. If TRUE, new UID are generated and replace the old ones.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must

keep. This dataframe is used to parse DICOM files.

Value

Returns an anonymyzed raw vector. See Note.

Note

The raw data is anonymized as follows:

- Each date of the DICOM file will be shifted by offset expressed in days.
- Each patient's name, and patient'ID are remplaced by new.PIN
- All other patient data are deleted, except age, weight, height, gender and shifted birthday.
- All address, phone, physician, operator, author, reviewer, service.
- If reset.private.tag = TRUE, the values of the tags not contained in the tag.dictionary are deleted.

dicom.raw.data.loader 29

Examples

dicom.raw.data.loader DICOM file loading in raw data

Description

the dicom.raw.data.loader function loads a DICOM file as raw data.

Usage

```
dicom.raw.data.loader(dcm.filename)
```

Arguments

dcm. filename Character string, representing the full name of a DICOM file.

Value

Returns a vector of raw data from dcm. filename.

See Also

dicom.browser, dicom.tag.parser

30 dicom.set.tag.value

Examples

dicom.set.tag.value

Change TAG value in DICOM raw data

Description

The dicom.set.tag.value function changes, in the DICOM raw data, the values of the TAG whose VR is a string of characters.

Usage

```
dicom.set.tag.value(
   dicom.raw.data,
   tag,
   tag.value,
   tag.dictionary = dicom.tag.dictionary(),
   ...
)
```

Arguments

dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

tag String vector, representing the list of tags whose value is to be changed. See

note 1.

tag.value String vector, representing the list of new tag values.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must

keep. This dataframe is used to parse DICOM files.

.. Additional arguments dicom. browser when previously calculated by dicom.browser

with argument full.info = TRUE.

dicom.tag.dictionary 31

Value

Returns a raw vector, with new tag values.

Note

- 1- The list of tags included in the DICOM file are given by the first columns of the dataframe provided by the functions dicom.browser and dicom.parser.
- 2- The dicom. set. tag. value function may take some processing time. To minimize this time, it is recommended to prepare in advance all the tags to be modified, and use the dicom. set. tag. value function only once, as shown in the example.

Examples

dicom.tag.dictionary DICOM TAG dictionary

Description

The dicom.tag.dictionary function gives the dictionary of tags used by default in the **espadon** package.

Usage

```
dicom.tag.dictionary(add.dict = c("raysearch.tag"))
```

Arguments

add.dict

Vector of the list of additional dictionaries. Put to NULL, if no additional dictionary is requested.

32 dicom.tag.parser

Value

Returns a 3-column dataframe, describing the VR (value representation) and the name of each DICOM TAG.

This dataframe is the fusion of the "nema.tag" dictionary, provided by *nema* [1], with the dictionaries defined in the add.dict vector:

• "raysearch.tag" dictionary is provided by RaySearch laboratories [2]

References

[1] DICOM nema (Online; accessed 2022-02-16). "Current Edition." https://www.dicomstandard.org/current.

[2] Raysearch Laboratories (Online; accessed 2022-04-25). "RAYPLAN 11A, DICOM Conformance Statement." https://www.raysearchlabs.com/siteassets/raystation-landing-page/dicom-conformance-statements/raystation-pdfs/rsl-d-rs-11a-dcs-en-1.0-2021-05-07-raystation-11a-dicpdf.

Examples

```
str (dicom.tag.dictionary ())
str (dicom.tag.dictionary (NULL))
```

dicom.tag.parser

DICOM TAG parser

Description

the dicom.tag.parser function decodes the content between two DICOM raw data addresses.

Usage

```
dicom.tag.parser(start, stop, VR, endian, dicom.raw.data, try.parse = FALSE)
```

start	Positive integer. Index of the first raw data to parse in the dicom.raw.data.
stop	Positive integer. Index of the last raw data to parse in the dicom.raw.data.
VR	Character string, representing the value representation of DICOM data. See DICOM standard.
endian	Character string, equal to "little" or "big".
dicom.raw.data	Raw vector, representing the binary extraction of the DICOM file.
try.parse	Boolean. If TRUE, the value, with an undocumented VR, is considered, as far as possible, as a string.

dicom.to.Rdcm.converter 33

Value

Returns the dicom.raw.data content between the addresses start and stop. Depending on the representation of the value (VR), it can be a character string or a numerical vector.

Examples

dicom.to.Rdcm.converter

Conversion of DICOM object into files that can be interpreted by the **espadon** package

Description

The dicom.to.Rdcm.converter function creates, for each DICOM object, a *.Rdcm file usefull for using **espadon** package. Each Rdcm file created is referenced by the date of acquisition of the object (if it is not available, its creation date), the patient's PIN, a reference number, an object number in this reference system, and the object modality (mr, ct, rtstruct...).

Usage

```
dicom.to.Rdcm.converter(
  dcm.files,
  pat.dest.dir,
  update = TRUE,
  ignore.duplicates = FALSE,
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

dcm.files	String vector, representing the list of the full names of the DICOM files of the same patient, or its directory.
pat.dest.dir	Character string representing the full name of the patient's directory, which will contain files converted into Rdcm objects.
update	Boolean. If set to TRUE, and if pat.dest.dir contains previously converted files, these files are updated, even if they are duplicated. They retain the same espadon reference frame assignment.

34 dicom.to.Rdcm.converter

```
ignore.duplicates
```

Boolean. If TRUE, the function ignores duplicated objects.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

verbose Boolean. If TRUE, a progress bar indicates the progress of the conversion.

Value

Returns the list of basenames of the created files.

Returns NULL if there are no DICOM files in dcm. files

Note

For each DICOM object, dicom.to.Rdcm.converter creates a *.Rdcm file whose basename is made up of the date of the acquisition (or creation date if previous not found), the patient's PIN, the pseudonym of the frame of reference ("ref1", "ref2"...), the number of the volume object in the directory in this frame of reference ("do1", "do2"...), and the object modality ("mr", "ct", "rtdose", "rtstruct"...).

For example: BASE = "20160514_a008e9ac_ref2_do1_mr"

```
# First, save toy.dicom.raw () raw data to a temporary file for testing.
pat.src.dir <- file.path (tempdir(), "PM_dcm")</pre>
dir.create (pat.src.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PM_rtplan", tmpdir = pat.src.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
# Create a temporary destination directory where the *.Rdcm file will be saved
pat.dest.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dicom.to.Rdcm.converter (pat.src.dir, pat.dest.dir, update = TRUE)
dicom.to.Rdcm.converter (dcm.filename, pat.dest.dir, update = TRUE)
list.files (pat.dest.dir)
# Cleaning temporary directories
unlink (pat.src.dir, recursive = TRUE)
unlink (pat.dest.dir, recursive = TRUE)
```

dicom.viewer 35

dicom.viewer	DICOM content viewer
--------------	----------------------

Description

the dicom.viewer function displays the data of a DICOM file.

Usage

```
dicom.viewer(
  dcm,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary(),
  height = 600,
  width = 900,
  ...
)
```

Arguments

dcm	espadon object of class "volume", "rtplan", "struct" provided by DICOM files, or DICOM filename, or Rdcm filename, or raw vector representing the binary extraction of the DICOM file.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
height, width	Height and width in pixel of the DICOM table.
•••	Additional argument dicom. browser when previously calculated by dicom.browser. Argument nb or dicom.nb representing the number of DICOM file, when dcm contains multiple DICOM files.

Value

Returns the DICOM file description in a browser window.

See Also

xlsx.from.dcm, xlsx.from.Rdcm, dicom.parser

```
if (interactive ()) dicom.viewer (toy.dicom.raw ())
```

36 display.2D.histo

display.2D.histo

Display of a 2D histogram

Description

The display. 2D. histo function displays the density map of a "histo2D" class object.

Usage

```
display.2D.histo(
  histo.2D,
  add = TRUE,
 main = NULL,
 x.lab = NULL,
 y.lab = NULL,
 x.lim = NULL,
 y.lim = NULL,
 bg = "#000000",
  i.rng = NULL,
  display.mode = c("mono.color", "rainbow.color", "line"),
  col = "#ffffff",
  alpha = 1,
  line.pc.levels = c(1, 100),
 line.lwd = 2,
  line.lty = 1
)
```

histo.2D	"histo2D" class object.
add	Boolean indicating whether to display the background image.
main	Title of the background image. If main = NULL, the title indicates "2D histogram".
x.lab	Label of the x-axis of the background image. If $x.lab = NULL$, this label is $histo.2D$x.file.src$
y.lab	Label of the y-axis of the background image. If y.lab = NULL, this label is histo.2D\$y.file.src.
x.lim	Vector, representing the range of the x-axis.
y.lim	Vector, representing the range of the y-axis.
bg	Background color of the image. By default, this color is black.
i.rng	Vector of 2 elements giving the minimum and the maximum intensity of the image. If i.rng = NULL, then the minimum is 0 and the maximum the maximum density.
display.mode	function display mode. See Details.

display.2D.histo 37

col	Color of the color gradient in display.mode = "mono.color", or of the level lines in display.mode = "line". By default, this color is white.
alpha	A number from 0 to 1, indicating the opacity of the image in "rainbow.color" mode. Default alpha = 1.
line.pc.levels	Vector of level lines in percent of maximum density in display.mode = "line". By default lines 1% and 100% are displayed.
line.lwd	Line thickness of the level lines in display.mode = "line".
line.lty	Type of lines for level lines in display.mode = "line".

Details

The display.mode argument can be set to three values: "mono.color", "rainbow.color", or "line". The 2D histogram graph is displayed by default in "mono.color" mode.

- The "mono.color" mode displays a gradient of the color defined by the col argument, depending on the intensity of \$density.map 2-dimensional array.
- The "rainbow. color" mode makes a display according to the "rainbow" palette, while managing the opacity of the colors.
- The "line" mode draws level lines defined in percent by the line.pc.levels argument.

Value

Returns a display of the density map of histo. 2D. This one must be an object of class "histo2D". See espadon.class for class definitions.

See Also

histo.2D.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct", "mr", "rtstruct"),</pre>
                              roi.name = "brain",
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]</pre>
MR <- patient$mr[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
T.MAT <- patient$T.MAT
# restriction of the volume around the RoI
CT.on.roi <- nesting.roi (CT, S, roi.name = "brain", vol.restrict = TRUE,
                          xyz.margin = c (1, 1, 1), alias = CT$description)
MR.on.CT <- vol.regrid (vol = MR, back.vol = CT.on.roi, interpolate = TRUE,
                        T.MAT = T.MAT, alias = CT$description,
                         description = NULL)
# selection of voxels included in the RoI.
roi.bin <- bin.from.roi (vol = CT.on.roi, struct = S, roi.sname = "brain",
```

38 display.3D.contour

display.3D.contour

Display the 3D contours of the RoI

Description

The display. 3D. contour function performs a 3D display of the selected RoI in the chosen coordinate system.

Usage

```
display.3D.contour(
   struct,
   roi.name = NULL,
   roi.sname = NULL,
   roi.idx = NULL,
   roi.col = NULL,
   roi.print = FALSE,
   roi.lwd = 1,
   roi.cex = 1,
   display.ref = struct$ref.pseudo,
   T.MAT = NULL,
   FoR.axis = FALSE,
   FoR.col = "black"
)
```

Arguments

struct	"struct" class object. See espadon.class for class definitions.
roi.name	Vector of exact names of the RoI in the struct object. By default ${\tt roi.name}$ = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.

display.3D.mesh

roi.col	Color of the RoI. If roi.col = NULL (default), the RoI colors are specified in the struct\$roi.info.
roi.print	Boolean vector indicating whether to display the pseudonym of the RoI.
roi.lwd	Line width of the RoI, by default at 1.
roi.cex	Numeric character expansion factor of RoI name if ${\tt roi.print}$ = TRUE, defaults to 1.
display.ref	Pseudonym of frame of reference of the display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, display.ref must be equal to NULL or to struct\$ref.pseudo.
FoR.axis	Boolean or numeric, by default set to FALSE. If FoR.axis = TRUE, the function displays 200 mm length director vectors of the frame of reference. If FoR.axis is numeric, it represent the length in mm of the director vectors.
FoR.col	Color of the frame of reference.

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all of the RoI are selected.

Value

If the concerned regions of interest (RoI) struct exist, it displays the 3D contours of these RoI in the current **RGL** window if it exists, in a new window otherwise.

Examples

display.3D.mesh

3D display of a mesh

Description

The display. 3D. mesh function performs a 3D display of a mesh.

Usage

```
display.3D.mesh(mesh, display.ref = mesh$ref.pseudo, T.MAT = NULL, ...)
```

40 display.3D.sections

Arguments

mesh	"mesh" class object, created by the mesh.from.bin function. See espadon.class for class definitions.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, mesh must be displayed in display.ref = mesh\$ref.pseudo.
	Additional arguments passed to shade3d as color, specular, alpha

Value

Returns a display of mesh in the current RGL window if it exists, in a new window otherwise.

See Also

mesh.from.bin.

Examples

display.3D.sections Display 3D sections of a patient

Description

The display. 3D. sections function displays transverse, sagittal and frontal views at a point in 3D.

Usage

```
display.3D.sections(
  obj,
  cross.pt = c(0, 0, 0),
  display.ref = obj$ref.pseudo,
```

display.3D.sections 41

```
T.MAT = NULL,
col = grey.colors(10, start = 0, end = 1, alpha = c(rep(0, 1), rep(1, 9))),
breaks = NULL,
trans = TRUE,
sagi = TRUE,
front = TRUE,
border = TRUE,
border.col = "#379DA2",
...
)
```

Arguments

obj	"volume" class object to display. See espadon.class for class definitions.
cross.pt	Vector of x, y, z coordinates, representing the cross point of the 3 planes in display.plane.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, obj must be displayed in display.ref = obj\$ref.pseudo.
col	Vector, representing the color palette of the image. Transparent colors are not represented.
breaks	One of:
	NULL: the minimum and the maximum value of the obj define the range.Vector giving the breakpoints of each color.
trans	Boolean. If TRUE (default), the transverse view is displayed.
sagi	Boolean. If TRUE (default), the sagittal view is displayed.
front	Boolean. If TRUE (default), the frontal view is displayed.
border	Boolean. If TRUE (default), the borders of the planes are displayed
border.col	Color of planes borders
	Argument for deprecated arguments

Value

Returns a display of transverse, sagittal and frontal views of obj at cross.pt in the current **RGL** window if it exists, in a new window otherwise. Palette colors are managed by col and breaks.

Examples

42 display.3D.stack

display.3D.stack

Display in 3D the selected planes of an espadon class volume

Description

The display. 3D. stack function displays in 3D the requested cutting planes of a "volume" class object.

Usage

```
display.3D.stack(
  vol,
  k.idx = unique(vol$k.idx[seq(1, vol$n.ijk[3], length.out = 10)]),
  display.ref = vol$ref.pseudo,
  T.MAT = NULL,
  col = grey.colors(10, start = 0, end = 1, alpha = c(rep(0, 1), rep(1, 9))),
  breaks = NULL,
  cube = TRUE,
  border = TRUE,
  ktext = TRUE,
  line.col = "#379DA2",
  line.lwd = 1,
  cex = 1
)
```

Arguments

vol	"volume" class object to display.
k.idx	vector of cutting plane numbers to be displayed, to be chosen in vol\$k.idx. By default k.idx is a vector of 10 uniformly distributed cutting planes in the volume.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, vol must be displayed in display.ref = vol\$ref.pseudo.
col	Vector, representing the color palette of the image. Transparent colors are not represented.
breaks	One of:
	 NULL: The minimum and the maximum value of the vol define the range. Vector giving the breakpoints of each color.
cube	Boolean. If TRUE (default), the "volume" edges are displayed.
border	Boolean. If TRUE (default), the borders of the planes defined in k.idx are displayed.
ktext	Boolean. If TRUE (default), the selected cutting plane numbers are displayed.
line.col	Color of cube, planes and texts displayed.

display.DVH 43

line.lwd Line width of the border and cube, by default at 1.cex Numeric character expansion factor of displayed plan numbers.

Value

Returns a display of the k.idx cutting planes of vol, in the current **RGL** window if it exists, in a new window otherwise. The colors of the palettes are managed by col and breaks.

Examples

display.DVH

Display of a DVH

Description

The display. DVH function displays the Dose Volume Histogram of a "dvh" class object. Y-units are cm^3 .

Usage

```
display.DVH(
   dvh,
   add = FALSE,
   xgrid = TRUE,
   ygrid = TRUE,
   MC.plot = FALSE,
   MC.col = grey.colors(4, rev = TRUE),
   ...
)
```

Arguments

dvh "dvh" class object.

add Boolean indicating whether to display the background image.

xgrid Boolean indicating the display of the x grid.

ygrid Boolean indicating the display of the y grid.

44 display.DVH

MC.plot	Boolean. If MC.plot = TRUE, then display. DVH displays, if they exist, the quantile zones (Prob = 0 , .025, .25, .5, .75, .975, 1) of MC DVH variations.
MC.col	Character string, a valid palette with 4 colours corresponding to 100% , 95% , 50% and median of MC data.
•••	Additional arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the cumulative histogram included in dvh, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the dvh\$vol variations, if they exist.

See Also

```
display.DVH.pc
```

Examples

```
# DVH without MCMC
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "",</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
H <- histo.from.roi (patient$rtdose[[1]], patient$rtstruct[[1]], roi.name = "ptv",</pre>
                     breaks = seq (0, 60, by = 2))
DVH <- histo.DVH (H)
display.DVH (DVH)
## Not run:
# loading of toy-patient objects
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"),</pre>
                              roi.name = "gizzard",
                              dxyz = c (2, 2, 2), beam.nb = 3)
# Calculation of the histogram
H <- histo.from.roi (patient$rtdose[[1]], patient$rtstruct[[1]],</pre>
                     roi.name = "gizzard",
                     breaks = seq (0, 60, by = 1), MC = 100)
# DVH
DVH <- histo.DVH (H)
display.DVH (DVH, MC.plot = TRUE, ylim = c (0, 40))
## End(Not run)
```

display.DVH.pc 45

display.DVH.pc Display of a cumulative DVH in percent of total volume	.ay.DVH.pc	Display of a cumulative DVH in percent of total volume	
---	------------	--	--

Description

The display.DVH.pc function displays the Dose Volume Histogram of "dvh" class object. Y-units are percents of total volume.

Usage

```
display.DVH.pc(
  dvh,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

Arguments

dvh	"dvh" class object. See espadon.class for class definitions.
add	Boolean indicating whether to display the background image.
xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display. DVH.pc displays, if they exist, the quantile zones (Prob = $0, .025, .25, .5, .75, .975, 1$) of MC DVH variations.
MC.col	Character string, a valid palette with 4 colours corresponding to 100% , 95% , 50% and median of MC data.
	Arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot in percent of total volume of the cumulative histogram included in dvh, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the dvh\$pcv variations, if they exist.

See Also

```
display.DVH
```

 $display.dV_dx$

Examples

display.dV_dx

Display of the volume density of a histogram

Description

The display.dV_dx function displays the volume density of a "histo" class object. Y-units are $cm^3.Gy^{-1}$.

Usage

```
display.dV_dx(
  histo,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

Arguments

histo "histo" class object. See espadon.class for class definitions.

Boolean indicating whether to display the background image.

xgrid Boolean indicating the display of the x grid.

Boolean indicating the display of the y grid.

MC.plot Boolean. If MC.plot = TRUE, then display.dV_dx displays, if they exist, the quantile zones (Prob = 0, .025, .25, .5, .75, .975, 1) of variations in volume density.

display.histo 47

MC.col	Character string, a valid palette with 4 colours corresponding to 100%, 95%, 50% and median of MC data.
	Additional arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the differential histogram included in histo, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the histo\$dv_dx variations, if they exist.

See Also

display.histo.

Examples

display.histo

Display of the counts of a histogram

Description

The display.histo function displays the counts of "histo" class object.

Usage

```
display.histo(
  histo,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

48 display.kplane

Arguments

histo	"histo" class object.
add	Boolean indicating whether to display the background image.
xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display.histo displays, if they exist, the quantile zones (Prob = $0, .05, .25, .5, .75, .95, 1$) of variations in counts.
MC.col	Character string, a valid palette with 4 colours corresponding to 100% , 95% , 50% and median of MC data.
	Additional arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the counts included in histo, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the histo\$counts variations, if they exist.

See Also

```
display.dV_dx.
```

Examples

display.kplane

Display of a plane of a volume

Description

The display.kplane function displays the requested plane of a "volume" class object. This function is low-level, used for example in the function display.plane with more intuitive arguments.

display.kplane 49

Usage

```
display.kplane(
  vol,
  k = vol$k.idx[ceiling(length(vol$k.idx)/2)],
 pt00 = c(0, 0),
  dxy = c(1, 1),
  col = grey.colors(255, start = 0, end = 1),
  breaks = NULL,
  sat.transp = FALSE,
 add = FALSE,
 main = NULL,
 abs.lab = "i",
 ord.lab = "j",
  flip = FALSE,
  flop = FALSE,
  bg = "#000000",
  abs.rng = NULL,
 ord.rng = NULL,
  interpolate = FALSE,
)
```

Arguments

vol	"volume" class object to display. See espadon.class for class definitions.
k	Number of the plane to display. By default k is located at mid-plane of the volume.
pt00	Origin point of the displayed plane. By default $pt00 = c (0, 0)$, corresponding to the bottom left of the displayed non-flipped image.
dxy	width and height of a pixel in the plane. If $dxy = c$ (1, 1) (default) abcissa and ordinate correspond to pixel number in the plane.
col	Vector, representing the color palette of the image.
breaks	One of:
	• NULL : the minimum and the maximum value of the vol define the range.
	• Vector giving the breakpoints of each color. Outside values are transparent, leaving the background visible, depending on sat.transp.
sat.transp	Boolean. If TRUE, outside values are transparent, else set to breaks limits colors.
add	Boolean indicating whether to display the background image.
main	Title of the background image. If main = NULL, the title just indicates the value of k.
abs.lab	Label of the image abcissa. By default abs.lab = 'i'.
ord.lab	Label of the image ordinate. By default ord.lab = 'j'.
flip	Boolean defaults to FALSE flipping the horizontal axis of the background image.
flop	Boolean defaults to FALSE flipping the vertical axis of the background image.

50 display.kplane

bg	Background color of the image. By default, this color is black.
abs.rng	Vector of 2 elements indicating the minimum and maximum background image abscissa to display.
ord.rng	Vector of 2 elements indicating the minimum and maximum background image ordinate to display.
interpolate	Boolean, indicating whether to apply linear interpolation to the image.
	others argument of plot function like xaxt, yaxt

Value

Returns a display of the k^{th} image plane of vol.

See Also

display.plane.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("ct","mr", "rtdose"),</pre>
                             dxyz = rep (step, 3), beam.nb = 3)
MR <- patient$mr[[1]]</pre>
CT <- patient$ct[[1]]
D <- patient$rtdose[[1]]</pre>
# display
display.kplane (CT)
display.kplane (MR, k = floor (length(MR$k.idx)*5/8),
                col = grey.colors (256, start = 0, end = 1),
                breaks = seq (0, 500, length.out = 257), bg = "darkblue")
display.kplane (D, k = floor (length(D$k.idx)*3/8),
                col = rainbow (256, s = seq (1, 0, length.out = 256),
                                start = 0, end = 4/6,
                                alpha = seq (0.8, 0, length.out=256),
                                rev = TRUE),
                bg = "darkblue", ord.flip = TRUE, sat.transp = FALSE,
                interpolate = FALSE)
display.kplane (CT, k = floor (length(CT$k.idx)/3), col = pal.RVV (1000),
                bg = "darkblue", ord.flip = TRUE, interpolate = FALSE)
```

display.legend 51

display.legend

Display of the RoI legend

Description

The display.legend function displays in an image the list of requested RoI and their associated color.

Usage

```
display.legend(
   struct = NULL,
   roi.name = NULL,
   roi.sname = NULL,
   roi.idx = NULL,
   lwd = 1,
   cex = 1,
   displayed.roi.name = NULL,
   bg = "black",
   text.col = "white"
)
```

Arguments

struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
lwd	Line thickness, defaults to 1
cex	Font size, default to 1.
displayed.roi.name	
	Vector. If different from NULL, it represents the replacement names of selected RoI if needed.
bg	color of the background.
text.col	color of the legend text.

Details

roi.name, roi.sname, and roi.idx indicates the RoI to display. If all three are set to NULL, all RoI are selected.

52 display.obj.links

Value

Returns display of the RoI names and their associated color in the active graphics window.

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = c("rtstruct"), dxyz = c (5, 5, 5))
S <- patient$rtstruct[[1]]
display.legend (struct = S, roi.idx = 2:10, lwd = 2)</pre>
```

display.obj.links

Display patient objects links

Description

The display.obj.links function displays a graph of connections between objects of a patient. The name of the objects corresponds to their modality (ct, mr, rtdose...) followed by their position in their respective lists in the patient list objects. Connected objects are linked by arrows. Objects sharing the same frame of reference have the same color except for objects with warnings, errors or missing planes which are all in grey. Approved objects are circled in green. By default, objects shapes are circles, except rtdose represented as squares.

Usage

```
display.obj.links(
  pat,
  obj.selected = NULL,
  exclusion = NULL,
  square = "rtdose",
  group.by.connected.FoR = TRUE,
  interactive = FALSE,
  random.seed = 314
)
```

Arguments

pat	"patient" class object, as loaded using load.patient.from.dicom, load.patient.from.Rdcm or toy.load.patient.
obj.selected	Dataframe (default to NULL) containing the objects already selected, created by a previous call of display.obj.links for example.
exclusion	Vector of patient file modalities that should not be displayed, as for instance "mr"
square	Vector of patient file modalities that should be enclosed by a square, as for instance c ("ct", "mr") If NULL no object name is squared.

display.palette 53

group.by.connected.FoR

Boolean. If TRUE (default), all objects sharing the same frame of reference or connected by a registration matrix have the same color. If group.by.connected.FoR

=FALSE, only objects sharing the same FoR have the same color.

interactive Boolean. If interactive = TRUE, buttons are available on the graph to get in-

formation about the objects and select or remove them from the data frame of the selected objects. Then simply click on the name of the object on which to apply the chosen action. If interactive = FALSE no interaction possible with

the plot.

random. seed Positive Integer or NULL. If random. seed = NULL, the objects are laid out ran-

domly. The layout is otherwise fixed.

Value

The function displays all patient objects, linked by an arrow when they are connected or a line when they belongs to the same DICOM object, and with a color and a shape depending on square, group.by.connected.FoR.

When interactive = TRUE, it returns a dataframe of the selected objects, or NULL if no object is selected.

Items are circled in green when the DICOM file has been approved. They are circled in red, when the DICOM series is incomplete (e.g. missing plan).

See Also

get.obj.connection

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (dxyz = c (5, 5, 5), beam.nb = 1)
display.obj.links (patient, group.by.connected.FoR = FALSE)
display.obj.links (patient, group.by.connected.FoR = TRUE)
display.obj.links (patient, group.by.connected.FoR = TRUE, random.seed=NULL)</pre>
```

display.palette

Display of the color scale of a color palette

Description

The Display.palette function displays the color scale as it is used for representations in espadon functions

54 display.palette

Usage

```
display.palette(
  col,
  breaks = NULL,
  factors = NULL,
  override.breaks = FALSE,
  new.window = FALSE,
  ...
)
```

Arguments

col Vector of colors like the ones generated by rainbow, heat.colors, etc.

breaks Vector of breaks for the color palette. It is the usual option for images or dose,

for instance. Its length must be one unit more than col length.

factors Vector containing the labels associated to each col. It should be used for tissue

identification or image segment labelling. Its length must be col length.

override.breaks

Boolean. When FALSE (by default) ordinates are set to breaks. when TRUE colors are uniformely displayed, and associated breaks set to the correct ordinates for

the given colors.

new.window Boolean. If TRUE, it opens a new window for displaying the palette.

... others parameters of plot or axis functions

Value

Returns in a new device (if new.window = TRUE), or in the active graphics window (if new.window = FALSE), the palette color defined by col and breaks in priority, or by col and factors.

Note

the breaks are not necessarily evenly spaced. In this case, the colour palette can be represented as the breaks are defined (default option) or by choosing a constant spacing for each colour and displaying the associated abscissa calculated from the breaks (override.breaks = TRUE).

Examples

display.plane 55

```
# for RVV palette, the function computes breaks between -1000 and 1000
display.palette (pal.RVV (255), new.window = TRUE)
# a palette for dose, for instance
display.palette (rainbow (255, start = 0, end = 4/6, rev = TRUE),
                 breaks = seq (0, 60, length.out = 256), new.window = TRUE)
# black & white palette for CTs or MRs
display.palette (grey.colors (255, start = 0, end = 1),
                 breaks = seq (0, 60, length.out = 256), new.window = TRUE)
# transparency affects colors depending on background (black in first exemple,
# white in the second one)
display.palette (pal.rainbow(255), breaks = seq (0, 60, length.out=256))
display.palette (pal.rainbow(255), breaks = seq (0, 60, length.out=256),
                 bg = "white", new.window = TRUE)
## End(Not run)
# colors contracted range using non uniform breaks in the plot window
display.palette (pal.rainbow(255),
                 breaks = seq (0, 1, length.out = 256)^0.25 * 60, bg="grey",
                 new.window = FALSE)
# the same using breaks override
display.palette (pal.rainbow(255),
                 breaks = seq (0, 1, length.out = 256)^0.25 * 60, bg="grey",
                 override.breaks = TRUE, new.window = FALSE)
```

display.plane

Display the transverse frontal or sagittal view in the patient reference system

Description

The display plane function displays an overlay of images and RoI closed planar contours on a plane defined by the equations x = constant (sagittal view), or y = constant (frontal view) or z = constant (transverse view) in a frame of reference chosen by the user.

Usage

```
display.plane(
  bottom = NULL,
  top = NULL,
  struct = NULL,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  struct.dxyz = c(0.5, 0.5, struct$thickness),
```

56 display.plane

```
display.ref = NULL,
 T.MAT = NULL,
 interpolate = TRUE,
 view.type = c("trans", "front", "sagi"),
 view.coord = 0,
 bg = "#000000",
 abs.rng = NULL,
 ord.rng = NULL,
 bottom.col = grey.colors(255, start = 0, end = 1),
 top.col = pal.rainbow(255),
 bottom.breaks = NULL,
 top.breaks = NULL,
 sat.transp = FALSE,
 struct.lwd = 2,
 main = NULL,
 legend.plot = TRUE,
 legend.shift = 0,
 legend.roi.pseudo = TRUE,
)
```

Arguments

bottom	"volume" class object, displayed using bottom.col palette. If bottom = NULL, no bottom image is displayed.
top	"volume" class object, displayed as an overlay, using top.col palette. If top = NULL, no overlay image is displayed.
struct	"struct" class object. If NULL, no RoI is displayed. Only RoI of closed planar or point type are displayed.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
struct.dxyz	3D vector. Used in case of bottom and top are set to NULL. It represents the voxel size in the display.ref frame of reference, used to calculate contours in frontal or sagittal view.
display.ref	Character string. Pseudonym of the frame of reference used for display. If NULL (default), the bottom image FoR, or top image FoR (when no bottom image), or struct FoR (when no volume displayed).
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, bottom, top and struct must have the same frame of reference.
interpolate	Boolean, indicating whether to apply linear interpolation, when calculating the bottom and top cuts, and then when displaying them. If interpolate = FALSE, the values of the nearest voxels are used. When TRUE (by delfault), trilinear interpolation is used.

display.plane 57

view. type Character string, defining the view to display. It must be set to

- "trans" for a transverse view,
- "front" for a frontal view or.
- "sagi" for a sagittal view.

view. coord Numeric vector of the coordinates along the normal vector of the selected view.

bg Background color of the image. By default, this color is black.

abs.rng Vector of 2 elements indicating the minimum and maximum abscissa to display

on the background image.

ord.rng Vector of 2 elements indicating the minimum and maximum ordinate to display

on the background image.

bottom.col, top.col

Vectors, representing the palette color of bottom and top.

bottom.breaks, top.breaks

One of:

• NULL : the minimum and the maximum value of bottom or top define the range.

• Vector giving the breakpoints of each color. Outside values are transparent, leaving the background visible, depending on sat.transp.

When breaks are specified, the number of breaks must be one unit more then the number of colors.

 ${\tt sat.transp} \qquad \qquad {\tt Boolean.} \ \ {\tt If} \ {\tt TRUE, outside } \ {\tt values} \ {\tt are transparent, else } \ {\tt set to bottom.breaks } \ {\tt or} \ \\$

top.breaks limits.

struct.lwd Line thickness of the RoI contours.

main Character string. When main different from NULL, it replaces the title, and re-

moves the subtitle and the maximum dose indication if top is of modality rtdose.

legend.plot Boolean, that indicates whether the RoI legend should be displayed on the im-

age. It is displayed by default.

legend. shift Numeric. It shifts (in mm) the display of the RoI legend on x-axis.

legend.roi.pseudo

Boolean. If TRUE, the name used for a RoI in the legend comes from the struct\$roi.info\$roi.pseudo

column, otherwise the struct\$roi.info\$name column.

... others parameters of plot function

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all closed planar or point RoI are selected. If a RoI is not present in the requested plane, the RoI legend won't mention it.

Value

Returns a display of the transverse, sagittal or frontal plane. This plane has the coordinate z = view.coord (transverse), y = view.coord (sagittal) pr x = view.coord (frontal). The display is an overlay of:

- a background image of uniform color bg
- the bottom image if it exists
- the top image if it exists
- the contours of the regions of interest if they exist in the plane considered.

Note

- 1- The main title is given by bottom, the subtitle by top.
- 2- When top is in the "rtdose" modality, the maximum dose is written on the image.

See Also

display.kplane, plot.volume, plot.struct, plot.mesh.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 4
patient <- toy.load.patient (modality = c("ct", "mr", "rtstruct", "rtdose"),</pre>
                              roi.name = "",
                              dxyz = rep (step, 3), beam.nb = 3)
CT <- patient$ct[[1]]
MR <- patient$mr[[1]]</pre>
D <- patient$rtdose[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
display.plane (bottom = CT, top = D, struct = S, view.coord = -30,
               interpolate = FALSE, legend.shift = -80)
# Display of CT in reference frame "ref1" and MR in "ref2"
display.plane (bottom = CT, top = MR, interpolate = FALSE)
# Display of CT and MR in reference frame "ref2"
display.plane (bottom = CT, top = MR, interpolate = FALSE, display.ref ="ref2",
               T.MAT = patient$T.MAT)
```

espadon.class

ESPADON class

Description

ESPADON class

Usage

```
espadon.class()
```

Value

Returns a vector of **espadon** class names.

Note

Each object of a class has specific features that are used to display or process that object.

- the "patient" class includes:
 - \$patient: dataframe providing patient's information as PIN, birth date and gender.
 - \$pat.pseudo: patient's pseudonym, initialized to the patient's PIN of \$patient dataframe.
 - \$description: dataframe describing the patient's DICOM objects: their modality (rtstruct, ct, mr, rtplan ...), the base name of the relevant source file in the patient's directory, the pseudonym of their frame of reference (ref1, ref2 ...), their number of sub-objects, their description if any, their numbers of slices/RoI for all sub-objects, their maximum voxels (for volume sub-objects), and finally the aliases of the sub-objects.
 - \$description.by.reg: list of DICOM objects descriptions that are linked by a transfer matrix.
 - \$T.MAT: list of class "t.mat" containing the information of the transfer matrices to move from one frame of reference to another. See load.T.MAT.
 - \$ct: list of CT, if any. They are named by their \$object.alias See load.obj.from.Rdcm.
 - \$mr: list of MRI, if any. They are formatted like the \$ct.
 - \$rtdose: list of dose matrices. They are formatted like the \$ct.
 - \$rtstruct: list of struct objects.
 - ...any DICOM objects other than the reg files, and those previously mentioned, or any modalities created by **espadon**.
 - \$dicom.dvh: if any, list of DVH computed in rt-dose DICOM files.
- the "t.mat" class includes:
 - \$ref.info: dataframe giving the correspondence between the frame of reference (column \$ref) of the DICOM object (TAG (0020,0052)) and its pseudonym (column \$ref_pseudo).
 - \$reg.info: list of dataframes: the first one gives the PID, birthday, and sex of the patient, the second one gives the name of the source file of transfer matrices.
 - \$matrix.description: dataframe giving the transfer matrix names (column \$t), its source frame of reference (column \$src), the destination frame of reference (column \$dest), and its type (\$type). Note that only the RIGID type is supported.
 - \$matrix.list: list of 4X4 transfer matrices. This list contains at least as many Identity matrices as there are ref.pseudo.

A **espadon** object of class "dvh", "fan", "histo", "histo2D", "mesh", "reg", "struct", "t.mat", "undef", "volume" is a list containing at least:

- \$patient: patient's PIN.
- \$patient.name: patient's name.
- \$patient.bd: patient's birthday.

- \$patient.sex: patient's sex
- \$file.basename: vector of .Rdcm or .dcm file basenames of the object, if it exists
- \$file.dirname: directory including the .Rdcm or .dcm file, if it exists
- \$object.name: name of the object.
- \$object.alias: alias of the object.
- \$frame.of.reference: value of TAG (0020,0052).
- \$ref.pseudo: pseudonym of the \$frame.of.reference
- \$modality: modality of the object (e.g. ct, mr, bin, rtplan..)
- \$description: description of the object.
- \$acq.date: date of the content (TAG (0008,0023) for ct and mr and rtimage, TAG (300A,0006) for rtplan, TAG (3006,0008) for rtstruct)
- \$creation.date: creation date of the object.

If the object was generated from a DICOM file, the list also contains:

- \$object.info: Information of the object. It includes:
- the SOP ID (value of TAG (0008,0016)),
- the transfer syntax UID (value of TAG (0002,0010)),
- the SOP implementation ID (value of TAG (0002,0012)),
- the SOP type (value of TAG (0008,0008)),
- the study ID (value of TAG (0020,0010)),
- the study UID (value of TAG (0020,000D)),
- the serie UID (value of TAG (0020,000E)),
- the scanning sequence (value of TAG (0018,0020)),
- the list of SOP labels (values of TAG (0008,0018)),
- the dicom source files,
- the encoding of the content of text tags (values of TAG (0008,0005))
- the description of the study (values of TAG (0008,1030)) in which all letters whose ASCCI code is not between 32 and 123
- the description of the serie (values of TAG (0008,103E)) in which all letters whose ASCCI code is not between 32 and 125
- the number of sub-objects.

if the object is linked to another DICOM object, the list also contains:

- \$ref.object.alias: Alias of the reference object.
- \$ref.object.info: Information of the reference object (not available for mr and ct). It includes:
 - the SOP ID of the reference object,
 - the list of SOP names of the reference object.
- the "dvh" class also includes:
 - \$nb.MC: set to histo\$nb.MC.
 - \$breaks: vector breakpoints.
 - \$mids: vector of cell centers.

• \$mids.unit: Character string, representing the unit of the abcissa of the DVH. For instance, "Gy".

- \$vol: cumulative volume receiving at least the doses defined by \$mids.
- \$pcv: percentage of the total volume receiving at least the doses defined by \$mids.
- if \$nb.MC is different from 0, the arrays MC.vol, MC.pcv and MC.dxyz are added. See histo.DVH.
- the "fan" class also includes:
 - \$origin: the xyz-coordinates of the source point.
 - \$direction: the xyz-coordinates of the main direction of the fan.
 - \$orientation: the xyz-coordinates of the two unit vectors of the plane orthogonal to the \$direction.
 - \$xyz: the xyz-coordinates of the unit vectors of the fan rays
 - \$local: depending on the generation of the fan rays, it can be the spherical coordinates, the deflection angles, the voxel coordinates...
- the "histo" class also includes:
 - \$nb.MC: number of Monte-Carlo simulations
 - \$breaks: vector breakpoints
 - \$mids: vector of cell centers.
 - \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy".
 - counts: count of voxels whose value is included in the limits defined by \$breaks.
 - dV_dx: differential histogram, expressed in cm³ by voxel units, at each \$mids.
 - if \$nb.MC is different from 0, the arrays MC.counts, MC.dV_dx and MC.dxyz are added. See histo.from.roi.
- the "histo2D" class also includes:
 - \$nb.pixels: number of elements in the density.map.
 - \$x.file.src: x label. See histo.2D.
 - \$y.file.src: y label. See histo.2D.
 - x.breaks: vector of x-axis breakpoints.
 - y.breaks: vector of y-axis breakpoints.
 - x.mids: vector of x-axis cell centers.
 - y.mids: vector of y-axis cell centers.
 - density.map: array of densities.
 - total.counts: number of counted voxels.
- the "mesh" class also includes:
 - \$nb. faces: set to the number of faces of the mesh.
 - \$mesh: list of 3 elements defining the mesh: \$vb, \$it and \$normals. See mesh.from.bin.

- the "reg" class also includes:
 - \$nb. of . ref: number of transfer matrices.
 - \$ref.data: list including the lists of information on transfer matrices, namely: the source frame of reference (\$src), the matrix type (\$type, for example 'RIGID') and the transfer matrix (\$matrix).
- the "rtplan" class also includes:
 - \$approval.status: value of TAG (300E,0002).
 - \$number: sub-object number.
 - \$plan.info: dataframe containing, if they exist,
 - \$label the label for the treatment plan,
 - \$plan. name the name for the treatment plan,
 - \$plan. description description of treatment plan,
 - \$tt.protocol the treatment protocol,
 - \$plan. intent the intent of this plan,
 - \$tt.site describing the anatomical treatment site,
 - \$geometry describing whether RT Plan is based on patient or treatment device geometry.
 - \$presc.dose: dataframe containing, if they exist,
 - \$ref.roi.nb value of TAG (3006,0084),
 - \$dose.ref.nb value of TAG (300A,0012),
 - \$dose.ref.id value of TAG (300A,0013),
 - \$struct.type value of TAG (300A,0014),
 - \$description value of TAG (300A,0016),
 - \$pt. coord value of TAG (300A,0018),
 - \$nominal.prior.dose value of TAG (300A,001A),
 - \$dose. type value of TAG (300A,0020),
 - \$constraint.weight value of TAG (300A,0021),
 - \$deliv.warn.dose value of TAG (300A,0022),
 - \$deliv.max.dose value of TAG (300A,0023),
 - \$targ.min.dose value of TAG (300A,0025),

 - \$targ.max.dose value of TAG (300A,0027),
 - \$targ.underdose.vol.frac value of TAG (300A,0028),
 - \$org.risk.full.vol.dose value of TAG (300A,002A),
 - \$org.risk.lim.dose value of TAG (300A,002B),
 - \$org.risk.max.dose value of TAG (300A,002C),
 - \$org.risk.overdose.vol.frac value of TAG (300A,002D)
 - \$fraction.info: dataframe containing, if they exist,
- \$fraction.id the id of the fraction group,
- \$description its description,
- \$planned.frac.nb the total number of treatments (Fractions) prescribed for current fraction group,

- \$frac.pattern.digit.per.day.nb the number of digits in \$frac.pattern used to represent one day,

- \$repeat.frac.cycle.le the number of weeks needed to describe treatment pattern,
- \$frac.pattern the value of TAG (300A,007B) describing treatment pattern every day,
- \$nb. of . beam the number of beams in current fraction group,
- \$beam.dose.meaning the value of TAG (300A,008B) indicating the meaning of Beam Dose,
- \$nb.of.brachy.app the number of brachy application setups in current fraction group.
 - \$fraction.beam (in case of beam treatment): dataframe containing, if they exist,
 - \$fraction.id,
 - \$nb.of.frac.planned,
 - \$beam. dose the value of TAG (00A,0084),
 - \$beam. specif. pt the value of TAG (300A,0082),
 - \$beam.meterset the value of TAG (300A,0086),
 - \$beam. type the value of TAG (300A,0090,
 - \$alt.dose the value of TAG (300A,0091),
 - \$alt. type the value of TAG (300A,0092,
 - \$duration.lim the value of TAG (300A,00C5),
 - \$beam. nb the value of TAG (300C,0006) or (300A,00C0),
 - \$beam. info (in case of beam treatment): dataframe containing, if they exist,
 - \$beam. nb the value of TAG (300C,0006) or (300A,00C0),
 - \$beam. name the value of TAG (300A,00C2),
 - \$beam. description the value of TAG (300A,00C3),
 - \$beam. type the value of TAG (300A,00C4),
 - \$radiation.type the value of TAG (300A,00C6),
 - \$high.dose.technique.type the value of TAG (300A,00C7),
 - \$treatment.machine.name the value of TAG (300A,00B2),
 - \$device.serial.nb the value of TAG (0018,1000),
 - \$primary.dosimeter.unit the value of TAG (300A,00B3),
 - \$referenced. tolerance. table. nb the value of TAG (300C,00A0),
 - \$src.axis.dist the value of TAG (300A,00B4),
 - \$referenced.patient.setup.nb the value of TAG (300C,006A),
 - \$treatment.delivery.type the value of TAG (300A,00CE),
 - \$wedges.nb the value of TAG (300A,00D0),
 - \$compensators.nb the value of TAG (300A,00E0),
 - \$total.compensator.tray.factor the value of TAG (300A,00E2),
 - \$boli.nb the value of TAG (300A,00ED),
 - \$blocks.nb the value of TAG (300A,00F0),
 - \$total.block.tray.factor the value of TAG (300A,00F2),
 - \$final.cumul.meterset.weight the value of TAG (300A,010E),
 - \$ctl.pts.nb the value of TAG (300A,0110),
 - \$radiation.mass.nb the value of TAG (300A,0302),
 - \$radiation.atomic.nb the value of TAG (300A,0304),
 - \$radiation.charge.state the value of TAG (300A,0306),
 - \$scan.mode the value of TAG (300A,0308),

- \$modulated.scan.mode.type the value of TAG (300A,0309),
- \$virtual.src.axis.dist the value of TAG (300A,030A),
- \$total.wedge.tray.water.equ.thickness the value of TAG (300A,00D7),
- \$total.compensator.tray.water.equ.thickness the value of TAG (300A,02E3),
- \$total.block.tray.water.equ.thickness the value of TAG (300A,00F3),
- \$range.shifters.nb the value of TAG (300A,0312),
- \$lateral.spreading.devices.nb the value of TAG (300A,0330),
- \$range.modulators.nb the value of TAG (300A,0340),
- \$fixation.light.azimuthal.angle the value of TAG (300A,0356),
- \$fixation.light.polar.angle the value of TAG (300A,0358).
- \$beam.ctl.pt (in case of beam treatment): list containing, for each beam,
- \$info a data.frame of control points information from DICOM
- \$beam. source the coordinates of the source in the patient frame of reference
- \$beam.direction the coordinates of the beam direction in the patient frame of reference
- \$beam.direction the coordinates of the beam orientation in the patient frame of reference
- \$beam. i socenter the coordinates of the isocenter in the patient frame of reference
- \$spot.map, if they exist, the coordinates of the spots in the patient frame of reference

For the moment, only the rotations of the gantry and the patient support, and the position of the isocenter are taken into account in the calculation of these coordinates.

- \$brachy.info (in case of brachy treatment): dataframe containing, if they exist,
 - \$fraction.id
 - \$nb.of.frac.planned,
 - \$brachy. dose the value of TAG (300A,00A4),
 - \$brachy.nb the value of TAG (300C,000C),
 - \$brachy.specif.pt the value of TAG (300A,00A).
- the "struct" class also includes:
 - \$nb.of.roi: number of regions of interest (RoI).
 - \$thickness: thickness between two consecutive planes of a contour.
 - \$ref.from.contour: reference frame change matrix, from the contour reference frame to the ref.pseudo reference frame
 - \$roi.info: dataframe. Information on RoI contours. It includes the followings columns:
 - -\$number: value of TAG (3006,0084) for the concerned RoI.
 - -\$name: value of TAG (3006,0026) for the concerned RoI.
 - -\$description: value of TAG (3006,0028) for the concerned RoI.
 - -\$generation.algorithm: value of TAG (3006,0036) for the concerned RoI.
 - -\$color: value of TAG (3006,002A) for the concerned RoI.
 - -\$dz: z step between planes for the concerned RoI.
 - -\$roi.pseudo: pseudonym of the RoI \$name. It can be changed by the user.
 - -\$min.x: minimum value x in mm of the RoI. Absent when data = FALSE.

- -\$max.x: maximum value x in mm of the RoI. Absent when data = FALSE.
- -\$min.y: minimum value y in mm of the RoI. Absent when data = FALSE.
- -\$max.y: maximum value y in mm of the RoI. Absent when data = FALSE.
- -\$min.z: minimum value z in mm of the RoI. Absent when data = FALSE.
- -\$max.z: maximum value z in mm of the RoI. Absent when data = FALSE.
- -\$vol: volume in mm³ of the RoI. Absent when data = FALSE.
- -\$Gx: position x in mm of the RoI center of gravity. Absent when data = FALSE.
- -\$Gy: position y in mm of the RoI center of gravity. Absent when data = FALSE.
- -\$Gz: position z in mm of the RoI center of gravity. Absent when data = FALSE.
- -\$continue: boolean, indicating whether the contours are on adjacent planes.
- \$roi.obs: dataframe. RT RoI observations (sequence TAG (3006,0080)). It includes the followings columns:
 - -\$nb: value of TAG (3006,0082) for the concerned RoI.
 - -\$roi.nb: value of TAG (3006,0084) for the concerned RoI.
 - -\$label: value of TAG (3006,0085) for the concerned RoI.
 - -\$code.value: value of TAG (0008,0100) in the Identification code sequence.
 - -\$code.scheme: value of TAG (0008,0102) in the Identification code sequence.
 - -\$code.scheme.v: value of TAG (0008,0103) in the Identification code sequence.
 - -\$code.meaning: value of TAG (0008,0104) in the Identification code sequence.
 - -\$type: value of TAG (3006,00A4) for the concerned RoI.
 - -\$interpreter:value of TAG (3006,00A6) for the concerned RoI.
- \$roi.data: exists only if the data is loaded. Contains the list of contour coordinates. The RoI of list number i is that of line i of roi.info. Each element of the list is a list giving the contour information for each plane, namely:
 - \$type: value of TAG (3006,0042).
 - \$pt: dataframe of the coordinates of the contour points.

If the contour is closed (i.e.\$type = "CLOSED_PLANAR"),

then the first point is repeated at the end.

- \$level: contour inclusion level. If this number is even,
- the inside of the closed contour belongs to the RoI.

Otherwise, if odd, the inside of the closed contour is excluded from the RoI.

- the "undef" class: is used for DICOM objects that will not be processed further by **espadon** functions. It can also include what the user wants.
- the "volume" class also includes:
 - \$number: sub-object number.
 - \$unit: unit of the voxel
 - \$rtdose.info: for rtdose only, dataframe with the type, the dose summation type and referenced beam number.
 - \$n.ijk: vector defining the number of indices i, j, k. The product prod(...\$n.ijk) represents the number of voxels in the 3D volume.

66 export

- \$slice. thickness: thickness in mm of a plane.
- \$min.pixel: minimum value of voxels in the volume.
- \$max pixel: maximum value of voxels in the volume.
- \$dxyz: x, y, z steps in mm.
- \$orientation: value of TAG (0020,0037). Vector, comprising the vectors i and j defining the orientation of the patient with respect to the volume planes.
- \$xyz0: in the patient frame of reference, position of the first voxel of each plane.
- \$xyz.from.ijk: transfer matrix of the voxels i, j, k indices to the position x, y, z in mm in the patient's frame of reference.
- \$k.idx: index of planes in the 3D volume.
- \$missing.k.idx: Boolean, indicating if k is a continuous sequence of integers.
- \$cube.idx: 3D volume vertices indices.
- \$vol3D. data: exists only if the data is loaded. 3D array of the voxel values of the 3D volume.

See Also

toy.load.patient, load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT histo.DVH, histo.vol, histo.from.roi, histo.from.bin, histo.2D, mesh.from.bin, load.obj.from.Rdcm

Examples

```
cat ("espadon class names are:", paste (espadon.class(), collapse = ", "))
```

export

Export espadon objects in dicom format

Description

The export function exports struct class objects and volume class objects with CT or RTDOSE modality in DICOM format.

Usage

```
export(
  obj,
  format = "dcm",
  ref.obj.list = NULL,
  use.original.UIs = FALSE,
  use.original.ref.UIs = TRUE,
  file.prefix = "",
  file.dirname = ".",
  file.name = NULL,
  tag.dictionary = dicom.tag.dictionary(),
  ...
)
```

export 67

Arguments

obj	espadon object of rtstruct, ct or rtdose modality to be exported.
format	Format of the export. For next use.
ref.obj.list	List of espadon objects which are referenced objects of obj. See Details.
use.original.UI	s,use.original.ref.UIs
	Booleans. If TRUE, study instance UID, serie instance UID and image type attribute are those indicated in \$object.info item. Otherwise, They are regenerated. See Details.
file.prefix	String. Prefix added to the generated filename, in case of file.name is NULL.
file.dirname	String. Name of the directory in which files are generated.
file.name	String. Base name of the generated files. in CT modality, a slice number is added as a suffix.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
•••	$Additional\ settings\ such\ as\ NAvalue\ (for\ "volume"\ data),\ '(0020,000D)',\ '(0020,000E)',\ '(0008,0008)'$

Details

The object you want to export may be known in a TPS thanks to these Unique Identifiers (UIs). If you want to create a DICOM object that is different and recognised by your TPS, it is important that the DICOM files you want to create have new UIs: in this case, set the use.original.UIs argument to FALSE. Otherwise, UIs in \$object.info item of your object will be used.

Your object may have been created from another DICOM object (i.e. a reference object). You can, for example, see these links with the function display.obj.links.

If you want to keep this links, you must indicate which objects are references in the ref.obj.list argument, in the form of a list of espadon objects. If these reference objects have their own identifier and you wish to keep them, you must set the use.original.ref.UIs argument to TRUE.

Otherwise (use.original.ref.UIs=FALSE), the UIs of the reference objects will be regenerated. It is therefore important that the reference objects contain all their data such as vol3D.data or roi.data.

It may be useful to impose a study number (tag '(0020,000D)'), serial number (tag '(0020,000E)'), or your Image Type Attribute (tag '(0008,0008)'). In this case, you need to add the arguments '(0020,000D)' = your_study_UID, '(0020,000E)' = your_serial_UID, '(0008,0008)' = your_image_type_attributer.

Value

Returns nothing, but generate DICOM files if conditions are required, and indicates the name or number of files created

Examples

```
# First, save toy patient objects to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "PM_Rdcm")
dir.create (pat.dir, recursive = TRUE)</pre>
```

68 fan.beam

fan.beam

Creation of pyramid fan object with constant angle step.

Description

The fan. beam function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of rectangular pyramid fan. Rays are uniformly distributed by angle.

Usage

```
fan.beam(
  alpha,
  dalpha,
  orientation = c(0, 0, 1, 1, 0, 0),
  origin = c(0, 0, 0),
  ref.pseudo = "ref1",
  frame.of.reference = "",
  alias = "",
  description = "beam fan"
)
```

Arguments

alnha

атриа	restrict number specifying the num tangle of the cometar count.
dalpha	Positive number specifying the step of the angle between the rays of the cone beam.
orientation	Vector orientation of the pyramid base composed by the 2 orthonormal vectors coordinates.

origin Numeric vector, giving the xyz coordinates of the fan origin. By default c (0,

Positive number specifying the half-angle of the conical beam

0,0).

fan.planar 69

```
ref.pseudo Character string, frame of reference pseudonym of the created object.

frame.of.reference Character string, frame of reference of the created object.

alias Character string, $alias of the created object.

description Character string, describing the the created object.
```

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz: a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local: a matrix of 2 columns indicating the deflection angle (in rad) in the main directions defined by orientation.

See Also

fan.planar, fan.sphere, fan.to.voxel.

Examples

```
fan <- fan.beam (alpha = 30, dalpha = 1)
head (fan$xyz)
library (rgl)
open3d ()
points3d (fan$xyz)</pre>
```

fan.planar

Creation of pyramid fan object passing through pixels of a plane.

Description

The fan.planar function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of rectangular pyramid fan. Rays are passing through all pixels of a plane, which represent the pyramid basis. It is for instance useful to compute rt-image.

Usage

```
fan.planar(
  vol,
  k = vol$k.idx[ceiling(length(vol$k.idx)/2)],
  origin = c(0, 0, 0),
  alias = "",
  description = "planar fan"
)
```

70 fan.sphere

Arguments

vol	"volume" class object.
k	Positive number specifying the plane index that the rays of the fan must cross. By default, ${\bf k}$ is the central plane.
origin	Numeric vector, giving the xyz coordinates of the fan origin. By default c (0, 0, 0).
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object.

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz : a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local.coords: a list of the ijkt DICOM coordinates of the crossed plane, and the transfer matrix to xyz.from.ijk to compute xyz coordinates in \$ref.pseudo.

See Also

fan.sphere, fan.beam, fan.to.voxel

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("ct"), dxyz = rep (step, 3))
fan <- fan.planar (patient$ct[[1]], origin = patient$ct[[1]]$xyz0[1,])
head (fan$xyz)
library (rgl)
open3d ()
points3d (fan$xyz)</pre>
```

fan. sphere Creation of spherical fan object.

Description

The fan. sphere function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of a spherical fan.

fan.sphere 71

Usage

```
fan.sphere(
  angle = 1,
  method = c("regular", "random"),
  origin = c(0, 0, 0),
  ref.pseudo = "ref1",
  frame.of.reference = "",
  alias = "",
  description = "spherical fan"
)
```

Arguments

angle	Positive number specifying the angle (or mean angle in case of method = "random") between two nearest vectors.
method	Requested method of fan calculation from among 'regular' and 'random'. By default, method = regular. See details.
origin	Numeric vector, giving the xyz coordinates of the fan origin. By default c (0, 0, 0).
ref.pseudo	Character string, frame of reference pseudonym of the created object.
frame.of.reference	
	Character string, frame of reference of the created object.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object.

Details

The "regular" and "random" method are explained by Deserno [1].

- If method = "regular", the returned vectors composing \$xyz matrix are regularly equidistributed at the specified angle.
- If method = "random", the returned vectors composing \$xyz matrix are randomly equidistributed at the specified angle.

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz: a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local: a matrix of 2 columns indicating the polar angle theta (rad) and the azimuthal angle phi (rad) of each ray are added.

References

[1] Deserno, Markus (Online; accessed 2022-08-24). "How to generate equidistributed points on the surface of a sphere." https://www.cmu.edu/biolphys/deserno/pdf/sphere_equi.pdf.

72 fan.to.voxel

See Also

fan.beam, fan.planar, fan.to.voxel

Examples

```
regular.fan <- fan.sphere (angle = 30)
head (regular.fan$xyz)
random.fan <- fan.sphere (angle = 30, method = "random")
head (random.fan$xyz)
library (rgl)
open3d ()
points3d (regular.fan$xyz)
open3d ()
points3d (random.fan$xyz)</pre>
```

fan.to.voxel

Indices of voxels crossed by a fan

Description

The fan.to.voxel function computes the indices of voxels crossed by a fan. It is useful for retrieving voxel values and voxel indices of a volume (dose or ct) along the fan rays.

Usage

```
fan.to.voxel(vol, fan, restrict = FALSE, vol.value = 1)
```

Arguments

vol "volume" class object.

fan "fan" class object created by fan.sphere for example.

restrict Boolean. If TRUE, only the voxels with a value equal to vol.value are taken

into account.

vol.value Value of the voxels taken into account, in case of restrict = TRUE

Value

Returns a dataframe of 4 columns. Each line gives:

- column "ray.index": the index (i.e. the row number) of the ray concerned in fan\$dxyz,
- column "vol.index": the index of the voxel crossed in vol\$vol.3Ddata,
- column "l.in": the distance between fan source (i.e. fan\$origin) and the first face of the voxel crossed by the ray,
- column "dl": the distance crossed by the ray in the voxel.

If the rays do not cross any voxel, the dataframe has no row.

get.extreme.pt 73

See Also

fan.beam, fan.planar, fan.sphere.

Examples

get.extreme.pt

Coordinates of the extreme points

Description

The get.extreme.pt function returns the x, y, z coordinates of the 2 extreme voxels of the rectangular parallelepiped, containing the objet obj of class volume, struct or mesh. These coordinates are given in the ref.pseudo frame of reference.

Usage

```
get.extreme.pt(obj, ref.pseudo = obj$ref.pseudo, T.MAT = NULL, ...)
```

Arguments

obj	object of class volume or struct or mesh.
ref.pseudo	Pseudonym of the frame of reference in which you want the coordinates.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, ref.pseudo must be equal to obj\$ref.pseudo.
•••	Additional arguments min, max (of voxel) if obj is of class 'volume'. Arguments roi.name, roi.sname, roi.idx if obj is of class 'struct'. Arguments vol (depracated), replaced by obj.

74 get.ijk.from.index

Value

Returns a dataframe of min and max columns, and x, y and z rows.

• If obj is a member of the class volume: the returned dataframe represents the coordinates of the 2 extreme points of the rectangle parallelepiped including all the voxels such as min <= obj\$vol3D.data <= max, if the arguments min or max exist, or including all the voxels otherwise.

- If obj is a member of the class struct: the returned dataframe represents the coordinates of the 2 extreme points of the rectangular parallelepiped including all the selected RoI.
- if obj is a member of the class mesh: the returned dataframe represents the coordinates of the 2 extreme points of the rectangular parallelepiped including all the mesh.

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = "ct", roi.name = "", dxyz = c (5, 5, 5))
CT <- patient$ct[[1]]

# xyz extreme coordinate
get.extreme.pt (CT)
get.extreme.pt (CT, min = 0)</pre>
```

get.ijk.from.index

Conversion of the indices of a point into ijk vector

Description

The get.ijk.from.index function converts the voxel indices of vol\$vol3D.data (for example, obtained with the function which) into a vector or matrix of DICOM indices i, j, k.

Usage

```
get.ijk.from.index(idx, vol)
```

Arguments

idx Index, or matrix of voxel indices of the array vol\$vol3D.data.

vol "volume" class object.

Value

Returns an i, j, k column matrix of the DICOM indices of the points of vol\$vol3D.data.

See Also

```
get.value.from.ijk, display.kplane
```

get.ijk.from.xyz 75

Examples

get.ijk.from.xyz

Indices relating to the coordinates of the points

Description

The get.ijk.from.xyz function calculates the i, j, k DICOM indices of the points given in the patient's reference frame.

Usage

```
get.ijk.from.xyz(xyz = matrix(c(0, 0, 0), ncol = 3), vol, verbose = FALSE)
```

Arguments

xyz	Vector of length 3, corresponding to the x, y, z coordinates (in mm) of a point in the patient's frame of reference, or 3-column matrix of x, y, z coordinates of several points.
vol	"volume" class object.
verbose	Boolean, default to FALSE. If verbose = TRUE, then the xyz coordinates are

Value

printed.

Returns a vector or a matrix of the i, j, k DICOM indices of the x, y, z coordinate points in the patient's frame of reference.

Note

The vector or matrix is made up of real numbers. It is up to the user to make the indices as integer. The indices of the first voxel vol are 0, 0, 0. WARNING: As i,j,k are DICOM indices, they are not directly related to array indices. To get the value of the vol\$vol3D.data, use the function get.value.from.ijk.

76 get.line

Examples

get.line

Image value along an axis

Description

The get.line function calculates the value of the points of a volume vol along an axis in any direction.

Usage

```
get.line(
  vol,
  origin = c(0, 0, 0),
  orientation = c(1, 0, 0),
  grid = seq(-100, 100, 1),
  interpolate = TRUE
)
```

Arguments

vol "volume" class object.

origin Vector of x, y, z coordinates belonging to the line to extract. If interpolate =

FALSE, these coordinates are replaced by the coordinates of the voxel closest to

origin.

orientation Directing vector of the line in the vol frame of reference. This vector is inter-

nally normalized.

grid Vector, representing the curvilinear coordinates on the line to extract.

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the

value of the voxels, relative to the values of adjacent voxels, is performed.

get.obj.connection 77

Value

Returns a dataframe, composed of the columns \$x, \$y, \$z, representing the coordinates of the points where the values are taken in vol volume, the column \$s representing the curvilinear abcissa, and the column \$value representing values along \$s.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 4
patient <- toy.load.patient (modality = "rtdose", roi.name = "",</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
# Dose at maximum dose
origin <- get.xyz.from.index (which.max (D$vol3D.data), D)</pre>
display.plane (bottom = D, view.coord = origin[3],
               bg = "#0000ff")
# Dose profile at x = origin[1] and z = origin[3].
1 <- get.line (vol = D, origin = origin,</pre>
               orientation = c(0, 1, 0), interpolate = FALSE)
plot (1$y, 1$value, type = "1")
grid ()
# Dose profile at y = origin[2] and z = origin[3].
1 <- get.line (D, origin = origin,</pre>
               orientation = c(1, 0, 0), interpolate = FALSE)
plot (1$s, 1$value, type = "1")
grid ()
```

get.obj.connection

List of connections between objects

Description

The get.obj.connection function describes with 4 matrices the different connections between the DICOM objects of the patient.

Usage

```
get.obj.connection(pat)
```

Arguments

pat

"patient" class object, as loaded using load.patient.from.dicom, load.patient.from.Rdcm or toy.load.patient.

78 get.plane

Value

Returns a list of 4 named matrices:

• the \$adjacency matrix matrix specifies the source objects that generated the destination objects: the column names correspond to the destinations, and the row names to the sources.

- the \$same.object matrix specifies the elements belonging to the same DICOM object.
- the \$components matrix specifies the objects belonging to the same study.
- the \$same.ref matrix specifies the objects that share the same frame of reference, or with frames of reference linked in T.MAT (by a DICOM reg file for instance)

See Also

display.obj.links

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (dxyz = c (5, 5, 5), beam.nb = 1)
get.obj.connection (patient)
display.obj.links (patient)</pre>
```

get.plane

Extracting a plane from a volume

Description

The get.plane function extracts a plane from a "volume" class object.

Usage

```
get.plane(
   vol,
   origin = c(0, 0, 0),
   plane.orientation = c(1, 0, 0, 0, 1, 0),
   alias = "plane.n",
   xgrid = NULL,
   ygrid = NULL,
   interpolate = TRUE,
   method = c("NN", "Av")
)
```

get.plane 79

Arguments

vol "volume" class object.

origin Vector of x, y, z coordinates, representing the origin of the plane to extract. If

interpolate = FALSE, these coordinates are replaced by the coordinates of the

voxel closest to origin.

plane.orientation

Vector orientation of the plane in the vol frame of reference, composed by the 2 vectors coordinates of the orthonormal basis of the plane. First vector is x-axis,

and second one is y-axis.

alias \$object.alias of the created object.

xgrid Vector, representing the grid of the plane abscissa. See note.

ygrid Vector, representing the grid of the plane ordinates. See note. If ygrid = NULL,

the ordinate is the line intercepting the volume and the step is set to the projec-

tion of vol\$dxyz onto the ordinate orientation.

interpolate Boolean, default to TRUE.

method method of interpolation, defaulted to 'NN' ie 'Nearest Neighbor'. See Details.

Details

The interpolation method is chosen from:

- 'NN': the value of a voxel is calculated from its nearest adajcent neighbors.
- 'Av': the value of a voxel is the weighted average of the voxels contained in a box, whose sides are automatically calculated from the back.bol\$dxyz steps.

Value

Returns a "volume" class object, containing only a single plane, at k = 0, in the same frame of reference as vol. This returned object has 2 new fields local.xgrid, and local.ygrid, representing the local grids of the abscissa (columns) and ordinate (rows) of the plane.

Returns NULL if plane doesn't exist.

Note

Determination of axes:

- the x-axis has plane.orientation[1:3] as unit vector.
- the y-axis has plane.orientation[4:6] as unit vector.
- If xgrid is not NULL, origin + x.grid * plane.orientation[1:3] are the coordinates of the points on the x axis.
- If ygrid is not NULL, origin + y.grid * plane.orientation[4:6] are the coordinates of the points on the y axis.
- If xgrid or ygrid are NULL, they are determined to represent as closely as possible the initial volume in the required cut.

80 get.rigid.M

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = "mr", dxyz = rep (step, 3))</pre>
MR <- patient$mr[[1]]
# mid-volume point
mid.point <- apply (get.extreme.pt (MR),1,mean)</pre>
plane <- get.plane (MR, origin = mid.point, interpolate = TRUE)</pre>
display.kplane (plane, interpolate = FALSE)
plane <- get.plane (MR, origin = mid.point, xgrid = seq (-50, 50, 1),
                    ygrid = seq (-50, 50, 1), interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
# 3 points on the inclined plane
pts <- t ((MR$xyz.from.ijk %*% MR$cube.idx) [1:3 , c (1, 2, 7)])
orientation <- orientation.create (A = pts[1,], B = pts[2,], C = pts[3,])
origin <- apply (pts, 2, mean)
plane <- get.plane (MR, origin = origin,</pre>
                     plane.orientation = orientation, interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
orientation <- orientation.create (A = c (0, 0, 0) , B = c (1, 1, 0),
                                    C = c (-1, 1, 0)
plane <- get.plane (MR, origin = origin,</pre>
                    plane.orientation = orientation, interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
```

get.rigid.M

Transfer matrix between two frames of reference

Description

The function get.rigid.M provides, from the T.MAT list created by the functions load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT, the 4x4 transfer matrix from the FoR (frame o reference) pseudonym src.ref to the FoR pseudonym dest.ref.

Usage

```
get.rigid.M(T.MAT, src.ref, dest.ref)
```

Arguments

T.MAT	"t.mat" class object, created by the functions load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT
src.ref	Pseudonym of the source frame of reference
dest.ref	Pseudonym of the destination frame of reference

get.roi.connection 81

Value

Returns the 4x4 transfer matrix dest.ref from src.ref.

Examples

get.roi.connection

List of inter-connections between RoI

Description

The get.roi.connection function describes the interconnections between Regions of Interest (RoI), from an imaging volume of "cluster" modality, created by struct.clustering.

Usage

```
get.roi.connection(vol)
```

Arguments

vol

"volume" class object of "cluster" modality, created by struct.clustering

Value

Returns the list of regions of interest (RoI), where each item in the list gives the inter-connections with other RoI.

See Also

struct.clustering

get.value.from.ijk

Description

The get.value.from.ijk function calculates the value of a "volume" class object at DICOM indices i, j, k, whether they are integers or not.

Usage

```
get.value.from.ijk(ijk, vol, interpolate = TRUE, s.ijk = c(1, 1, 1))
```

Arguments

ijk	Vector or 3-column matrix of DICOM indices.
vol	"volume" class object.
interpolate	Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the value of the voxels, relative to the values of adjacent voxels, is performed.
s.ijk	Vector of 3 positive numbers greater than or equal to 1, representing the new voxel ijk-sizes in which averaging is calculated.

Value

Returns a vector of the values of the volume at the requested DICOM indices.

See Also

```
get.ijk.from.index.
```

get.value.from.mesh 83

get.value.from.mesh

Voxel value at a given depth of a mesh

Description

The get.value.from.mesh function is used to retrieve the values of an object of class "volume" at the desired depth of a surface described by a mesh. If the mesh corresponds to the "patient" contour, the zero depth is the skin, positive depths enter the patient and negative depths exit to the outside.

Usage

```
get.value.from.mesh(
  mesh,
  vol,
  method = c("point", "disk", "sphere"),
  depth = 0,
  radius = 5,
  spacing = 1,
  T.MAT = NULL,
  FUN = median,
  ...
)
```

Arguments

mesh	espadon "mesh" class object, or rgl/Rvcg "mesh3d" class object. "mesh3d" class object shall an additional field "ref.pseudo" specifying the mesh frame of reference.
vol	"volume" class object.
method	string specifying the desired method for retrieving measurements in vol. by default "point". Other methods exist, for example "disk" or "sphere. See details.
depth	Numeric value, representing the depth, relative to the surface of the mesh, at which values are retrieved. 0 corresponds to the surface, positive values enter the volume used to define the mesh and negative values leave it.
radius	Positive number, defining the radius of the disk or sphere, depending on the desired method.
spacing	spacing of the measurement points on the disk or sphere.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom, load.T.MAT or ref.add.
FUN	function to be applied to reduce the data ("disk" or "sphere" method) to a single value. Default, median value.
	Additional arguments passed to FUN if needed.

84 get.value.from.mesh

Details

The get.value.from.mesh function works at each vertex of the mesh. It moves along the normal at that point to the desired depth.

- When the method is "point", it simply retrieves the value of the volume vol specified at that point.
- When the method is "disk", the values are retrieved on the disk orthogonal to the normal, with radius radius.
- When the method is "sphere", the values are retrieved inside the sphere of radius radius.

For "disk" or "sphere", the measurement points are spaced by spacing. For radius=5 and spacing=1, "disk" and "sphere" perform 78 and 523 measurements respectively. In both cases, the measured values must be reduced to a single result using the FUN function. By default, espadon uses the median, but it can be provided with more complex functions to filter the data efficiently (see example below).

Value

Returns a vector of values measured at the requested depth, with the desired method, filtered by FUN, at each vertex of the mesh.

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct", "rtstruct"), roi.name = "",</pre>
                               dxyz = rep (step, 3))
CT <- patient$ct[[1]]
S <- patient$rtstruct[[1]]</pre>
# creation of the patient mesh
bin <- bin.from.roi (CT, struct = S, roi.name = "patient",</pre>
                      verbose = FALSE)
mesh.patient <- mesh.from.bin (bin, alias = "patient", verbose = FALSE)</pre>
# density value on the skin contour, extracted from CT
density <- get.value.from.mesh (mesh.patient, CT ,depth = 0)</pre>
# Display of mesh, with RVV pal
density[density < -1000] <- -1000
density[density > 1000] <- 1000</pre>
col <- pal.RVV(255)[cut (density, seq (-1000, 1000, length.out = 256),</pre>
                         include.lowest=TRUE)]
library (rgl)
open3d ()
display.3D.mesh (mesh.patient, col = col)
```

get.value.from.xyz 85

get.value.from.xyz Voxel values

Voxel values on a selection of points

Description

The get.value.from.xyz function calculates the voxel values at the x, y, z coordinate points in the chosen frame of reference.

Usage

```
get.value.from.xyz(
  xyz,
  vol,
  xyz.ref.pseudo = NULL,
  T.MAT = NULL,
  interpolate = TRUE,
  verbose = FALSE
)
```

Arguments

xyz	Vector of length 3, corresponding to the x, y, z coordinates (in mm) of a point in xyz.ref.pseudo frame of reference, or 3-column matrix or dataframe of x, y, z coordinates of several points.
vol	"volume" class object.
xyz.ref.pseudo	ref.pseudo in which the xyz coordinate points are given. This ref.pseudo must exist in the T.MAT list. If ref.pseudo is NULL then the point with coordinates xyz is considered to be in the patient frame of reference vol\$ref.pseudo.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT. If T.MAT = NULL, xyz.ref.pseudo must be equal to vol\$ref.pseudo or NULL.
interpolate	Boolean, default to FALSE. If interpolate = TRUE, a trilinear interpolation of the value of the voxels, relative to the values of adjacent voxels, is performed.
verbose	Boolean, default to FALSE. If verbose = TRUE, then the xyz coordinates are printed.

Value

Returns a vector of the voxel values at the requested coordinates.

See Also

get.xyz.from.index

86 get.volume.from.bin

Examples

get.volume.from.bin

Volume selected by binary volume

Description

The get.volume.from.bin function calculates the volume in cm^3 of the selection specified by a "volume" class object of "binary" modality.

Usage

```
get.volume.from.bin(bin)
```

Arguments

bin

"volume" class object, of "binary" modality.

Value

Returns the volume of the binary selection, in cm^3 .

See Also

get.volume.from.roi

get.volume.from.roi 87

get.volume.from.roi Volume of a region of interest (RoI)

Description

The get.volume.from.roi function extracts the volume in cm^3 of one or more RoI, from the \$roi.info of the "struct" class object.

Usage

```
get.volume.from.roi(struct, roi.name = NULL, roi.sname = NULL, roi.idx = NULL)
```

Arguments

struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all RoI are selected.

Value

Returns a vector of the volumes in cm^3 of the requested RoI.

See Also

get.volume.from.bin, select.names

88 get.xyz.from.index

Examples

get.xyz.from.index

Conversion of the indices of a point, into xyz coordinate vector in the patient's frame of reference

Description

The get.xyz.from.index function converts the indices of a voxel of vol\$vol3D.data (for example, obtained with the function which) into a vector or matrix of x, y, z coordinates in the patient's frame of reference.

Usage

```
get.xyz.from.index(idx, vol)
```

Arguments

idx Index, or matrix of voxel indices in the array vol\$vol3D.data. The first index of the array is 1.

vol "volume" class object.

Value

Returns a column-matrix of coordinates in the patient's reference frame, corresponding to the indices idx.

grid.equal 89

```
get.xyz.from.index (which (D$vol3D.data >= 0.99 * Dmax), D)
# or
get.xyz.from.index (which (D$vol3D.data >= 0.99 * Dmax, arr.ind = TRUE), D)
```

grid.equal

Comparison of the grids of two volume objects

Description

The grid.equal function checks that two volumes share the same grid, i.e. the same frame of reference, the same origin point, and the same dx, dy, dz steps.

Usage

```
grid.equal(vol1, vol2)
```

Arguments

```
vol1, vol2 "volume" class objects
```

Value

Returns TRUE if the 2 volumes share the same grid.

Examples

histo.2D

2D histograms of 2 volumes

Description

The histo.2D function creates a "histo2D" class object, containing the two-dimensional array of histograms of two "volume" class objects that have the same grid.

90 histo.2D

Usage

```
histo.2D(
  x.vol,
  y.vol,
  x.breaks = NULL,
  y.breaks = NULL,
  include.outer = TRUE,
  alias = "",
  description = ""
)
```

Arguments

x.vol, y.vol "volume" class objects. The 2 volumes must have the grid (i.e. share the same voxels location).

x.breaks, y.breaks

Vectors giving the breakpoints of x and y axes. See Details.

include.outer Boolean. If include.outer = TRUE, the values out the x.breaks and y.breaks

of each volume are counted in the first and the last cell of the histograms. They

are not taken into account otherwise.

alias Character string, \$alias of the created object

description Character string, describing the created object.

Details

The arguments x.breaks and y.breaks represent the scales of the x and y axes of 2D-histogram graph. If they are NULL, the histo. 2D function defaults to 256 cells between the values vol\$min.pixel and vol\$max.pixel.

Value

Returns a "histo2D" class object. This is a list including:

- \$patient: set to x.vol\$patient.
- \$patient.name: set to x.vol\$patient.name.
- \$patient.bd: set to x.vol\$patient.bd.
- \$patient.sex: set to x.vol\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the histo2D object.
- \$frame.of.reference: set to x.vol\$frame.of.reference.
- \$ref.pseudo: set to x.vol\$ref.pseudo.
- \$modality: set to "histo2D".

histo.2D 91

- \$description: description of the histo2D object.
- \$creation.date: set to Sys.Date.
- \$nb.pixels: number of elements in the density.map.
- \$x.file.src: set to x.vol\$object.alias
- \$y.file.src: set to y.vol\$object.alias
- x.breaks: vector of x-axis breakpoints.
- y.breaks: vector of y-axis breakpoints.
- x.mids: vector of x-axis cell centers.
- y.mids: vector of y-axis cell centers.
- density.map: array of densities.
- total.counts: number of counted voxels.

See Also

display.2D.histo.

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct", "mr", "rtstruct"),</pre>
                              roi.name = "brain",
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]
MR <- patient$mr[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
T.MAT <- patient$T.MAT
# restriction of the volume around the RoI
CT.on.roi <- nesting.roi (CT, S, roi.name = "brain", vol.restrict = TRUE,
                          xyz.margin = c (1, 1, 1), alias = CT$description)
MR.on.CT <- vol.regrid (vol = MR, back.vol = CT.on.roi, interpolate = TRUE,
                         T.MAT = T.MAT, alias = CT$description,
                         description = NULL)
# selection of voxels included in the RoI.
roi.bin <- bin.from.roi (vol = CT.on.roi, struct = S, roi.sname = "brain",</pre>
                         verbose = FALSE)
MR.select <- vol.from.bin (MR.on.CT, roi.bin, alias = MR$description)</pre>
CT.select <- vol.from.bin (CT.on.roi, roi.bin, alias = CT$description)</pre>
# 2D histogram
H2D <- histo.2D (MR.select, CT.select, x.breaks = seq (50, 400, 10),
  y.breaks = seq (50, 400, 10), alias = "H2D MR CT")
str (H2D)
```

92 histo.DVH

histo.DVH Cumulative Dose Volume Histogram

Description

The histo. DVH function calculates, for each dose, the volume receiving at least this dose.

Usage

```
histo.DVH(histo, alias = "", description = histo$description)
```

Arguments

histo "histo" class object.

alias Character string, \$alias of the created object.

description Character string, describing the the created object. If the description = NULL

(default value), it will be set to histo\$description.

Value

Returns a "dvh" class object. This is a list including:

- \$patient: set to histo\$patient.
- \$patient.name: set to histo\$patient.name.
- \$patient.bd: set to histo\$patient.bd.
- \$patient.sex: set to histo\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the dvh object..
- \$frame.of.reference: set to histo\$frame.of.reference.
- \$ref.pseudo: set to histo\$ref.pseudo.
- \$modality: set to "dvh".
- \$description: description of the dvh object. By default, set to histo\$description.
- \$creation.date: set to Sys.Date.
- \$nb.MC: set to histo\$nb.MC.
- \$breaks: vector breakpoints.
- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the DVH. For instance, "Gy", when vol is a rtdose.
- \$vol: cumulative volume receiving at least the doses defined by \$mids.

histo.from.bin 93

- \$pcv: percentage of the total volume receiving at least the doses defined by \$mids.
- \$MC.vol: cumulative volume associated with histo\$MC.dV_dx, if it exists.
- \$MC.pcv: percentage of the total volume associated with histo\$MC.dV_dx, if it exists.
- \$MC.dxyz: set to histo\$MC.dxyz, if it exists.

See Also

histo.from.roi, histo.from.bin, histo.vol, display.DVH, display.DVH.pc

Examples

histo.from.bin

Histogram according to a binary

Description

The histo.from.bin function computes the voxels histogram of the selection defined by the binary object sel.bin of a "volume" class object.

Usage

```
histo.from.bin(vol, sel.bin, breaks = NULL, alias = "", description = NULL)
```

Arguments

vol	"volume" class object
sel.bin	"volume" class object, of binary modality
breaks	Vector giving the breakpoints between histogram cells. If breaks = NULL, the chosen breakpoints are those used by the hist function by default. If breaks are specified, outside values of vol\$vol3D.data are not taken into account.
alias	Character string, \$alias of the created object
description	Character string, describing the the created object. If the description = NULL (default value), it will be set to vol\$description

94 histo.from.roi

Value

Returns a "histo" class object. See histo.vol.

See Also

histo.from.roi, histo.vol, display.histo, display.dV_dx

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct","rtstruct"), roi.name = "",</pre>
                              dxyz = rep (step, 3))
bin.patient <- bin.from.roi (patient$ct[[1]], struct = patient$rtstruct[[1]],</pre>
                              roi.name = "patient", verbose = FALSE)
# ct histogram in patient volume
H <- histo.from.bin (patient$ct[[1]], sel.bin = bin.patient, breaks = NULL,
                      alias = "patient_hist")
str(H)
## Not run:
# Skin dose histogram
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "",</pre>
                              dxyz = c (2, 2, 2), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
# Creation of the skin contour of 3 mm
bin.patient <- bin.from.roi (D, struct = S, roi.name = "patient",</pre>
                              alias = "patient", verbose = FALSE)
inverse.patient <- bin.inversion (bin.patient, alias = "inv (patient)")</pre>
expansion <- bin.dilation (inverse.patient, radius = 3,
                            alias = "inv (patient) + 3")
contour.3mm <- bin.intersection (bin.patient, expansion,</pre>
                                  alias = "contour 3 mm")
# Dose histogram in this volume
H <- histo.from.bin (D, sel.bin = contour.3mm, breaks = NULL,
                      alias = "Skin dose")
str(H)
## End(Not run)
```

histo.from.roi

Histogram according to a RoI

Description

The histo.from.roi function calculates the histogram of the volume voxels belonging to a RoI.

histo.from.roi 95

Usage

```
histo.from.roi(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  breaks = NULL,
  MC = NULL,
  sd = c(1, 1, 1),
  offset = c(0, 0, 0),
  over.sampling.factor = 1,
  alias = "",
  description = NULL
)
```

"volume" class object

Arguments

vol

	· · · · · · · · · · · · · · · · · · ·
struct	"struct" class object.
roi.name	Exact name of a RoI in struct object. By default roi.name = NULL. See Details.
roi.sname	Name or part of name of a RoI in struct object. By default roi.sname = NULL. See Details.
roi.idx	Value of the index of a RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.
breaks	Vector giving the breakpoints between histogram cells. If breaks = NULL, the chosen breakpoints are those used by the hist function by default. If breaks are specified, outside values of vol\$vol3D.data are not taken into account.
MC	If different from NULL (default value), number of calculations that will be performed, by Monte-Carlo, by randomly moving the chosen RoI over a random distance, generated according to a normal distribution with mean translation defined by offset and standard deviation sd.
sd	Vector representing the standard deviation of distances in the 3 directions x, y and z.
offset	Vector representing the translation of the RoI in the 3 directions x, y and z.
over.sampling.factor	
	Strictly positive integer, or a vector of 3 strictly positive integers, default to 1. Defined to oversample grids of vol. Oversampling can be very time consuming.
alias	Character string, \$alias of the created object
description	Character string, describing the the created object. If the description = NULL (default value), it will be set to struct\$roi.info\$roi.pseudo[roi.idx]

96 histo.from.roi

Details

roi.name, roi.sname, and roi.idx must select only one RoI.

Value

Returns "histo" class object. This is a list including:

- \$alias: alias of the histo object.
- \$description: description of the histo object.
- \$breaks: vector breakpoints
- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy", when vol is a rtdose.
- counts: count of voxels whose value is included in the limits defined by \$breaks.
- dV_dx: differential histogram, expressed in cm3 by voxel units, at each \$mids.
- MC. counts: array of MC rows. Each row i represents the histogram of the voxels contained in the RoI, whose points have been shifted by \$MC.dxyz[i,].
- MC.dV_dx: array of MC rows. Each row i represents the differential histogram of the voxels contained in the RoI, the points of which have been shifted by \$MC.dxyz[i,].
- MC. dxyz: array of MC rows, representing the offset applied to the RoI.

Note

Using Monte-Carlo can be time consuming for large RoI.

If you only want the result just for a translation, use the arguments MC = 1, sd = 0 and offset = desired translation vector.

See Also

histo.vol, histo.from.bin, display.histo, display.dV_dx

histo.vol 97

histo.vol Histogram of a volume

Description

The histo.vol function calculates the voxel values histogram of "volume" class object.

Usage

```
histo.vol(vol, breaks = NULL, alias = "", description = NULL, weight = NULL)
```

Arguments

vol "volume" class object.

breaks Vector giving the breakpoints between histogram cells. If breaks = NULL, the

chosen breakpoints are those used by the hist function by default. If breaks are

specified, outside values of vol\$vol3D. data are not taken into account.

alias Character string, \$alias of the created object.

description Character string, describing the the created object. If the

weight espadon object of the volume class, where weight\$vol3D.data represents the

weight of vol\$vol3D.data. description = NULL (default value), it will be set

to vol\$description.

Value

Returns a "histo" class object. This is a list including:

- \$patient: set to vol\$patient.
- \$patient.name: set to vol\$patient.name.
- \$patient.bd: set to vol\$patient.bd.
- \$patient.sex: set to vol\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the histo object.
- \$frame.of.reference: set to vol\$frame.of.reference.
- \$ref.pseudo: set to vol\$ref.pseudo.
- \$modality: set to "histo".
- \$description: description of the histo object.
- \$creation.date: set to Sys.Date.
- \$nb.MC: set to 0.
- · \$breaks: vector breakpoints

98 load.obj.data

- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy", when vol is a rtdose.
- counts: count of voxels whose value is included in the limits defined by \$breaks.
- dV_dx : differential histogram, expressed in cm^3 by voxel units, at each \$mids.

See Also

histo.from.roi, histo.from.bin, display.histo, display.dV_dx

Examples

load.obj.data

Load data of an espadon class object

Description

The load.obj.data function loads all the data of an **espadon** object of class 'struct' or 'volume'.

Usage

```
load.obj.data(obj, tag.dictionary = dicom.tag.dictionary())
```

Arguments

```
obj struct or "volume" class object
```

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must

keep. This dataframe is used to parse DICOM files in case obj was extracted

from DICOM files.

Value

Returns the **espadon** object with data \$vol3D.data or \$roi.data

See Also

load.obj.from.dicom and load.obj.from.Rdcm

load.obj.from.dicom

Examples

load.obj.from.dicom

Loading an espadon object from DICOM files or folder

Description

Loading an espadon object from DICOM files or folder.

Usage

```
load.obj.from.dicom(
  dcm.files,
  data = TRUE,
  ref.pseudo = "ref1",
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

Arguments

dcm.files	String vector, representing the list of the full names of the DICOM files of the same DICOM object, or its directory.
data	Boolean. Only valid for objects usable by the espadon package, namely ct, mr, rtdose, rtstruct, pt If data = TRUE, either the values of the voxels when modality is (ct, mr, rtdose, pt), or the coordinates of the RoI when modality is rtstruct, are loaded into memory.
ref.pseudo	String, $ref.pseudo$ (i.e. pseudonym of the frame of reference) to assign to the loaded object.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
verbose	Boolean. If TRUE, a progress bar indicates the progress of the conversion.

100 load.obj.from.Rdcm

Value

Returns an **espadon** object of class "dvh", "histo", "histo2D", "mesh", "rtplan", "struct", "undef" or "volume" depending on the object modality. See espadon.class for class definitions.

See Also

load.obj.data and load.obj.from.Rdcm

Examples

load.obj.from.Rdcm

Loading an espadon object from *.Rdcm file

Description

The load.obj.from.Rdcm function loads a DICOM object into memory, creating a list containing the information necessary for its subsequent use with the **espadon** package.

Usage

```
load.obj.from.Rdcm(
  Rdcm.filename,
  data = TRUE,
  nb = NULL,
  upgrade.to.latest.version = FALSE
)
```

Arguments

Rdcm.filename Character string, representing the full name of a *.Rdcm file created by dicom.to.Rdcm.converter. load.obj.from.Rdcm

data

Boolean. Only works for objects usable by the **espadon** package, namely ct, mr, rtdose, rtstruct, pt... If data = TRUE, either the values of the voxels when modality is (ct, mr, rtdose), or the coordinates of the RoI when modality is rtstruct, are loaded into memory.

nh

Vector of integers, active only if data = TRUE, and only operating on rtstruct. If nb = NULL, all the RoI of rtstruct are loaded into memory. Otherwise only data of the RoI indices defined by the vector nb are loaded.

```
upgrade.to.latest.version
```

Boolean. If TRUE, the function attempts to upgrade to the latest version, parsing the DICOM data. It may take longer to load the data. Consider using the Rdcm.upgrade function.

Value

Returns an **espadon** object of class "dvh", "histo", "histo2D", "mesh", "rtplan", "struct", "undef" or "volume" depending on the object modality. See espadon.class for class definitions.

See Also

load.obj.data and load.obj.from.dicom

```
# First, save toy patient objects to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dir.create (pat.dir, recursive = TRUE)
patient <- toy.load.patient (modality = c("ct", "mr"), roi.name = "",</pre>
                              dxyz = c (4, 4, 4))
save.to.Rdcm (patient$ct[[1]], dirname = pat.dir)
save.to.Rdcm (patient$mr[[1]], dirname = pat.dir)
save.T.MAT (patient$T.MAT, dirname = pat.dir)
# Rdcm files in pat.dir
list.files(pat.dir)
CT <- load.obj.from.Rdcm (file.path (pat.dir,</pre>
                                      list.files(pat.dir, pattern="ct1[.]Rdcm")[1]),
                           data=TRUE)
MR <- load.obj.from.Rdcm (file.path (pat.dir,
                                      list.files(pat.dir, pattern="mr1[.]Rdcm")[1]),
                           data=TRUE)
Reg <-load.obj.from.Rdcm (file.path (pat.dir,"ref1_from_ref2.Rdcm"), data=TRUE)</pre>
str(Reg)
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

load.patient.from.dicom

Loading patient data from DICOM files

Description

The load.patient.from.dicom function is used to load or pre-load in memory all patient objects from DICOM files.

Usage

```
load.patient.from.dicom(
  dcm.files,
  data = FALSE,
  dvh = FALSE,
  ignore.duplicates = FALSE,
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

Arguments

dcm.files	String vector, representing the list of the full names of the DICOM files of the same patient, or its directories.
data	Boolean. If data = TRUE, the voxels value of the "volume" class objects, or the coordinates of the RoI (region of interest) of the struct class objects, are loaded into memory.
dvh	Boolean. if dvh = TRUE and if they exist, patient DVH are loaded, for convenience. They are not used as is in espadon package.

ignore.duplicates

Boolean. If TRUE, the function ignores duplicated objects.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must

keep. This dataframe is used to parse DICOM files.

verbose Boolean. If TRUE, a progress bar indicates loading progress.

Value

Returns an **espadon** object of class "patient", describing the information from dcm.files. See espadon.class for a description of the "patient" class.

See Also

dicom.to.Rdcm.converter, load.patient.from.Rdcm, load.obj.data, load.obj.from.dicom, load.obj.from.Rdcm and load.T.MAT.

Examples

```
# First, save toy.dicom.raw () raw data to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "toy_dcm")</pre>
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "toyrtplan", tmpdir = pat.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
# loading patient. Here the toy patient ha only a unique rt-plan object
patient <- load.patient.from.dicom (pat.dir, data = FALSE)</pre>
str (patient, max = 2)
# description of object
patient$description
# transfer matrices :
patient$T.MAT
# rt-plan object
str (patient$rtplan[[1]])
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

load.patient.from.Rdcm

Loading patient data from *.Rdcm files

Description

The load.patient.from.Rdcm function is used to load or pre-load in memory all patient objects converted in *.Rdcm files.

Usage

```
load.patient.from.Rdcm(
  dirname,
  data = FALSE,
  dvh = FALSE,
  upgrade.to.latest.version = FALSE,
  ignore.duplicates = FALSE
)
```

Arguments

dirname

Full paths of the directories of a single patient, or vector of full.path of Rdcm.files.

data

Boolean. If data = TRUE, the voxels value of the "volume" class objects, or the coordinates of the RoI (region of interest) of the struct class objects, are loaded into memory.

104 load.Rdcm.raw.data

dvh Boolean. if dvh = TRUE and if they exist, patient DVH are loaded, for convenience. They are not used as is in **espadon** package.

upgrade.to.latest.version

Boolean. If TRUE, the function attempts to upgrade to the latest version, parsing the DICOM data. It may take longer to load the data. Consider using the Rdcm.upgrade function.

ignore.duplicates

Boolean. If TRUE, the function ignores duplicated objects.

Value

Returns an **espadon** object of class "patient", describing the information contained in dirname. See **espadon.class** for a description of the "patient" class.

See Also

dicom.to.Rdcm.converter, load.patient.from.dicom, load.obj.data, load.obj.from.dicom, load.obj.from.Rdcm and load.T.MAT.

Examples

load.Rdcm.raw.data

Loading a *.Rdcm file

Description

the load.Rdcm.raw.data function loads the content of a *.Rdcm file.

Ioad.Rdcm.raw.data

Usage

```
load.Rdcm.raw.data(
  Rdcm.filename,
  address = TRUE,
  data = TRUE,
  upgrade.to.latest.version = FALSE
)
```

Arguments

Rdcm.filename Character string, representing the full name of a *.Rdcm file created by dicom.to.Rdcm.converter.

address Boolean. If TRUE, a dataframe with the address of the tags in the raw DICOM

data is returned.

data Boolean. If TRUE, the DICOM information are returned as an R list.

upgrade.to.latest.version

Boolean. If TRUE, the function attempts to upgrade to the latest version, parsing the DICOM data. It may take longer to load the data. Consider using the Rdcm.upgrade function.

Value

Returns a list containing the information, converted by espadon, of a DICOM object...

See Also

dicom.to.Rdcm.converter, load.obj.from.Rdcm.

```
# For testing, save first toy.dicom.raw () raw data to a temporary file, and
# convert it in Rdcm fie
pat.src.dir <- file.path (tempdir(), "PM_dcm")</pre>
dir.create (pat.src.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PM_rtplan", tmpdir = pat.src.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
pat.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dicom.to.Rdcm.converter (pat.src.dir, pat.dir, update = TRUE)
lf <- list.files (pat.dir, pattern = "[.]Rdcm$", full.names = TRUE)</pre>
1f
# Inspect Rdcm raw data
L <- load.Rdcm.raw.data (lf[1])
str (L, max.level =3)
# Cleaning temporary directory
unlink (pat.src.dir, recursive = TRUE)
```

106 load.T.MAT

load.T.MAT	Loading of information about transfer matrices between frames of ref-
	erence of patient Rdcm objects.

Description

The load.T.MAT function lists all the frames of reference of the objects included in the patient directory. It concatenates all the information of the reg matrices of a directory, creating, among other things, a list of 4x4 transfer matrices between frames of reference.

Usage

```
load.T.MAT(dirname, upgrade.to.latest.version = FALSE)
```

Arguments

dirname Full paths of the directories of a single patient, or vector of full.path of Rdcm.files. upgrade.to.latest.version

Boolean. If TRUE, the function attempts to upgrade to the latest version, parsing the DICOM data. It may take longer to load the data. Consider using the Rdcm.upgrade function.

Value

Returns a "t.mat" class object. It is a list that includes:

- \$ref.info: dataframe giving the correspondence between the frame of reference (column \$ref) of the DICOM object (TAG (0020,0052)) and its pseudonym (column \$ref_pseudo).
- \$reg.info:list of dataframes: the first one gives the PID, birthday, and sex of the patient, the second one gives the name of the source file of transfer matrices.
- \$matrix.description: dataframe giving the transfer matrix names (column \$t), its source frame of reference (column \$src), the destination frame of reference (column \$dest), and its type (\$type). Note: only the RIGID type is supported.
- \$matrix.list: list of 4X4 transfer matrices. This list contains at least as many Identity matrices as there are ref.pseudo.

mesh.from.bin

```
list.files(pat.dir)
T.MAT <- load.T.MAT (pat.dir)
T.MAT
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)</pre>
```

mesh.from.bin

Creation of a mesh according to a binary volume

Description

The mesh. from. bin function creates a mesh class object from a volume object of "binary" modality.

Usage

Arguments

bin	"volume" class object of "binary" modality.
alias	Character string, \$alias of the mesh defining the \$alias of the created object.
tol	Tolerance in mm, applied for mesh simplification. See vcgClean. The default value, equal to half the smallest voxel edge, limits meshing errors.
smooth.iteration	
	Number of iterations applied in the smoothing algorithm. See vcgSmooth.
smooth.type	character: select smoothing algorithm. Available are "taubin", "laplace", "HClaplace", "fujiLaplace", "angWeight" (and any sensible abbreviations). By default, set to "taubin". See vcgSmooth.
smooth.lambda	numeric: parameter for Taubin smooth. See vcgSmooth.
smooth.mu	numeric: parameter for Taubin smooth. See vcgSmooth.
smooth.delta	numeric: parameter for Scale dependent laplacian smoothing (see reference below).and maximum allowed angle (in radians) for deviation between normals Laplacian (surface preserving). See vcgSmooth.
verbose	Boolean, by default set to FALSE. Allows you to inhibit comments.

108 mesh.from.bin

Value

Returns a "mesh" class object. This is a list including the following 6 elements:

- \$patient: set to bin\$patient.
- \$patient.bd: set to bin\$patient.bd.
- \$patient.name: set to bin\$patient.name.
- \$patient.sex: set to bin\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: set to the alias argument of the function.
- \$frame.of.reference: set to bin\$frame.of.reference.
- \$ref.pseudo: set to bin\$ref.pseudo.
- \$modality: set to "mesh".
- \$description: By default, set to paste (bin\$object.alias, "mesh").
- \$creation.date: set to Sys.Date.
- \$nb. faces: set to the number of faces of the mesh.
- \$mesh: list of 3 elements defining the mesh:
- \$vb: array made up of the generalized coordinates (x, y, z, 1) of the vertices of the triangles.

There are as many columns as there are vertices.

- \$it: array of the 3 indices of the vertices forming a triangle, arranged by column.

There are as many columns as there are triangles in the mesh.

- \$normals: array made up of the generalized coordinates (x, y, z, 1) of the normal vectors of each triangle.

There are as many columns as there are vertices.

Note

To compute the mesh, all NA voxels of the binary volume bin are set to FALSE. If all voxels are equal to FALSE, the function returns the code NULL.

See Also

vcgSmooth

mesh.in.new.ref

```
# creation of the patient mesh
bin <- bin.from.roi (CT, struct = S, roi.name = "patient", verbose = FALSE)
mesh.patient <- mesh.from.bin (bin, alias = "patient", verbose = FALSE)
str (mesh.patient)</pre>
```

mesh.in.new.ref

Change of frame of reference of a mesh

Description

The mesh.in.new.ref function allows you to change the frame of reference of a mesh.

Usage

```
mesh.in.new.ref(
  mesh,
  new.ref.pseudo,
  T.MAT = NULL,
  alias = "",
  description = NULL)
```

Arguments

mesh "volume" class object.

new.ref.pseudo pseudonym of the frame of reference in which the mesh should be located. This new.ref.pseudo must exist in the T.MAT list.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom, load.T.MAT or ref.add.

alias Character string, \$alias of the created object.

Character string, describing the created object. If description = NULL (default value), it will be that of the mesh.

Value

Returns "mesh" class object in the new frame of reference new.ref.pseudo.

110 mesh.repair

mesh.repair

Repair of a mesh

Description

The mesh.repair function repairs holes in a mesh class object.

Usage

```
mesh.repair(mesh, verbose = TRUE)
```

Arguments

mesh "mesh" class object.

verbose Boolean, by default set to FALSE. Allows you to inhibit comments.

Value

Returns a mesh, repaired by removing degenerated triangles and filling holes.

mesh.spheric.proj 111

mesh.spheric.proj

Adding spherical coordinates to a mesh

Description

The mesh.spheric.proj function adds latitude and longitude coordinates to a mesh. These features map the mesh surface to a sphere. Latitude and longitude are computed using the heat diffusion approach explained by *Brechbühler and al* [1].

Usage

```
mesh.spheric.proj(mesh, verbose = TRUE)
```

Arguments

mesh

"mesh" class object.

verbose

Boolean, by default set to FALSE. Allows you to inhibit comments.

Value

returns a "mesh" class object in which mesh contains Lat and 1on evaluated at vertices. The function allows to have a parameterized surface for later computations as curvature or shape index, hence, nor the surface, nor the angles are preserved. In the DICOM frame of reference, latitude goes along Z axis (from feet = -1 to head = +1) and longitude turns counter clockwise (from -1 to +1).

Note

This funtion is time consuming.

References

[1] Brechbuhler C, Gerig G, Kubler O (1995). "Parametrization of Closed Surfaces for 3-D Shape Description." *Computer Vision and Image Understanding*, **61**(2), 154-170. ISSN 1077-3142, doi:10.1006/cviu.1995.1013.

nesting.bin

nesting.bin

Restrict volume to a binary selection

Description

The nesting.bin function restricts a "volume" class object to the rectangular parallelepiped circumscribed to the selected voxels.

Usage

```
nesting.bin(
  obj,
  sel.bin,
  alias = "",
  description = NULL,
  xyz.margin = c(0, 0, 0),
  obj.restrict = FALSE,
  ...
)
```

Arguments

```
obj
                  "volume" class object, containing data to restrict.
                  "volume" class object, of "binary" modality, specifying the selected voxels.
sel.bin
                  Character string, $alias of the created object.
alias
                  Character string, describing the tre created object. If description = NULL, it
description
                  will be paste (obj$description, "restricted to", sel.bin$description).
xyz.margin
                  Numeric vector of length 3, by default set to c(0, 0, 0). See details.
obj.restrict
                  Boolean. Used if obj is of class "volume". If obj.restrict = TRUE, the rectan-
                  gular parallelepiped circumscribed to the selected voxels, enlarged by xyz.margin
                  cannot exceed the initial volume.
                  Argument such as T.MAT, or for deprecated arguments
```

nesting.cube 113

Details

If obj is of class "volume", xyz.margin represents the distances in mm to be added to the x, y and z directions of the rectangular parallelepiped circumscribed to the voxels selected in sel.bin, in the obj frame of reference.

If obj is of class "mesh", sel.bin will undergo a margin expansion xyz.margin before the mesh points are selected..

Value

Returns a "volume" class object, in which 3D volume is limited to the rectangular parallelepiped circumscribed to the voxels selected by sel.bin, increased by the requested margins.

See Also

add.margin, nesting.cube and nesting.roi.

Examples

nesting.cube

Restriction of a volume to a rectangular parallelepiped

Description

The nesting. cube function restricts or increases a volume to the rectangular parallelepiped defined by its 2 extreme vertices.

Usage

```
nesting.cube(obj, pt.min, pt.max, alias = "", description = NULL, ...)
```

114 nesting.roi

Arguments

obj	object of class volume or mesh.
pt.min	minimum x, y, z coordinates of the rectangular parallelepiped vertex.
pt.max	maximum x, y, z coordinates of the rectangular parallelepiped vertex.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object. If the description = NULL (default value), it will be set to obj\$description.
	Additional arguments vol (depracated), replaced by obj.

Value

Returns a "volume" class object, in which 3D volume is restricted or increased to be circumscribed to the requested rectangular parallelepiped. If the created volume exceeds the initial volume, new voxels are set to NA.

See Also

add.margin, nesting.roi and nesting.bin.

Examples

nesting.roi

Restrict volume to RoI

Description

The nesting.roi function restricts a "volume" class object to the rectangular parallelepiped circumscribed to the chosen RoI.

nesting.roi 115

Usage

```
nesting.roi(
  obj,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  xyz.margin = c(0, 0, 0),
  obj.restrict = FALSE,
  T.MAT = NULL,
  alias = "",
  description = NULL,
  ...
)
```

Arguments

obj	object of class volume or mesh.
struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Index of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
xyz.margin	Vector of length 3, representing the distances in mm to be added to the x, y and z directions of the rectangular parallelepiped circumscribed to the chosen RoI, in the cutting planes frame of reference. By default xyz.margin = $c(0, 0, 0)$.
obj.restrict	Boolean. Used if obj is of class"volume". If obj.restrict = TRUE, the rectangular parallelepiped circumscribed to the selected voxels, enlarged by xyz.margin cannot exceed the initial volume.
T.MAT	"t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to obj\$ref.pseudo.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object. If description = NULL, it will be that of the obj, plus "restricted to" the selected RoI.
	Additional arguments such as vol (depracated), replaced by obj.

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all RoI are selected.

Value

Returns a "volume" class object, in which 3D volume is limited to the rectangular parallelepiped circumscribed to the chosen RoI, increased by the requested margins.

obj.create

See Also

add.margin, nesting.cube and nesting.bin.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct", "rtstruct"),</pre>
                              roi.name = "brain", dxyz = rep (step, 3))
CT <- patient$ct[[1]]
S <- patient$rtstruct[[1]]</pre>
CT.brain <- nesting.roi (CT, S, roi.sname = "brain")</pre>
CT.brain.with.margin <- nesting.roi (CT, S, roi.sname = "brain",
                                          xyz.margin = c (10,10,10))
\mbox{\tt\#} display at the center of gravity of the cerebellum Gz
Gz <- S$roi.info [grep("^brain",S$roi.info$roi.pseudo),]$Gz</pre>
display.plane (bottom = CT.brain, view.coord = Gz,
               struct = S, bottom.col = pal.RVV (1000),
               bottom.breaks = seq (-1000, 1000, length.out = 1001),
               bg = "#00ff00", interpolate = FALSE, legend.shift = -20)
display.plane (bottom = CT.brain.with.margin,view.coord = Gz,
               struct = S, bottom.col = pal.RVV (1000),
               bottom.breaks = seq(-1000, 1000, length.out = 1001),
               bg = "#00ff00", interpolate = FALSE, legend.shift = -20)
```

obj.create

Espadon object creating

Description

The obj. create function creates an espadon object with the essential properties it must have.

Usage

```
obj.create(class = c("", "volume", "struct", "mesh"), alias = "")
```

Arguments

class Character string, representing an espadon class from among "volume", "struct" or "mesh".

alias Character string, \$alias of the created object.

Value

Returns a espadon class object (see espadon.class for class definitions).

orientation.create 117

Examples

```
# Creation of an espadon mesh of a cube
M <- obj.create (class = "mesh")
M$mesh <- Rvcg::vcgIsotropicRemeshing (Rvcg::vcgBox(),0.5)
M$nb.faces <- ncol (M$mesh$it)
rgl::wire3d (M$mesh)</pre>
```

orientation.create

Creation of orientation

Description

The orientation.create function creates the orientation vectors of a plane:

- from 3 points A, B and C (see details),
- or from 2 vectors B and C, resp. defining x and y-axis (see details),
- or from 2 points A, B defining x-axis, and the normal vector to the plane (see details),
- or from a vector B defining x-axis, and the normal vector to the plane (see details).

Usage

```
orientation.create(A = c(0, 0, 0), B = NULL, C = NULL, normal = NULL)
```

Arguments

A	Vector of the x, y and z coordinates of point A, by default equal to $c(0,0,0)$ in the case where B and C are vectors.
В	Vector of x, y and z coordinates of point or vector B.
С	Vector of x, y and z coordinates of point or vector C.
normal	Vector of x, y and z coordinates of normal vector.
-	4

Details

When using B and C, B-A define the x-axis unit vector. The unit vector of the y-axis is orthonormal to the x-axis, coplanar with A, B and C, and in the direction of A to C.

When using B and normal, the unit vector of the x-axis is orthonormal to the normal vector, in the direction of A to B. The unit vector of the y-axis is defined so as to constitute a direct orthonormal basis with the unit vector of the x-axis and the normal vector of the plane.

Value

Returns the orientation of the plane. That means the concatenation of 2 vectors, defining an orthonormal basis of the plane.

pal.rainbow

Examples

```
A <- c (-29.93, 18.85, 4.34)
B <- c (28.73, 15.36, 4.46)
C <- c (1.53, 75.21, 13.51)
orientation.create (A, B, C)
```

pal.rainbow

Rainbow palette

Description

The pal.rainbow function produces a color palette adapted to dose repesentation.

Usage

```
pal.rainbow(n, alpha = seq(0.8, 0, length.out = n))
```

Arguments

n Integer, number of colors to be in the palette

alpha Vector representing the opacity, in the range of 0 (transparent) to 1 (opaque). If

alpha = NULL, all colors are opaque.

Value

Returns a color-labeled vector of size n.

```
pal <- pal.rainbow (255)

image (x = seq (0, 70, length.out = 255), y = 1,
    z = matrix (seq (0, 70, length.out = 255), ncol = 1),
    col = pal,
    main = "Rainbow colors")</pre>
```

pal.RVV

pal.RVV	Conversion of Hounsfied Units to Realistic Volume Vizualization col-
	ors

Description

The RVV.pal function produces a color palette where Hounsfield Units in the range -1000 HU to 1000 HU are converted into realistically colorized virtual anatomy (for use with CT), developed by *J.C. Silverstein and al* [1]

Usage

```
pal.RVV(n, alpha = NULL, min.col = "#000000", max.col = "#FFFFFF")
```

Arguments

n	Integer, number of colors to be in the palette
alpha	Vector representing the opacity, in the range of 0 (transparent) to 1 (opaque). If alpha = NULL (default), all colors are opaque, and no alpha channel is added to the colors.
min.col, max.col	l
	respectively the color below -1000HU (by default, black, i.e. "#000000") and above +1000HU (by default, white, i.e. "#FFFFFF")

Value

Returns a color-labeled vector of size n.

References

[1] Silverstein JC, Parsad NM, Tsirline V (2008). "Automatic perceptual color map generation for realistic volume visualization." *Journal of Biomedical Informatics*, **41**(6), 927-935. ISSN 1532-0464, doi:10.1016/j.jbi.2008.02.008.

```
pal <- pal.RVV (256)

image (x = seq (-1000, 1000, length.out = 1024), y = 1,
    z = matrix (seq (-1100, 1100, length.out = 1024), ncol = 1),
    col = pal,
    main = "Realistic Volume Vizualization colors")</pre>
```

plot plot

plot

plot a 2D cut of a 3D object

Description

The plot function displays the requested map of espadon objects of class "volume", "struct", "mesh".

Usage

```
## S3 method for class 'volume'
plot(
 Х,
  . . . ,
  view.type = "trans",
  view.coord = NULL,
  flip = FALSE,
  flop = FALSE,
  cut.interpolate = TRUE,
  display.interpolate = FALSE,
  col = grey.colors(255, start = 0, end = 1),
 breaks = NULL,
  sat.transp = FALSE
)
## S3 method for class 'struct'
plot(
 Х,
 view.type = "trans",
  view.coord = NULL,
  flip = FALSE,
  flop = FALSE,
  roi.name = NULL,
  roi.sname = NULL,
 roi.idx = NULL,
 back.dxyz = c(0.5, 0.5, x$thickness)
)
## S3 method for class 'mesh'
plot(
 Х,
 view.type = "trans",
  view.coord = NULL,
  flip = FALSE,
  flop = FALSE,
```

plot 121

```
col = "#ff0000"
)
```

Arguments

X	object of class "volume", "struct" or "mesh". See espadon.class
	others parameters of plot functions. See details
view.type	character string among the values 'ij', 'ji', 'ik', 'ki', 'jk', 'kj', 'yx', 'xz', 'zx', 'yz', 'zy', 'trans', 'front' or 'sagi" representing the map to be displayed.
view.coord	value representing the coordinate where the map is displayed. This parameter can be a 3D-vector, representing the coordinate of the point on the displayed map. If NULL, the display is located in the center of the object.
flip	Boolean defaults to FALSE flipping the horizontal axis of the background image.
flop	Boolean defaults to FALSE flipping the vertical axis of the background image.
cut.interpolat	e
	Boolean, indicating whether to calculate the volume cut using linear interpolation.
display.interp	
	Boolean, indicating whether to apply linear interpolation when displaying the cut.
col	Vector, representing the color palette of the image, if x is of class 'volume'. Color of the mesh outline if object x is of class 'mesh'.
breaks	One of:
	• NULL : the minimum and the maximum value of the object x define the range.
	• Vector giving the breakpoints of each color. Outside values are transparent, leaving the background visible, depending on sat.transp.
sat.transp	Boolean. If TRUE, outside values are transparent, else set to breaks limits colors.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
back.dxyz	3D vector of voxel size, used to calculate contours in frontal or sagittal view.

Details

 \dots can be xlim, ylim, add, bg etc. If view.type is egal to 'trans' or 'front' or 'sagi', the direction of xlim and ylim is ignored.

If roi.name, roi.sname, and roi.idx are all set to NULL, all closed planar or point RoI, present at view.coord are selected.

plot plot

Value

Returns a display of the k^{th} image plane of x.

See Also

display.plane, display.kplane, display.palette, pal.RVV, pal.rainbow.

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("ct","mr", "rtdose", "rtstruct"),</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
MR <- patient$mr[[1]]</pre>
CT <- patient$ct[[1]]
D <- patient$rtdose[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
# display 1
layout (matrix(c(1,1,2,3), ncol=2), widths=c(1,0.2))
plot (CT, view.coord = 0, col = pal.RVV(255))
S_plot <- plot (S, view.coord = 0, add = TRUE, lwd = 2)
display.palette (col = pal.RVV(255), main="HU")
display.legend (S_plot, bg="white", text.col="black", lwd = 2, cex = 1.1)
# display 2
layout (matrix(c (1, 1, 2, 3), ncol = 2), widths = c <math>(1, 0.2))
# Coordinates of the PTV barycenter in CT frame of reference
G <- as.numeric (S$roi.info[S$roi.info$roi.pseudo == "ptv",</pre>
                             c ("Gx", "Gy", "Gz")])
# Coordinates of the PTV baricenter in MR frame of reference
G_MR \leftarrow as.numeric (c(G, 1) %*%
                   t(get.rigid.M(CT$ref.pseudo, MR$ref.pseudo,
                                T.MAT = patient$T.MAT)))[1:3]
plot (MR, view.type = "sagi", view.coord = G_MR,
      col = grey.colors (255, start = 0, end = 1),
      breaks = seq (0, 500, length.out = 256), bg = "darkblue")
plot (vol.in.new.ref(D, MR$ref.pseudo, T.MAT = patient$T.MAT),
      view.type = "sagi", view.coord = G_MR,
      col = pal.rainbow(255), add = TRUE)
display.palette (col = grey.colors (255, start = 0, end = 1),
                 breaks = seq (0, 500, length.out = 256), main="MR")
display.palette (col = pal.rainbow(255),
                 breaks = seq (D$min.pixel, D$max.pixel, length.out = 256),
                 main="Gy")
layout(1)
```

Rdcm.inventory 123

Rdcm.inventory

Inventory of espadon objects from Rdcm files

Description

The Rdcm. inventory function creates, from Rdcm files in a patient's directory, a dataframe describing objects.

Usage

```
Rdcm.inventory(dirname, upgrade.to.latest.version = FALSE)
```

Arguments

dirname

Character string, representing the full name of patient directory, including Rdcm files.

upgrade.to.latest.version

Boolean. If TRUE, the function attempts to upgrade to the latest version, parsing the DICOM data. It may take longer to load the data. Consider using the Rdcm.upgrade function.

Value

Returns a dataframe, providing information of DICOM objects.

Rdcm.upgrade

Rdcm.upgrade

Updating Rdcm files.

Description

The Rdcm. upgrade function updates Rdcm files that were created with a previous version.

Usage

```
Rdcm.upgrade(Rdcm.files)
```

Arguments

Rdcm.files

String vector, representing the list of the full names of the Rdcm files, or its directories.

Value

Saves the updated Rdcm files. If the Rdcm files were generated from the dicom files, the data is updated from the DICOM fields.

```
# First, save toy patient objects to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dir.create (pat.dir, recursive = TRUE)
patient <- toy.load.patient (modality = c("ct", "mr"), roi.name = "",</pre>
                             dxyz = c (4, 4, 4)
save.to.Rdcm (patient$ct[[1]], dirname = pat.dir)
save.to.Rdcm (patient$mr[[1]], dirname = pat.dir)
save.T.MAT (patient$T.MAT, dirname = pat.dir)
# Rdcm files in pat.dir
list.files(pat.dir)
# test of Rdcm.upgrade
Rdcm.upgrade (pat.dir)
# or
Rdcm.upgrade (list.files (pat.dir, full.names = TRUE))
# Cleaning temporary directories
unlink (pat.dir, recursive = TRUE)
```

ref.add 125

ref.add

Adding a frame of reference in T.MAT

Description

The ref. add function adds the transfer matrices from or to a new frame of reference defined from 2 unit vectors and an origin point.

Usage

```
ref.add(
    src.ref,
    orientation = c(1, 0, 0, 0, 1, 0),
    origin = c(0, 0, 0),
    new.ref.pseudo = "newref",
    T.MAT = NULL
)
```

Arguments

src.ref	Character string, pseudonym of the frame of reference in which the orientation vector and the origin point origin are defined.
orientation	Vector of 6 or 9 elements, composed of the coordinates of the 2 orthonormal vectors (i, j), or of the 3 orthonormal vectors (i, j, k) of the new coordinate system, in the src.ref frame of reference.
origin	Vector of the x, y, z coordinates of the origin point of the new frame of reference in the src.ref frame of reference. Default to $c(0, 0, 0)$.
new.ref.pseudo	Character string, pseudonym of the new frame of reference to add.
T.MAT	"t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, then only the link between src.ref and new.ref.pseudo is computed.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to new.ref.pseudo pseudonym of the new frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-new.ref.pseudo", "new.ref.pseudo<- new.ref.pseudo", "new.ref.pseudo<- src.ref".

Returns a NULL if orientation is not well defined.

See Also

ref.cutplane.add, ref.remove, ref.srctodest.add.

126 ref.cutplane.add

Examples

ref.cutplane.add

Adding volume's cutting planes frame of reference in T.MAT

Description

The ref.cutplane.add function adds in T.MAT the transfer matrices from or to volume's cutting planes frame of reference.

Usage

```
ref.cutplane.add(
  vol,
  origin = c(0, 0, 0),
  ref.cutplane = paste0(vol$ref.pseudo, "m"),
  T.MAT = NULL
)
```

Arguments

vol "volume" class object.

origin Vector of the x, y, z coordinates of the origin point of the cutting planes frame

of reference.

ref.cutplane Name of the volume's cutting planes frame of reference. By default ref.cutplane

= paste0 (vol\$ref.pseudo,".m").

T.MAT "t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm

or load.T.MAT. If T.MAT = NULL, then only the link between vol\$ref.pseudo

and ref. name is established.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to volume's cutting planes frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-ref.cutplane", "ref.cutplane", "ref.cutplane", "ref.cutplane".

See Also

ref.add, ref.srctodest.add, ref.remove.

ref.remove 127

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = "mr", roi.name = "", dxyz = c (4, 4, 4))
MR <- patient$mr[[1]]
MR$xyz.from.ijk

# creation of t.mat, containing the transfer matrix to the frame of reference
# of the MR cutting planes
t.mat <- ref.cutplane.add (MR)

# Change of frame of reference
MR.m <- vol.in.new.ref (MR, paste0 (MR$ref.pseudo, "m"), t.mat)

MR.m$xyz.from.ijk</pre>
```

ref.remove

Deletion of a frame of reference in T.MAT

Description

The ref. remove function removes the management of a frame of reference in T.MAT.

Usage

```
ref.remove(ref.name, T.MAT)
```

Arguments

ref.name Character string, pseudonym of the frame of reference to delete.

T.MAT "t.mat" class object in which the ref. name frame of reference is to be deleted.

Value

Returns a "t.mat" class object, which no longer contains transfer matrices from or to the ref.pseudo ref.name. ref.cutplane.add.

128 ref.srctodest.add

ref.srctodest.add

Linking two existing frames of reference in T.MAT

Description

The ref.srctodest.add function links the source frame of reference with the destination frame of reference.

Usage

```
ref.srctodest.add(src.ref, dest.ref, TM = diag(4), T.MAT = NULL)
```

Arguments

src.ref	Character string, pseudonym of the source frame of reference.
dest.ref	Character string, pseudonym of the destination frame of reference.
TM	4x4 tansfert matrix for moving from src.ref to dest.ref.
T.MAT	"t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, then only the link between src.ref and dest.ref is established.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to dest.ref pseudonym of the new frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-dest.ref", "dest.ref<- dest.ref", "dest.ref<-.src.ref".

See Also

ref.add, ref.cutplane.add, ref.remove.

rt.chi.index 129

rt.chi.index

Chi index 2D - 3D

Description

The rt.chi.index function computes the local or global Chi index from a measurement and a reference. These latter are "volume" class objects containing one (2D) or several planes (3D).

Usage

```
rt.chi.index(
  vol,
  vol.ref,
  abs = TRUE,
  vol.max = vol.ref$max.pixel,
  dose.th = 0.02,
  delta.r = 3,
  analysis.th = 0.05,
  local = FALSE,
  local.th = 0.3,
  project.to.isocenter = TRUE,
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object, which represents the measured volume.
vol.ref	"volume" class object, which represents the reference volume.
abs	Boolean. If TRUE (default), the absolute value of Chi is computed.
vol.max	Positive number, by default equal to the maximum value of the reference volume. See Details.
dose.th	Number between 0 and 1, used to determine the dose difference criterion. See Details.
delta.r	Positive number, in mm. Distance difference criterion.
analysis.th	Number between 0 and 1. Only the voxels whose value are greater than or equal analyse.th \star vol.max are processed.
local	Boolean. If local = FALSE (default), a global Chi index is computed, and a local Chi index otherwise.
local.th	Number between 0 and 1. Local threshold, only used if local = TRUE. See Details.
project.to.isocenter	

Boolean. If TRUE, and if vol and vol.ref are of modality "rtimage", the size of the pixels is corrected to correspond to that found if the sensor was at the isocenter.

130 rt.chi.index

alias Character string, \$object.alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to Chi index setup.

Details

The Chi index of a voxel n was defined by Bakai and al [1]. It is computed from the formulae:

$$\chi_n = \frac{D_i - Dref_n}{\sqrt{\Delta D^2 + \Delta r^2 \cdot \|\nabla Dref_n\|^2}}$$

If abs = TRUE, the used formulae is:

$$\chi_n = \frac{|D_i - Dref_n|}{\sqrt{\Delta D^2 + \Delta r^2 \cdot \|\nabla Dref_n\|^2}}$$

with D_i the measured dose at voxel i, $Dref_n$ the reference dose at voxel n, $\nabla Dref_n$ the gradient of reference dose at voxel n, Δr the distance difference criterion equal to delta.r, and ΔD the distance difference criterion at voxel n defined as follows:

- If local = FALSE a global Chi index is computed and $\Delta D = dose.th \cdot vol.max$.
- If local = TRUE, then $\Delta D = dose.th \cdot Dref_n$ when $Dref_n \geq local.th \cdot vol.max$, and $\Delta D = dose.th \cdot local.th \cdot vol.max$ otherwise.

Value

Returns a "volume" class object (see espadon.class for class definitions). The \$vol3D.data field represents the Chi index. Two fields are added: the \$setup field recalls the calculation setup, and the \$chi.info field details the number of dose points, the number of evaluated dose points, the rate of evaluated dose points, the rate of absolute values of the Chi index below 1, above 1.2 and 1.5, the max and the mean Chi index.

References

[1] Bakai A, Alber A, Nüsslin F (2003). "A revision of the Gamma-evaluation concept for the comparison of dose distributions." *Physics in Medicine and Biologys*, **48**(21), 3543–3553.

See Also

rt.gamma.index

rt.gamma.index 131

rt.gamma.index

Gamma index 2D - 3D

Description

The rt.gamma.index function computes the local or global Gamma index from a measurement and a reference. These latter are "volume" class objects containing one (2D) or several planes (3D).

Usage

```
rt.gamma.index(
  vol,
  vol.ref,
  over.sampling.factor = 1,
  vol.max = vol.ref$max.pixel,
  dose.th = 0.02,
  delta.r = 3,
  analysis.th = 0.05,
  local = FALSE,
  local.th = 0.3,
  project.to.isocenter = TRUE,
  alias = "",
  description = NULL
)
```

Arguments

```
vol "volume" class object, which represents the measured volume.
vol.ref "volume" class object, which represents the reference volume.
over.sampling.factor
```

Strictly positive integer, or a vector of 3 strictly positive integers, default to 1. Defined to oversample grids of vol and vol.ref. Oversampling can be very time consuming.

rt.gamma.index

vol.max	Positive number, by default equal to the maximum value of the reference volume. See Details.	
dose.th	Number between 0 and 1, used to determine the dose difference criterion. See Details.	
delta.r	Positive number, in mm. Distance difference criterion.	
analysis.th	Number between 0 and 1. Only the voxels whose value is greater than or equal analyse.th*vol.max are processed.	
local	Boolean. If local = FALSE (default), a global Gamma index is computed, and a local Gamma index otherwise.	
local.th	Number between 0 and 1. Local threshold, only used if local = TRUE. See Details.	
project.to.isocenter		
	Boolean. If TRUE, and if vol and vol.ref are of modality "rtimage", the size of the pixels is corrected to correspond to that found if the sensor was at the isocenter.	
alias	Character string, \$object.alias of the created object.	
description	Character string, describing the created object. If description = NULL (default value), it will be set to Gamma index setup.	

Details

The Gamma index of a voxel n was defined by Low and al [1]. It is computed from the formulae:

$$\gamma_n = min\left(\sqrt{\frac{(D_i - Dref_n)^2}{\Delta D^2} + \frac{r_i^2}{\Delta r^2}}\right)$$

whith D_i the measured dose at voxel i, $Dref_n$ the reference dose at voxel n, r_i the distance between voxels i and n, Δr the distance difference criterion equal to delta.r, ΔD the distance difference criterion at voxel n defined as follows:

- If local = FALSE a global Gamma index is computed and $\Delta D = dose.th \cdot vol.max$.
- If local = TRUE, then $\Delta D = dose.th \cdot Dref_n$ when $Dref_n \geq local.th \cdot vol.max$, and $\Delta D = dose.th \cdot local.th \cdot vol.max$ otherwise.

Value

Returns a "volume" class object (see espadon.class for class definitions). The \$vol3D.data field represents the Gamma index. Two fields are added: the \$setup field recalls the calculation setup, and the \$gamma.info field details the number of dose points, the number of evaluated dose points, the rate of evaluated dose points, the rate of Gamma indices below 1, above 1.2 and 1.5, the max and the mean Gamma index.

References

[1] Low DA, Harms WB, Mutic S, Purdy JA (1998). "A technique for the quantitative evaluation of dose distributions." *Medical Physics*, **25**(5), 656–661.

rt.indices.from.bin

See Also

rt.chi.index

Examples

```
# Creation of a reference volume and measured volume
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c ("rtdose", "rtstruct"),</pre>
                              roi.name = "ptv", dxyz = c (3, 3, 3))
D.ref <- patient$rtdose[[1]]</pre>
# We will assume that the measured dose is equal to the reference dose shifted
# by 3 pixels on the x axis
D.meas <- vol.copy (D.ref, alias = "measured_dose")</pre>
D.meas$vol3D.data[1:(D.meas$n.ijk[1] - 3) ,,] <- D.ref$vol3D.data[4:D.ref$n.ijk[1],,]</pre>
D.max <- as.numeric(quantile(as.numeric(D.ref$vol3D.data),</pre>
                              probs = 99.99/100, na.rm = TRUE))
gamma <- rt.gamma.index (D.meas, D.ref, delta.r = 6, vol.max = D.max)</pre>
gamma$gamma.info
# Display gamma index at isocenter
G.iso <- patient$rtstruct[[1]]$roi.info$Gz[</pre>
  patient$rtstruct[[1]]$roi.info$name == "ptv"]
display.plane(gamma, view.coord = G.iso,
              bottom.col = c ("#00FF00", "#007F00", "#FF8000", "#FF0000",
                               "#AF0000"),
              bottom.breaks = c(0, 0.8, 1, 1.2, 1.5, gamma*max.pixel),
              bg = "blue", interpolate = FALSE)
```

rt.indices.from.bin Dosimetry, volume, conformity, homogeneity indices from binary selection

Description

The rt.indices.from.bin function calculates, from a "volume" class object of modality "rtdose", all the standard indicators of radiotherapy, as long as their options are transmitted, for the target and healthy "volume" object of modality "binary".

Usage

```
rt.indices.from.bin(
  vol,
  target.bin.list = NULL,
  healthy.bin.list = NULL,
  T.MAT = NULL,
  presc.dose = NA,
  healthy.tol.dose = NA,
  healthy.weight = 1,
  dosimetry = c("D.min", "D.max", "D.mean", "STD"),
```

134 rt.indices.from.bin

```
volume.indices = c("V.tot", "area", "V.prescdose"),
     conformity.indices = c("PITV", "CI.lomax2003", "CN", "NCI", "DSC", "CI.distance",
       "CI.abs_distance", "CDI", "CS3", "ULF", "OHTF", "gCI", "COIN", "G_COSI", "COSI"),
      homogeneity.indices = c("HI.RTOG.max_ref", "HI.RTOG.5_95", "HI.ICRU.max_min",
         "HI.ICRU.2.98_ref", "HI.ICRU.2.98_50", "HI.ICRU.5.95_ref", "HI.mayo2010",
         "HI.heufelder"),
      gradient.indices = c("GI.ratio.50"),
      D.xpc = NULL,
      D.xcc = NULL,
      V.xpc = NULL,
      V.xGy = NULL,
      DVH = TRUE,
      verbose = TRUE,
    )
Arguments
    vol
                      "volume" class object, of "rtdose" modality.
    target.bin.list
                     list of "volume" class objects, of "binary" modality. The $object.alias field
                      of each target.bin.list object represents the name of the selected region of
                      interest of the target volume.
    healthy.bin.list
                     list of "volume" class objects, of "binary" modality. The $object.alias field
                      of each healthy.bin.list object represents the name of the selected region of
                     interest of the healthy tissues.
    T.MAT
                      "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T. MAT
                     = NULL, all $ref.pseudo of bin.list elements must be equal to vol$ref.pseudo.
                     vector of prescription doses that serve as reference doses for the target RoI.
    presc.dose
    healthy.tol.dose
                     vector of tolerance dose of each healthy RoI.
    healthy.weight Vector of weight, indicating the importance of the healthy RoI.
                      Vector indicating the requested dose indicators from among 'D.min', 'D.max',
    dosimetry
                      'D.mean' and 'STD'. If D. xpc is different from NULL, it will be added.
    volume.indices Vector indicating the requested volume indices from among 'V.tot', 'V.prescdose'
                      (i.e. volume over presc.dose) and 'area'. If V.xGy is different from NULL, it
                      will be added.
    conformity.indices
                      Vector. Requested conformity indices from among 'PITV', 'PDS', 'CI.lomax2003',
                      'CN', 'NCI', 'DSC', 'CI.distance', 'CI.abs_distance', 'CDI', 'CS3', 'ULF',
                      'OHTF', 'gCI', 'COIN', 'COSI' and 'G_COSI'.
    homogeneity.indices
                      Vector. Requested homogeneity indices from among 'HI.RTOG.max_ref', 'HI.RTOG.5_95',
                      'HI.ICRU.max_min', 'HI.ICRU.2.98_ref', 'HI.ICRU.2.98_50', 'HI.ICRU.5.95_ref',
                      'HI.mayo2010' and 'HI.heufelder.'
```

rt.indices.from.bin

gradient.indices	
	Vector. Requested gradient indices from among 'GI.ratio.50', 'mGI'.
D.xpc	Vector of the percentage of the volume, for which the dose coverage is requested.
D.xcc	Vector of the volume in cm^3 , for which the dose coverage is requested.
V.xpc	Vector of the percentage of the reference dose, received by the volume to be calculated.
V.xGy	Vector of the minimum dose in Gy, received by the volume to be calculated.
DVH	boolean. If TRUE (default), the function returns the DVHs of the target and healthy ROIs.
verbose	Boolean. if TRUE (default) a progress bar is displayed.
	others parameters such as DVH.step.

Value

Return a list of indices dataframe. For details, see rt.indices.from.roi.

See Also

rt.indices.from.roi.

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for better
# result)
step <- 5
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "eye",</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
struct <- patient$rtstruct[[1]]</pre>
T.MAT <- patient$T.MAT
# creation of the list of target binary volumes
target.roi.idx <- select.names (struct$roi.info$roi.pseudo, roi.sname = "ptv")</pre>
healthy.roi.idx <- select.names (struct$roi.info$roi.pseudo, roi.sname = "eye")
target.bin.list <- lapply (target.roi.idx , function (idx) {</pre>
  vr <- nesting.roi (D, struct, roi.idx = idx, xyz.margin = c (5, 5, 5),</pre>
                     T.MAT = T.MAT, alias = struct$roi.info$name[idx])
  b <- bin.from.roi(vr, struct, roi.idx = idx, T.MAT = T.MAT,
                    alias = struct$roi.info$name[idx], verbose = FALSE)
  })
names (target.bin.list) <- struct$roi.info$name[target.roi.idx]</pre>
healthy.bin.list <- lapply (healthy.roi.idx , function (idx) {
  vr <- nesting.roi (D, struct, roi.idx = idx, xyz.margin = c (5, 5, 5),
                     T.MAT = T.MAT, alias = struct$roi.info$name[idx])
  b <- bin.from.roi(vr, struct, roi.idx = idx, T.MAT = T.MAT,
                    alias = struct$roi.info$name[idx], verbose = FALSE)
})
names (healthy.bin.list) <- struct$roi.info$name[healthy.roi.idx]</pre>
```

rt.indices.from.roi Dosimetry, volume, conformity, homogeneity indices from RoI

Description

The rt.indices.from.roi function calculates, from a "volume" class object of modality "rtdose", standard indicators of radiotherapy in relation to the target and healthy RoI, as long as their options are transmitted.

Usage

```
rt.indices.from.roi(
  vol,
  struct = NULL,
 T.MAT = NULL
  target.roi.name = NULL,
  target.roi.sname = NULL,
  target.roi.idx = NULL,
  healthy.roi.name = NULL,
  healthy.roi.sname = NULL,
  healthy.roi.idx = NULL,
  presc.dose = NA,
  healthy.tol.dose = NA,
  healthy.weight = 1,
  dosimetry = c("D.min", "D.max", "D.mean", "STD"),
  volume.indices = c("V.tot", "area", "V.prescdose"),
  conformity.indices = c("PITV", "PDS", "CI.lomax2003", "CN", "NCI", "DSC",
   "CI.distance", "CI.abs_distance", "CDI", "CS3", "ULF", "OHTF", "gCI", "COIN",
    "G_COSI", "COSI"),
  homogeneity.indices = c("HI.RTOG.max_ref", "HI.RTOG.5_95", "HI.ICRU.max_min",
    "HI.ICRU.2.98_ref", "HI.ICRU.2.98_50", "HI.ICRU.5.95_ref", "HI.mayo2010",
    "HI.heufelder"),
  gradient.indices = c("GI.ratio.50", "mGI"),
 D.xpc = NULL,
 D.xcc = NULL,
 V.xpc = NULL
 V.xGy = NULL,
 DVH = TRUE,
 verbose = TRUE,
)
```

Arguments

vol "volume" class object, of "rtdose" modality.

struct "struct" class object.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT

= NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.

target.roi.name

Exact name of target RoI in struct object. By default target.roi.name =

NULL. See Details.

target.roi.sname

Name or part of name of target RoI in struct object. By default target.roi.sname = NULL. See Details.

target.roi.idx Value of the index of target RoI that belong to the struct object. By default target.roi.idx = NULL. See Details.

healthy.roi.name

Exact name of healthy RoI in struct object. By default healthy.roi.name = NULL.

healthy.roi.sname

Name or part of name of healthy RoI in struct object. By default healthy.roi.sname - NIII I

healthy.roi.idx

Value of the index of healthy RoI that belong to the struct object. By default healthy.roi.idx = NULL.

presc.dose Vector of prescription doses that serve as reference doses for the target RoI.

healthy.tol.dose

Vector of tolerance doses of each healthy RoI.

healthy.weight Vector of weights, indicating the importance of the healthy RoI.

dosimetry Vector indicating the requested dose indicators from among 'D.min', 'D.max', 'D.mean' and 'STD'. If D.xpc is different from NULL, it will be added.

volume.indices Vector indicating the requested volume indices from among 'V.tot', 'V.prescdose' (i.e. volume over presc.dose) and 'area'. If V.xGy is different from NULL, it will be added.

conformity.indices

Vector. Requested conformity indices from among 'PITV', 'PDS', 'CI.lomax2003', 'CN', 'NCI', 'DSC', 'CI.distance', 'CI.abs_distance', 'CDI', 'CS3', 'ULF', 'OHTF', 'gCI', 'COIN', 'COSI' and 'G_COSI'.

homogeneity.indices

Vector. Requested homogeneity indices from among 'HI.RTOG.max_ref', 'HI.RTOG.5_95', 'HI.ICRU.max_min', 'HI.ICRU.2.98_ref', 'HI.ICRU.2.98_50', 'HI.ICRU.5.95_ref', 'HI.mayo2010' and 'HI.heufelder.'

gradient.indices

Vector. Requested gradient indices from among 'GI.ratio.50', 'mGI'.

D. xpc Vector of the percentage of the volume, for which the dose coverage is requested.

D. xcc Vector of the volume in cm^3 , for which the dose coverage is requested.

V.xpc	Vector of the percentage of the reference dose, received by the volume to be calculated.
V.xGy	Vector of the minimum dose in Gy, received by the volume to be calculated.
DVH	boolean. If TRUE (default), the function returns the DVHs of the target and healthy ROIs.
verbose	Boolean. if TRUE (default) a progress bar is displayed.
	others parameters such as DVH.step.

Details

If target.roi.name, target.roi.sname, and target.roi.idx are all set to NULL, all RoI containing 'ptv' (if they exist) are selected.

If target.roi.name, target.roi.sname, and target.roi.idx are all set to NULL,no target RoI are selected.

If healthy.roi.name, healthy.roi.sname, and healthy.roi.idx are all set to NULL, no healthy RoI are selected.

Value

Returns a list containing (if requested)

- -dosimetry: dataframe containing, for all target and healthy structures:
 - the requested dosimetry: D.min (Gy), D.max (Gy), D.mean (Gy) and STD (Gy), respectively the minimum, maximum, mean and standard deviation of the dose in the regions of interest.
 - the requested \$D.x%: (Gy) Dose covering x percent of structure volume.
 - the requested D.xcc: (Gy) Dose covering x (cm^3) of structure volume.

All numbers are displayed with a resolution depending on the DVH quantization step. The default resolution is 1e-3.

- -volume: dataframe containing, for all target and healthy structures, and isodoses:
 - the requested volume.indices: $V_{tot}(cm^3)$ (except for isodose) the total volume of the regions of interest, area (cm^2) (except for isodose) their surface areas, $V_{tot}(cm^3)$ the volumes receiving at least presc.dose Gy,
 - the requested V. xGy (cm^3) : volumes receiving at least x Gy.
 - the requested V. xpc (cm^3) Volume receiving at least x% of the reference dose.
- -conformity: dataframe containing, if requested,
 - PITV: (1) Prescription Isodose Target Volume, or conformity index defined by E.Shaw [1]

$$PITV = \frac{V_{presc.dose}}{V_{target}}$$

• PDS: (1) Prescription Dose Spillage defined by SABR UK Consortium 2019 [2]

$$PDS = \frac{V_{presc.dose}}{V_{target \ \geq \ presc.dose}} = \frac{V_{presc.dose}}{V_{target} \ \cap \ V_{presc.dose}}$$

• CI.lomax2003: (1) Conformity Index defined by Lomax and al [3]

$$CI_{lomax2003} = \frac{V_{target \geq presc.dose}}{V_{presc.dose}} = \frac{V_{target} \cap V_{presc.dose}}{V_{presc.dose}}$$

• CN: (1) Conformation Number defined by *Van't Riet and al* [4]. It corresponds to conformity index defined by *Paddick* [5]

$$CN = CI_{paddick2000} = \frac{V_{target \geq presc.dose}^2}{V_{target} \cdot V_{presc.dose}} = \frac{(V_{target} \cap V_{presc.dose})^2}{V_{target} \cdot V_{presc.dose}}$$

• NCI: (1) New conformity index, inverse of CN, defined by Nakamura and al [6]

$$NCI = \frac{1}{CN}$$

• DSC: (1) Dice Similarity Coefficient [7]

$$DSC = 2 \cdot \frac{V_{target \geq presc.dose}}{V_{target} + V_{presc.dose}} = 2 \cdot \frac{V_{target} \cap V_{presc.dose}}{V_{target} + V_{presc.dose}}$$

• CI.distance: (1) Conformity Index based on distance defined by Park and al [8]

$$CI.distance = \frac{100}{N} \sum_{i=1}^{N} \frac{dist_{S_{presc.dose}} \rightarrow G_{target} - dist_{S_{target}} \rightarrow G_{target}}{dist_{S_{target}} \rightarrow G_{target}}$$

where $dist_{S_{presc.dose} \to G_{target}}$ is the distance between the surface of the prescription dose volume and the centroid of the target, and $dist_{S_{target} \to G_{target}}$ the surface of the target volume and the centroid of the target. N is the number of directions where the distances are calculated. These directions are computed every 1° . If the centroid is not within the target volume, then CI.distance = NA.

• CI.abs_distance: (1) Conformity Index based on distance defined by Park and al [8]

$$CI.abs_distance = \frac{100}{N} \sum_{i=1}^{N} \frac{|dist_{S_{presc.dose}} \rightarrow G_{target} - dist_{S_{target}} \rightarrow G_{target}|}{dist_{S_{target}} \rightarrow G_{target}|}$$

• CDI: (1) Conformity Distance Index defined by Wu and al [9]

$$CDI = 2\frac{V_{presc.dose} + V_{target} - 2 \ V_{target} \ge presc.dose}{S_{target} + S_{presc.dose}} = \frac{V_{presc.dose} + V_{target} - 2 \ \cdot \ V_{target} \ \cap \ V_{presc.dose}}{S_{target} + S_{presc.dose}}$$

where S_{target} is the surface of the target volume and $S_{presc.dose}$ is the surface of the prescription dose volume.

• CS3: (1) Triple Point Conformity Scale defined by Ansari and al [10]

$$CS3 = \frac{V_{0.95 \cdot presc.dose} + V_{presc.dose} + V_{1.05 \cdot presc.dose}}{3 \cdot V_{target}}$$

• ULF: (1) Underdosed lesion factor defined by Lefkopoulos and al [11]

$$ULF = \frac{V_{target} < presc.dose}{V_{target}} = \frac{V_{target} \cap \overline{V_{presc.dose}}}{V_{target}}$$

• OHTF:(1) Overdosed healthy tissues factor defined by Lefkopoulos and al [11]

$$OHTF = \frac{\sum V_{healthy \ \geq \ presc.dose}}{V_{target}} = \frac{\sum V_{healthy} \ \cap \ V_{presc.dose}}{V_{target}}$$

• gCI: (1) Geometric Conformity Index defined by Lefkopoulos and al [11]

$$qCI = ULF + OHTF$$

• COIN: Conformity Index defined by Baltas and al [12]

$$COIN = \frac{V_{target \geq presc.dose}^2}{V_{target} \cdot V_{presc.dose}} \cdot \prod^{N_{healthy}} \left(1 - \frac{V_{healthy \geq presc.dose}}{V_{healthy}}\right)$$

• gCOSI : generalized COSI, defined by *Menhel and al* [13].

$$gCOSI = 1 - \sum_{healthy}^{N_{healthy}} healthy.weight \cdot \frac{\frac{V_{healthy} \geq healthy.tol.dose}{V_{healthy}}}{\frac{V_{target} \geq presc.dose}{V_{target}}}$$

-COSI: if "COSI" is requested in conformity.indices, it returns a dataframe of Critical Organ Scoring Index for each healthy organ, at each presc.dose, and for each target. COSI is defined by *Menhel and al* [13]

$$COSI = 1 - \frac{\frac{V_{healthy \geq healthy, tol.dose}}{V_{healthy}}}{\frac{V_{target} \geq presc.dose}{V_{target}}}$$

- homogeneity: dataframe containing
 - HI.RTOG.max_ref: (1) Homogeneity Index from RTOG defined by E.Shaw [1]

$$HI.RTOG.max_ref = \frac{D_{max}}{presc.dose}$$

where D_{max} is the maximum dose in the target volume.

• HI.RTOG.5_95: (1) Homogeneity Index from RTOG [1]

$$HI.RTOG.5_95 = \frac{D.5pc}{D.95pc}$$

where D.5pc and D.95pc are respectively the doses at 5% and 95% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.max_min: (1) Homogeneity Index from ICRU report 62 [14]

$$HI.ICRU.max_min = \frac{D_{max}}{D_{min}}$$

where D_{max} and D_{min} are respectively the maximum and the minimum dose in the target volume.

• HI.ICRU.2.98_ref: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.2.98_ref = 100 \frac{D.2pc - D.98pc}{presc.dose}$$

where D.2pc and D.98pc are respectively the doses at 2% and 98% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.2.98_50: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.2.98_50 = 100 \frac{D.2pc - D.98pc}{D.50pc}$$

where D.2pc, D.98pc and D.50pc are respectively the doses at 2%, 98% and 50% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.5.95_ref: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.5.95_ref = 100 \frac{D.5pc - D.95pc}{presc.dose}$$

where D.5pc and D.95pc are respectively the doses at 5% and 95% of the target volume in cumulative dose-volume histogram.

• HI.mayo2010: (1) Homogeneity Index defined by Mayo and al [16]

$$HI.mayo2010 = \sqrt{\frac{D_{max}}{presc.dose} \cdot (1 + \frac{\sigma_D}{presc.dose})}$$

where D_{max} is the maximum dose in the target volume, and σ_D the standard deviation of the dose in the target volume.

• HI. heufelder: (1) Homogeneity Index defined by Heufelder and al [17]

$$HI.heufelder = e^{-0.01 \cdot (1 - \frac{\mu_D}{presc.dose})^2} \, \cdot \, e^{-0.01 \cdot (\frac{\sigma_D}{presc.dose})^2}$$

where μ_D and σ_D are respectively the mean and the standard deviation of the dose in the target volume.

- gradient : dataframe containing
 - GI. ratio. 50: Gradient Index based on volumes ratio defined by *Paddick and Lippitz* [18]

$$GI.ratio.50 = \frac{V_{0.5 \cdot presc.dose}}{V_{presc.dose}}$$

• mGI: Modified Gradient Index defined by SABR UK Consortium 2019 [2]

$$mGI = \frac{V_{0.5 \ \cdot \ presc. dose}}{V_{target} \geq presc. dose} = \frac{V_{0.5 \ \cdot \ presc. dose}}{V_{target} \ \cap \ V_{presc. dose}}$$

- DVH: dataframe containing the Dose-Volume histograms of target and healthy ROIs, in percent.

References

- [1] Shaw E, Kline R, Gillin M, Souhami L, Hirschfeld A, Dinapoli R, Martin L (1993). "Radiation therapy oncology group: Radiosurgery quality assurance guidelines." *International journal of radiation oncology, biology, physics*, **27**(5), 1231-1239. ISSN 0360-3016, doi:10.1016/0360-3016(93)90548A.
- [2] UK SABR Consortium (Online; accessed 2022-04-01). "Stereotactic Ablative Radiation Therapy (SABR): a resource. v6.1, January 2019." https://www.sabr.org.uk/wp-content/uploads/2019/04/SABRconsortium-guidelines-2019-v6.1.0.pdf.

[3] Lomax NJ, Scheib SG (2003). "Quantifying the degree of conformity in radiosurgery treatment planning." *International journal of radiation oncology, biology, physics*, **55**(5), 1409-1419. ISSN 0360-3016, doi:10.1016/S03603016(02)045996.

- [4] Riet AV, Mak AC, Moerland MA, Elders LH, Van der Zee W (1997). "A conformation number to quantify the degree of conformality in brachytherapy and external beam irradiation: Application to the prostate." *International journal of radiation oncology, biology, physics*, **37**(3), 731-736. ISSN 0360-3016, doi:10.1016/S03603016(96)006013.
- [5] Paddick I (2000). "A simple scoring ratio to index the conformity of radiosurgical treatment plans. Technical note." *Journal of neurosurgery*, **93 Suppl 3**, 219-222.
- [6] Nakamura J, Verhey L, Smith V, Petti P, Lamborn K, Larson D, Wara W, Mcdermott M, Sneed P (2002). "Dose conformity of Gamma Knife radiosurgery and risk factors for complications." *International journal of radiation oncology, biology, physics*, **51**, 1313-9. doi:10.1016/S0360-3016(01)017576.
- [7] Dice LR (1945). "Measures of the Amount of Ecologic Association Between Species." *Ecology*, **26**(3), 297–302. ISSN 00129658, 19399170.
- [8] Park JM, Park S, Ye S, Kim J, Carlson J, Wu H (2014). "New conformity indices based on the calculation of distances between the target volume and the volume of reference isodose." *The British journal of radiology*, **87**, 20140342. doi:10.1259/bjr.20140342.
- [9] Wu Q, Wessels BW, Einstein DB, Maciunas RJ, Kim EY, Kinsella TJ (2003). "Quality of coverage: Conformity measures for stereotactic radiosurgery." *Journal of Applied Clinical Medical Physics*, **4**, 374-381.
- [10] Ansari S, Satpathy S, Singh P, Lad S, Thappa N, Singh B (2018). "A new index: Triple Point Conformity Scale (CS3) and its implication in qualitative evaluation of radiotherapy plan." *Journal of Radiotherapy in Practice*, 17, 1-4. doi:10.1017/S1460396917000772.
- [11] Lefkopoulos D, Dejean C, balaa ZE, Platoni K, Grandjean P, Foulquier J, Schlienger M (2000). "Determination of dose-volumes parameters to characterise the conformity of stereotactic treatment plans." In chapter XIII, 356-358. Springer Berlin Heidelberg. ISBN 978-3-540-67176-3, doi:10.1007/9783642597589 135.
- [12] Baltas D, Kolotas C, Geramani KN, Mould RF, Ioannidis G, Kekchidi M, Zamboglou N (1998). "A conformal index (COIN) to evaluate implant quality and dose specification in brachytherapy." *International journal of radiation oncology, biology, physics*, **40 2**, 515-24. doi:10.1016/s0360-3016(97)007323.
- [13] Menhel J, Levin D, Alezra D, Symon Z, Pfeffer R (2006). "Assessing the quality of conformal treatment planning: a new tool for quantitative comparison." *Physics in Medicine and Biology*, **51**(20), 5363–5375.
- [14] Landberg T, Chavaudra J, Dobbs J, Gerard J, Hanks G, Horiot J, Johansson K, Möller T, Purdy J, Suntharalingam N, Svensson H (1999). "ICRU Report 62: Prescribing, Recording and Reporting Photon Beam Therapy (Supplement to ICRU Report 50),3. Absorbed Doses." *Reports of the International Commission on Radiation Units and Measurements*, os-32(1), 21-25.
- [15] ICRU (2010). "Report 83: Prescribing, Recording, and Reporting Photon-Beam Intensity-Modulated Radiation Therapy (IMRT)." *Reports of the International Commission on Radiation Units and Measurements*, **10**(1), 1-3.
- [16] Mayo CS, Ding L, Addesa A, Kadish S, Fitzgerald TJ, Moser R (2010). "Initial Experience With Volumetric IMRT (RapidArc) for Intracranial Stereotactic Radiosurgery." *International Jour-*

save.T.MAT

nal of Radiation Oncology*Biology*Physics, **78**(5), 1457-1466. ISSN 0360-3016, doi:10.1016/j.ijrobp.2009.10.005.

[17] Heufelder J, Zink K, Scholz M, Kramer K, Welker K (2003). "Eine Methode zur automatisierten Bewertung von CT-basierten Bestrahlungsplanen in der perkutanen Strahlentherapie." *Zeitschrift fur Medizinische Physik*, 13(4), 231-239. ISSN 0939-3889, doi:10.1078/0939388900175.

[18] Paddick I, Lippitz BE (2006). "A simple dose gradient measurement tool to complement the conformity index." *Journal of neurosurgery*, 105 Suppl, 194-201.

All this references are compiled by

- Kaplan LP, Korreman SS (2021). "A systematically compiled set of quantitative metrics to describe spatial characteristics of radiotherapy dose distributions and aid in treatment planning." *Physica Medica*, **90**, 164-175. ISSN 1120-1797, doi:10.1016/j.ejmp.2021.09.014. and
- Patel G, Mandal A, Choudhary S, Mishra R, Shende R (2020). "Plan evaluation indices: A journey of evolution." *Reports of Practical Oncology & Radiotherapy*, 25. doi:10.1016/j.rpor.2020.03.002..

See Also

rt.indices.from.bin.

Examples

save.T.MAT

Save a T.MAT class object

Description

The save.T.MAT function saves the data required by load.T.MAT, load.patient.from.dicom or load.patient.from.Rdcm to generate T.MAT, as pre-formatted Rdcm files.

Usage

```
save.T.MAT(T.MAT, dirname)
```

144 save.to.Rdcm

Arguments

T.MAT "t.mat" class object to save.

dirname Directory where new reg .Rdcm files will be saved.

Details

Reg files from DICOM files cannot be updated with the save. T. MAT function. Only transfer matrices added with ref.add or ref.cutplane.add will be saved.

Value

Returns TRUE, if all reg files generating T.MAT are saved.

Examples

save.to.Rdcm

Save a espadon object in a pre-formatted *.Rdcm file

Description

The function save.to.Rdcm allows you to save an object created by **espadon** in a pre-formatted *.Rdcm file. This object will also be accessible by the load.patient.from.Rdcm function.

Usage

```
save.to.Rdcm(obj, object.name = obj$object.alias, dirname = obj$file.dirname)
```

Arguments

```
obj espadon object of class "volume", "struct", "mesh", "histo", "dvh", "histo2D".

Object.name Character string, representing the name of the object, default to obj$object.alias.

Directory where new files from obj will be saved.
```

select.names 145

Value

Returns TRUE, if paste0(object.name, ".Rdcm") exists in dirname.

Returns FALSE, if object.name is not a valid file name, or if the file that is created would replace a *.Rdcm file created by dicom.to.Rdcm.converter.

Note

save.to.Rdcm can not replace an *.Rdcm file created by dicom.to.Rdcm.converter.

Examples

select.names

Regions of Interest (RoI) indices

Description

The select.names function allows you to select words from a vector of words, according to several criteria, eliminating spaces and case.

Usage

```
select.names(names, roi.name = NULL, roi.sname = NULL, roi.idx = NULL)
```

Arguments

names	Words vector
roi.name	Vector of words to compare to names. By default roi.name = NULL. See Details
roi.sname	Vector of words or parts of words to compare. By default roi.sname = NULL. See Details
roi.idx	Index vector. By default roi.idx = NULL. See Details.

set.reference.obj

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all RoI are selected.

Value

Returns the indices of the elements of the word vector names satisfying one or more of the following conditions:

- ASCII // TRANSLIT transcriptions, without spaces, of names and roi.name, are identical.
- ASCII // TRANSLIT transcriptions, without spaces of roi.sname are identical to part of ASCII // TRANSLIT transcriptions, without spaces of names.
- names indices belong to the index vector roi.idx.

Examples

```
# loading patient objects
names <- c ("Eye left", "EyeR", "OPTICAL nerve L", "opical nervR", "chiasma")
# RoI selection.
select.names (names = names, roi.name = c("eye left", "eye right"))
select.names (names = names, roi.sname = c("eye", "ner"))
select.names (names = names, roi.idx = 4:9)</pre>
```

set.reference.obj

Set the reference objects of a espadon object

Description

The function set.reference.obj adds to an espadon object the information identifying the espadon objects from which it derives.

Usage

```
set.reference.obj(obj, ref.obj, add = TRUE)
```

Arguments

```
espadon object of class "dvh", "fan", "histo", "histo2D", "mesh", "rtplan", "struct", 
"undef" or "volume".

ref.obj espadon object of class "dvh", "fan", "histo", "histo2D", "mesh", "rtplan", "struct", 
"undef" or "volume". List of espadon objects.

add Boolean. If TRUE, the reference objects are added to those already contained 
by obj.
```

Value

Returns the espadon object obj, containing the ref.object.alias and ref.object.info fields identifying its reference objects

sp.similarity.from.bin 147

Examples

sp.similarity.from.bin

Volume-based spatial similarity metrics calculated from binary modality 3D volumes.

Description

The sp.similarity.from.bin function computes volumetric Dice similarity coefficient, Dice-Jaccard coefficient and Dice surface similarity coefficient.

Usage

```
sp.similarity.from.bin(
  vol.A,
  vol.B,
  coeff = c("DSC", "DJC", "MDC", "under.MDC", "over.MDC")
```

Arguments

vol.A, vol.B "volume" class objects, of "binary" modality. vol.B is the reference for MDC calculation.

coeff

Vector indicating the requested metrics from among 'DSC' (Dice similarity coefficient), 'DJC' (Dice-Jaccard coefficient), and 'MDC' (mean distance to conformity). Equal to NULL if not requested.

Value

returns a dataframe containing (if requested):

• volumetric Dice similarity coefficient DSC defined by :

$$DSC = 2\frac{V_A \cap V_B}{V_A + V_B}$$

• Dice-Jaccard coefficient DJC defined by :

$$DJC = \frac{V_A \ \cap \ V_B}{V_A \ \cup \ V_B}$$

• mean distance to conformity MDC, over-contouring mean distance over . MDC and under-contouring mean distance under . MDC, defined by *Jena et al* [1]

References

[1] Jena R, et al. (2010). "A novel algorithm for the morphometric assessment of radiotherapy treatment planning volumes." *Br J Radiol.*, **83**(985), 44-51. doi:10.1259/bjr/27674581.

See Also

sp.similarity.from.mesh

Examples

```
sp.similarity.from.mesh
```

Distance-based spatial similarity metrics calculated from the mesh.

Description

The sp.similarity.from.mesh function computes Hausdorff distances and surface Dice similarity coefficient.

Usage

```
sp.similarity.from.mesh(
  mesh1,
  mesh2,
  hausdorff.coeff = c("HD.max", "HD.mean"),
  hausdorff.quantile = c(0.5, 0.95),
  surface.tol = seq(0, 10, 1)
)
```

sp.similarity.from.mesh 149

Arguments

```
mesh1, mesh2 espadon mesh class objects

hausdorff.coeff

Vector indicating the requested Hausdorff distance metrics from among 'HD.max', 'HD.mean'.

Equal to NULL if not requested. NULL, it will be added.

hausdorff.quantile

numeric vector of probabilities with values between 0 and 1, representing the quantiles of the unsigned distances between mesh1 and mesh2. Equal to NULL if not requested.

surface.tol

numeric vector representing the maximum margins of deviation which may be tolerated without penalty. Equal to NULL if not requested.
```

Value

Returns a list containing (if requested):

- Hausdorff: dataframe including the maximum, mean and quantiles
- smetrics: dataframe with the columns:
 - tol: the requested tolerances
 - sDSC: the surface Dice similarity coefficients, defined by Nikolov et al [1]
 - sAPL: the surface Added Path Length in mm^2 , introduced (in pixels) by Vaassen et al [2]

References

- [1] Nikolov S, et al. (2018). "Deep learning to achieve clinically applicable segmentation of head and neck anatomy for radiotherapy." *ArXiv*, **abs/1809.04430**.
- [2] Vaassen S, et al. (2020). "Evaluation of measures for assessing time-saving of automatic organat-risk segmentation in radiotherapy." *Physics and Imaging in Radiation Oncology*, **13**, 1-6.

See Also

sp.similarity.from.bin

```
library (Rvcg)
# espadon mesh of two spheres of radius R1=10 and R2=11, separated by dR = 3
sph <- vcgSphere ()
mesh1 <- obj.create ("mesh")
mesh1$nb.faces <- ncol (sph$it)
mesh1$mesh <- sph
mesh2 <- mesh1</pre>
R1 <- 10
R2 <- 11
dR <- 3
mesh1$mesh$vb[1:3,] <- R1 * mesh1$mesh$normals[1:3,] + mesh1$mesh$vb[1:3,]
mesh2$mesh$vb[1:3,] <- R2 * mesh2$mesh$normals[1:3,] + mesh2$mesh$vb[1:3,] +</pre>
```

150 struct.clustering

struct.clustering

Clustering volumes by RoI

Description

The struct.clustering function creates a new volume in which voxels are clustered and labeled by region of interest defined in an rt-struct.

Usage

```
struct.clustering(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  alias = "",
  description = NULL,
  verbose = TRUE
)
```

Arguments

vol	"volume" class object.
struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (struct\$object.alias, "clustering")
verbose	Boolean. if TRUE (default), the RoI studied are listed.

struct.from.bin 151

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all RoI are selected.

Value

Returns a "volume" class object (see espadon.class for class definitions), of "cluster" modality. This object contains the custor = cust

See Also

get.roi.connection

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 5
patient <- toy.load.patient (modality = c ("mr", "rtstruct"),</pre>
                              dxyz = rep (step, 3))
MR <- patient$mr[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
cluster.vol <- struct.clustering (MR, S, T.MAT = patient$T.MAT, verbose = FALSE)</pre>
head (cluster.vol$cluster.info)
# Display
n = nrow(cluster.vol$cluster.info)
col = paste0(c ("#000000", rainbow (n-1)), "af")
breaks <- seq (cluster.vol$min.pixel - 0.5, cluster.vol$max.pixel + 0.5,</pre>
              length.out = n+1)
par0 <- par()
par(mfrow = c(1,2), mar = c(1,15,1,1))
display.palette(col, factors = cluster.vol$cluster.info$label)
par(mar = c(1, 1, 1, 1))
display.plane (MR, cluster.vol, main = "RoI clustering", view.coord = 0,
                top.col = col, top.breaks = breaks, interpolate = FALSE)
par(mfrow= par0$mfrow,mar=par0$mar)
```

struct.from.bin

Creation of struct class object from a binary volume

Description

The struct.from.bin function creates a struct object with a unique RoI, defined by the contours of binary volume.

152 struct.from.bin

Usage

Arguments

vol	"volume" class object, of binary modality.
roi.name	Character string, representing the name of created RoI.
roi.nb	Positive integer, representing the number of created RoI.
roi.color	Color of the created RoI, in hex code format ("#RRGGBB").
roi.type	Type of RoI, from among "", "EXTERNAL", "PTV", "CTV", "GTV", "TREATED_VOLUME", "IRRAD_VOLUME", "OAR", "BOLUS", "AVOIDANCE", "ORGAN", "MARKER", "REGISTRATION", "ISOCENTER", "CONTRAST_AGENT", "CAVITY", "BRACHY_CHANNEL" "BRACHY_ACCESSORY", "BRACHY_SRC_APP", "BRACHY_CHNL_SHLD", "SUPPORT", "FIXATION", "DOSE_REGION", "CONTROL" and "DOSE_MEASUREMENT"
external.only	Boolean. If TRUE, only external contours are kept.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object.

Value

Returns a "struct" class object (see espadon.class for class definition), including the unique roi. name as region of interest.

struct.from.mesh 153

```
S.sphere <- struct.from.bin (b.sphere, roi.name = "sphere", external.only = TRUE)
str (S.sphere$roi.info)</pre>
```

struct.from.mesh

Creation of struct class object from an espadon mesh

Description

The struct.from.mesh function creates a struct object with a unique RoI, defined by the contours of a mesh.

Usage

```
struct.from.mesh(
  mesh,
  z,
  thickness = NULL,
  roi.name = mesh$object.alias,
  roi.nb = 1,
  roi.color = "#ff0000",
  roi.type = "",
  alias = "",
  description = NULL,
  force.closed = TRUE,
  verbose = TRUE,
  ...
)
```

Others parameters

Arguments

espadon mesh class object.
z-coordinate vector where mesh contours are computed.
struct thickness between 2 adjacent contours. If NULL (default), it is deduced from z.
Character string, representing the name of created RoI.
Positive integer, representing the number of created RoI.
Color of the created RoI, in hex code format ("#RRGGBB").
Type of RoI, from among "", "EXTERNAL", "PTV", "CTV", "GTV", "TREATED_VOLUME", "IRRAD_VOLUME", "OAR", "BOLUS", "AVOIDANCE", "ORGAN", "MARKER", "REGISTRATION", "ISOCENTER", "CONTRAST_AGENT", "CAVITY", "BRACHY_CHANNEL", "BRACHY_ACCESSORY", "BRACHY_SRC_APP", "BRACHY_CHNL_SHLD", "SUPPORT", "FIXATION", "DOSE_REGION", "CONTROL" and "DOSE_MEASUREMENT".
Character string, \$alias of the created object.
Character string, describing the the created object.
Boolean. Set to TRUE when the mesh represents the surface of a closed volume.
Boolean. If TRUE (default), a progress bar indicates the state of calculation.

struct.in.new.ref

Value

Returns a "struct" class object (see espadon.class for class definition), including the unique roi.name as region of interest.

Examples

```
# Creation of an espadon mesh of a cube
M <- obj.create (class = "mesh")
M$mesh <- Rvcg::vcgIsotropicRemeshing (Rvcg::vcgBox(),0.5)
M$nb.faces <- ncol (M$mesh$it)

S <- struct.from.mesh (M, z = seq(-1,1,0.5))
display.3D.contour(S)</pre>
```

struct.in.new.ref

Change of frame of reference of a "struct" class object.

Description

The struct.in.new.ref function allows you to change the frame of reference of a struct.

Usage

```
struct.in.new.ref(struct, new.ref.pseudo, T.MAT, alias = "")
```

Arguments

struct	"struct" class object.
new.ref.pseudo	pseudonym of the frame of reference in which the struct should be located. This new.ref.pseudo must exist in the T.MAT list.
T.MAT	$"t.mat"\ class\ object,\ created\ by\ load.patient.from.dicom,\ load.patient.from.Rdcm,\ load.T.MAT\ or\ ref.add.$
alias	Character string, \$alias of the created object.

Value

Returns "struct" class object in the new frame of reference new.ref.pseudo.

See Also

vol.in.new.ref

struct.merge 155

Examples

struct.merge

Merging of structures into a new structure

Description

The struct.merge function merges two structures into a new one. It is useful for comparing contours, for example.

Usage

```
struct.merge(
  ref.struct,
  add.struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  suffix = "",
  alias = "",
  description = ""
```

Arguments

ref.struct	struct class object. All RoI of this structure are kept.
add.struct	struct class object. Only the selected RoI are kept for merging.
roi.name	$\label{thm:condition} Vector of exact names of the RoI in the \verb add.struct object. By default \verb roi.name = \verb NULL . See Details .$
roi.sname	Vector of names or parts of names of the RoI in the add.struct object. By default ${\tt roi.sname}$ = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the add.struct object. By default roi.idx = NULL. See Details.
suffix	Character string. '-suffix' is added to RoI name.
alias	Character string, \$alias of the resulted object.
description	Character string, describing the the resulted object.

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all RoI of add.struct are selected.

156 study.deployment

Value

Returns a struct class object. See espadon.class for class definitions.

Note

Beware that, when merging structures, some RoI may have same name or roi.info\$roi.pseudo. In this case struct.merge prints a warning message. Consider changing suffix to avoid the ambiguity.

See Also

struct.from.bin.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
patient <- toy.load.patient (modality = c("rtdose"),</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
# isodose 50% Dmax Gy and 90% Dmax
bin50 <- bin.from.vol (D, min = 0.5 * D$max.pixel)</pre>
bin90 <- bin.from.vol (D, min = 0.9 * D$max.pixel)</pre>
S.isodose50 <- struct.from.bin (bin50, roi.name = "50pc" ,</pre>
                                 roi.color = "#00FFFF")
S.isodose90 <- struct.from.bin (bin90, roi.name = "90pc" ,</pre>
                                 roi.color = "#FFFF00")
S <- struct.merge (S.isodose50, S.isodose90, alias = "isodose",
                    description = paste ("isodose of", D$object.alias))
# Dmax location :
z.dmax <- get.xyz.from.index(which (D$vol3D.data == D$max.pixel), D)[1,3]</pre>
display.plane(top = D, struct = S, view.coord = z.dmax, legend.shift = -50)
```

study.deployment

Deployment of DICOM files from multiple patients

Description

The study.deployment function duplicates DICOM data from multiple patients, so that it becomes data independent of the original data. This function simplifies the analysis of multi-center or multi-expert studies in dosimetry challenges, contouring consensus searches, etc.

study.deployment 157

Usage

```
study.deployment(
  pats.dir,
  deploy.dir,
  design.matrix = matrix(TRUE, nrow = length(dir(pats.dir)), ncol = 1, dimnames =
    list(basename(dir(pats.dir)), "expert_1")),
  pid.prefix = "",
  white.list = c("instance", "reference"),
  black.list = c("frame of reference", "class"),
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

pats.dir	Name of the directory in which all patient directories are stored, each containing the DICOM files to be duplicated.
deploy.dir	Name of the directory where all patient files will be duplicated.
design.matrix	Boolean matrix. See Details.
pid.prefix	string vector of length 1 or string vector of length ncol(design.matrix), representing the prefix added to the new unique identifier of the deployed patient (tag $(0010,0020)$).
white.list	Names vector, representing a part of the DICOM tag name UI value representation, other than those defined by the DICOM standard, which will be modified. By default, the UID name containing 'instance' or 'reference' will be modified.
black.list	Names vector, representing a part of the DICOM tag name UI value representation, other than those defined by the DICOM standard, which will not be modified. By default, the frame of reference UID will not be modified.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Details

The design.matrix argument defines how patients DICOM files will be deployed. The names of the lines must match the names of the directories contained in pats.dir. The names of the columns are for example the different experts or hospitals who will study the patient files. These experts will only review the patients files defined by rownames(design.matrix)[design.matrix[,"expert"]].

Value

Creates the deploy.dir directory, containing the expert directories defined by the design.matrix column names. Each expert directory contains as many patient directories as defined by the design.matrix row names. All patients will be independent of each other. The new created patients have the pats.dir as name, and expert name as first name, and an independent patient ID, with prefix pid.prefix.

158 toy.dicom.raw

Examples

```
# First, save toy.dicom.raw () raw data to a temporary file/pats.dir/toy_PM
# for testing.
toy_PM.dir <- file.path (tempdir(), "pats.dir","toy_PM")</pre>
dir.create (toy_PM.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "toyrtplan", tmpdir = toy_PM.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
# function test:
pats.dir <- dirname (toy_PM.dir)</pre>
deploy.dir <- file.path (tempdir(), "deploy.dir")</pre>
design.matrix <- matrix(TRUE, nrow = length (dir (pats.dir)), ncol=3,</pre>
                         dimnames = list (basename (dir (pats.dir)),
                                           c("Dr Quinn", "Dr Who", "Dr House")))
design.matrix
study.deployment (pats.dir, deploy.dir, design.matrix,
                  pid.prefix = c("zz_", "yy_", "xx_"))
# check result
list.files(deploy.dir, recursive = TRUE)
load.patient.from.dicom(deploy.dir)$patient
# Cleaning temporary directory
unlink (pats.dir, recursive = TRUE)
unlink (deploy.dir, recursive = TRUE)
```

toy.dicom.raw

toy DICOM raw data

Description

The toy.dicom.raw loads raw data from a dummy DICOM file. It is used for the test.

Usage

```
toy.dicom.raw()
```

Value

Returns the raw data of a dummy DICOM file of rtplan modality.

```
toy.dicom.raw ()
```

toy.load.patient 159

toy.load.patient

Load a toy patient for test

Description

The toy.load.patient creates a dummy "patient" class object. It is used for the test.

Usage

```
toy.load.patient(
  modality = c("ct", "mr", "rtdose", "rtstruct"),
  roi.name = c("eye", "optical nerve", "brain", "labyrinth processing unit",
        "energy unit", "gizzard", "ghost container", "exhaust valve"),
  dxyz = c(1, 1, 1),
  beam.nb = 7
)
```

Arguments

modality	String vector, whose elements are chosen among the modalities "ct", "mr", "rt-struct" and "rtdose".
roi.name	String vector, whose elements are chosen among the regions of interest (RoI) "eye", "optical nerve", "brain", "labyrinth processing unit", "energy unit", "gizzard", "ghost container" and "exhaust valve". Note that the RoI "couch", "patient" and "ptv" are still present.
dxyz	Vector of length 3, representing the x, y, z steps in mm, between ct, mr and rtdose voxels.
beam.nb	Positive integer. Number of radiotherapy beams in rtdose modality.

Value

Returns an toy object of "patient" class, containing the modalities defined in modality. See espadon.class for class definitions.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 5
pat <- toy.load.patient (dxyz = rep (step, 3), beam.nb = 2)
str (pat, max.level = 2)</pre>
```

160 vol.copy

vector	.product
vector	. broduct

Vector product of two vectors

Description

Vector product of two vectors

Usage

```
vector.product(v1, v2)
```

Arguments

```
v1 Vector of x, y, z coordinates
v2 Vector of x, y, z coordinates
```

Value

Returns the x, y, z coordinates of the vector product of v1 and v2

Examples

```
vector.product(c (1, 0, 0), c (0, 1, 0))
```

vol.copy

Creating a volume from another one

Description

The vol.copy function creates a "volume" class object, with the same grid as the vol volume object.

Usage

```
vol.copy(vol, alias = "", modality = NULL, description = NULL, number = NULL)
```

Arguments

vol	"volume" class object, template of the created object.
alias	Character string, \$object.alias of the created object.

modality Character string, modality of the created volume. If modality = NULL, then the

created object will have the modality of vol.

description Character string, description of the returned object. If description = NULL,

then the created object will have the description of vol.

number number of the returned volume. If number = NULL, then the returned object will

have the number of vol.

vol.create 161

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data is initialized to NA.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = "ct", roi.name = "",dxyz = c (4, 4, 4))
CT <- patient$ct[[1]]

# creating a volume
vol.from.CT <- vol.copy (CT, alias = "ct reference")
str (vol.from.CT)</pre>
```

vol.create

Volume creating

Description

The vol.create function creates a volume object from a user-defined grid.

Usage

```
vol.create(
  n.ijk,
  dxyz,
  mid.pt = NULL,
  pt000 = NULL,
  default.value = NA,
  value.sd = 0,
  ref.pseudo = "ref1",
  frame.of.reference = "",
  alias = "",
  modality = "",
  description = "",
  number = 0
)
```

Arguments

n.ijk	Vector of length 3, representing the number of elements on the 1, 1 and k axes.
dxyz	Vector of length 3, representing the x , y , z steps in mm, between voxels. See details.
mid.pt	Vector of length 3, representing the x, y, z coordinates of the midpoint of the volume. See details.
pt000	Vector of length 3, representing the x, y, z coordinates of the first voxel of the first plane.

162 vol.from.bin

default.value Numerical or boolean value, representing the default value of the voxels.

value.sd Positive integer. If default.value is a number, it represents standard deviation

of the volume default value.

ref. pseudo Character string, frame of reference pseudonym of the created object. By defaukt

equal to "ref1"

frame.of.reference

Character string, frame of reference of the created object.

alias Character string, \$alias of the created object.

modality Character string, \$modality of the created object.

description Character string, describing the the created object.

number Integer, by default set to 0, number of the created object.

Details

If mid.pt and pt000 are both equal to NULL, then mid.pt = c(0, 0, 0) by default. If mid.pt and pt000 are both different from NULL, then only mid.pt is taken into account.

Value

Returns a "volume" class object (see espadon.class for class definitions), in which the grid is defined by pt000 or mid.pt, dxyz n.ijk. If default.value are initialized to FALSE, then modality = "binary". The orientation of the patient is orthonormal to the grid.

Examples

vol.from.bin

Volume class object according to binary selection

Description

The vol.from.bin function selects a part of a "volume" class object of "binary" modality which has the same grid. It is especially useful to restrict voxel data in region of interest.

Usage

```
vol.from.bin(vol, sel.bin, alias = "", description = NULL)
```

vol.gradient 163

Arguments

vol	"volume" class object, containing data to restrict.
sel.bin	"volume" class object, of "binary" modality. vol and sel.bin must have the same grid.
alias	Character string, \$alias of the created object
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol \$object, alias, "from", sel, bin \$object, alias)

Value

Returns a "volume" class object (see espadon.class for class definitions), in which non-selected voxels have the value NA, and selected voxels have the original value of vol.

Examples

vol.gradient

Gradient of a volume

Description

The vol.gradient function calculates the 3D gradient of a "volume" class object

Usage

```
vol.gradient(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "gradient").

vol.in.new.ref

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol, gradient of vol.

Examples

vol.in.new.ref

Change of frame of reference of a volume

Description

The vol.in.new.ref function allows you to change the frame of reference of a volume.

Usage

```
vol.in.new.ref(vol, new.ref.pseudo, T.MAT, alias = "", description = NULL)
```

Arguments

vol	"volume" class object.
new.ref.pseudo	pseudonym of the frame of reference in which the volume should be located. This new.ref.pseudo must exist in the T.MAT list.
T.MAT	"t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT or ref.add.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be that of the vol volume.

Value

Returns "volume" class object in the new frame of reference new.ref.pseudo.

See Also

struct.in.new.ref

vol.median 165

Examples

vol.median

Median filter on a volume

Description

The vol.median function applies a 26-connectivity median filter on all the voxels of a "volume" class object.

Usage

```
vol.median(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "median").

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol, in which voxels are filtered by a 26-connectivity median filter.

166 vol.oversampling

vol.oversampling

Oversampling a volume

Description

The vol. oversampling function oversamples the grid of a "volume" class object.

Usage

```
vol.oversampling(
  vol,
  fact.ijk = 2,
  alias = "",
  interpolate = TRUE,
  description = NULL
)
```

Arguments

vol "volume" class object.

fact.ijk Strictly positive integer, or a vector of 3 strictly positive integers.

alias Character string, \$alias of the created object.

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the

value of the voxels, relative to the values of adjacent voxels, is performed.

description Character string, describing the the created object. If description = NULL, it

will be paste ("oversampling", vol\$description).

Value

Returns a "volume" class object, in which 3D volume grid is oversampled: the voxel size is divided by fact.ijk.

See Also

vol.subsampling.

vol.regrid 167

vol.regrid

Transform the grid of a volume class object into the grid of another

Description

The vol. regrid function transforms the grid of a volume according to the grid of another.

Usage

```
vol.regrid(
  vol,
  back.vol,
  T.MAT = NULL,
  interpolate = TRUE,
  method = c("NN", "Av"),
  alias = "",
  description = NULL,
  verbose = TRUE
)
```

Arguments

vol	"volume" class object to regrid.
back.vol	"volume" class object whose grid will be used for regriding. Its \$ref.pseudo must exist in the T.MAT list.
T.MAT	"t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT or ref.add. If T.MAT = NULL, back.vol\$ref.pseudo must be equal to vol\$ref.pseudo.
interpolate	Boolean, default to TRUE.
method	method of interpolation, defaulted to 'NN' ie 'Nearest Neighbor'. See Details.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object. If description = NULL, it will be that of vol.
verbose	Boolean. if TRUE (default) a progress bar is displayed.

168 vol.repair

Details

The interpolation method is chosen from:

- 'NN': the value of a voxel is calculated from its nearest adajcent neighbors.
- 'Av': the value of a voxel is the weighted average of the voxels contained in a box, whose sides are automatically calculated from the back.bol\$dxyz steps.

Value

Returns a copy of vol, in which grid is that of back.vol.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c ("mr", "rtdose"),</pre>
                              dxyz = rep (step, 3), beam.nb = 4)
MR <- patient$mr[[1]]</pre>
D <- patient$rtdose[[1]]</pre>
# change grid
D.on.MR <- vol.regrid (vol = D, back.vol = MR, interpolate = TRUE,
                       T.MAT = patient$T.MAT, alias = "",
                        description = NULL, verbose = FALSE)
# maximum dose location
max.dose.in.MR <- get.xyz.from.index (which.max (D.on.MR$vol3D.data), D.on.MR)
display.plane (bottom = MR, view.coord = max.dose.in.MR[3],
               top= D.on.MR, bottom.col = grey.colors(255, start = 0, end = 1),
               bottom.breaks = seq (0, 500, length.out = 256),
               bg = "#00ff00", interpolate = FALSE)
```

vol.repair

repairing missing planes of volumes

Description

The vol. repair function repairs missing planes in volumes.

Usage

```
vol.repair(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "repair").

vol.subsampling 169

Details

Missing planes at download can generate errors or unpredictible results in espadon processing. The vol.repair function detects such missing planes and recreates their value by interpolation.

Value

Returns a "volume" class object (see espadon.class for class definitions), with no missing plane, if vol is to be repaired. Returns vol otherwise.

Examples

```
step <- 4
patient <- toy.load.patient (modality = c("ct", "mr", "rtstruct", "rtdose"),</pre>
                               roi.name = "",
                               dxyz = rep (step, 3), beam.nb = 3)
CT <- patient$ct[[1]]
# this function removes a plane in a volume in order to simulate
# a dicom transfer issue
remove.plane <- function (vol, k) {
  idx <- which (vol$k.idx == k)</pre>
  vol$n.ijk[3] \leftarrow vol$n.ijk[3] - 1
  vol$xyz0 <- vol$xyz0[-idx, ]</pre>
  vol$k.idx <- vol$k.idx[-idx]</pre>
  vol\mbox{missing.k.idx} \leftarrow TRUE
  vol$vol3D.data <- vol$vol3D.data[, , -idx]</pre>
  return (vol)
}
# Creation of CT.damaged without the 29th slice.
CT.damaged<- remove.plane (CT, 29)
CT.fix <- vol.repair (CT.damaged)</pre>
# Display
par (mfrow=c(3, 3))
for (k in 28:30) {
display.kplane (CT, k, main = paste("CT @ k =",k),interpolate = FALSE)
display.kplane (CT.damaged, k, main = "damaged CT",interpolate = FALSE)
display.kplane (CT.fix, k, main = "fixed CT", interpolate = FALSE)
}
```

vol.subsampling

Subsampling a volume

Description

The vol. subsampling function sub-samples the grid of a "volume" class object.

vol.subsampling

Usage

```
vol.subsampling(
  vol,
  fact.ijk = 2,
  interpolate = TRUE,
  method = c("NN", "Av"),
  alias = "",
  description = NULL
)
```

Arguments

vol "volume" class object.

fact.ijk Strictly positive integer, or a vector of 3 strictly positive integers.

interpolate Boolean, default to TRUE.

method method of interpolation, defaulted to 'NN' ie 'Nearest Neighbor'. See Details.

alias Character string, \$alias of the created object.

description Character string, describing the the created object. If description = NULL, it will be paste ("subsampling", vol\$description).

Details

The interpolation method is chosen from:

- 'NN': the value of a voxel is calculated from its nearest adajcent neighbors.
- 'Av': the value of a voxel is the weighted average of the voxels contained in a box, whose sides are automatically calculated from the back.bol\$dxyz steps.

Value

Returns a "volume" class object, in which 3D volume grid is subsampled: the voxel size is multiplied by fact.ijk and the center location of the volume is invariant.

See Also

vol.oversampling.

vol.sum 171

vol.sum

Sum of 2 volumes

Description

The vol. sum function adds two "volume" class objects of the same grid and of the same modality.

Usage

```
vol.sum(vol1, vol2, alias = "", description = NULL)
```

Arguments

vol1, vol2

"volume" class objects. The 2 volumes must have the same modality, and the same grid (i.e. share the same position of the voxels).

Character string, \$alias of the created object.

Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1\$object.alias, "+", vol2\$object.alias).

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol1 and vol2, sum of vol1 and vol2.

172 xlsx.from.dcm

xlsx.from.dcm

Converting DICOM files to .xlsx files

Description

The xlsx.from.dcm function creates an Excel file from DICOM files.

Usage

```
xlsx.from.dcm(
  dcm.filenames,
  xlsx.filenames,
  multipage = TRUE,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dcm.filenames	String vector, representing the list of full names of DICOM files.
xlsx.filenames	String vector, representing the list of full names of created *.xlsx files. If multipage = TRUE, only the xlsx.filenames[1] is used.
multipage	Boolean. If TRUE, all dcm.filenames are converted into multiple pages of the same *.xlsx file.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Value

Returns a boolean vector, establishing the existence of the created Excel files.

```
# First, save toy.dicom.raw () raw data to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "PM_dcm")
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PMrtplan", tmpdir = pat.dir, fileext = ".dcm")
zz <- file (dcm.filename, "wb")
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
list.files (pat.dir)
# Creating an Excel file</pre>
```

xlsx.from.Rdcm 173

xlsx.from.Rdcm

Converting .Rdcm files to .xlsx files

Description

A *.Rdcm file contains the list of contents, in dataframe form, of the DICOM files of the same object. The xlsx.from.Rdcm function creates, from a *.Rdcm file, an Excel file, in which each page contains the dataframe representation of a DICOM file of the same object.

Usage

```
xlsx.from.Rdcm(
  Rdcm.filenames,
  dest.dirname = dirname(Rdcm.filenames),
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

Rdcm.filenames	String vector, representing the *.Rdcm filenames to be converted.
dest.dirname	String vector of the same length as Rdcm.filenames, indicating the directory where the $*.xlsx$ files will be created.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Value

Returns a boolean vector, establishing the existence of the created Excel files which have the same basenames as the *.Rdcm files.

174 xlsx.from.Rdcm

```
# First, create a Rdcm file from toy.dicom.raw () to a temporary file for testing.
pat.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PM_rtplan", tmpdir = pat.dir, fileext = ".dcm")</pre>
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
dicom.to.Rdcm.converter (dcm.filename, pat.dir, update = TRUE)
file.remove (dcm.filename)
list.files (pat.dir)
# Creating an Excel file
Rdcm.filenames <- list.files (pat.dir, pattern = "[.]Rdcm$",</pre>
                               recursive = TRUE, full.names = TRUE)
xlsx.from.Rdcm (Rdcm.filenames)
list.files (pat.dir)
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

Index

```
add.margin, 4, 7, 12, 13, 19, 113, 114, 116
                                                   display.3D.stack,42
add. shape, 5, 10
                                                   display.dV_dx, 46, 48, 94, 96, 98
                                                   display.DVH, 43, 45, 93
bin.closing, 6, 12, 13, 19
                                                   display.DVH.pc, 44, 45, 93
bin.clustering, 7
                                                   display.histo, 47, 47, 94, 96, 98
bin.cuboid, 6, 8
                                                   display.kplane, 48, 58, 74, 122
bin.cylinder, 6
                                                   display.legend, 51
bin.cylinder(bin.cuboid), 8
                                                   display.obj.links, 52, 67, 78
bin.dilation, 7, 11, 13, 19
                                                   display.palette, 53, 122
bin.ellipsoid, 6
                                                   display.plane, 48, 50, 55, 122
bin.ellipsoid(bin.cuboid), 8
bin.erosion, 7, 12, 12, 19
                                                   espadon.class, 5, 7, 8, 10, 12-14, 17-21, 37,
bin.from.roi, 6, 13
                                                            38, 40, 41, 45, 46, 49, 58, 69–71,
bin.from.vol, 14, 15
                                                            100-102, 104, 116, 121, 130, 132,
bin.intersection, 17
                                                            151, 152, 154, 156, 159, 161–165,
bin.inversion, 18
                                                            169, 171
bin.opening, 7, 12, 13, 19
                                                   export, 66
bin.subtraction, 20
bin.sum, 21
                                                   fan.beam, 68, 70, 72, 73
                                                   fan.planar, 69, 69, 72, 73
castlow.str, 22, 23
                                                   fan.sphere, 69, 70, 70, 72, 73
castup.str, 22, 22
                                                   fan.to.voxel, 69, 70, 72, 72
dicom.browser, 23, 25, 29–31, 35
                                                   get.extreme.pt, 73
dicom.parser, 24, 31, 35
                                                   get.ijk.from.index, 74, 82
dicom.patient.anonymiser, 26
                                                   get.ijk.from.xyz,75
dicom.raw.data.anonymizer, 28
                                                   get.line, 76
dicom.raw.data.loader, 24, 25, 29
                                                   get.obj.connection, 53, 77
dicom.set.tag.value, 30
                                                   get.plane, 78
dicom.tag.dictionary, 23, 25, 27, 28, 30,
                                                   get.rigid.M, 80
         31, 34, 35, 67, 98, 99, 102, 157, 172,
                                                   get.roi.connection, 81, 151
         173
                                                   get.value.from.ijk, 74, 75, 82
dicom.tag.parser, 24, 25, 29, 32
                                                   get.value.from.mesh, 83
dicom.to.Rdcm.converter, 33, 100, 102,
                                                   get.value.from.xyz, 85
         104, 105, 145
                                                   get.volume.from.bin, 86, 87
dicom.viewer, 25, 35
                                                   get.volume.from.roi, 86, 87
                                                   get.xyz.from.index, 85, 88
display. 2D. histo, 36, 91
display.3D.contour, 38
                                                   grid.equal, 89
display.3D.mesh, 39
                                                   hist, 93, 95, 97
display.3D.sections, 40
```

INDEX

histo.2D, <i>37</i> , <i>61</i> , <i>66</i> , 89	rt.gamma.index, <i>130</i> , 131
histo.DVH, <i>61</i> , <i>66</i> , 92	rt.indices.from.bin, 133, 143
histo.from.bin, 66, 93, 93, 96, 98	rt.indices.from.roi, <i>135</i> , 136
histo.from.roi, 61, 66, 93, 94, 94, 98	
histo.vol, 66, 93, 94, 96, 97	save.T.MAT, 143
	save.to.Rdcm, 144
load.obj.data, 98, <i>100-102</i> , <i>104</i>	select.names, 87, 145
load.obj.from.dicom, 98, 99, 101, 102, 104	set.reference.obj, 146
load.obj.from.Rdcm, 59, 66, 98, 100, 100,	shade3d, <i>40</i>
102, 104, 105	sp.similarity.from.bin, 147, 149
load.patient.from.dicom, 52, 66, 77, 80,	sp.similarity.from.mesh, 148, 148
83, 85, 95, 102, 104, 109, 115, 125,	struct.clustering, 81, 150
126, 128, 143, 154, 164, 167	struct.from.bin, 151, <i>156</i>
load.patient.from.Rdcm, 9, 14, 39-42, 52,	struct.from.mesh, 153
56, 66, 73, 77, 80, 83, 85, 95, 102,	struct.in.new.ref, 154, <i>164</i>
103, 109, 115, 125, 126, 128, 134,	struct.merge, 155
137, 143, 150, 154, 164, 167	study.deployment, 156
load.Rdcm.raw.data, 104	
load.T.MAT, 9, 14, 39–42, 56, 59, 66, 73, 80,	toy.dicom.raw, 158
83, 85, 95, 102, 104, 106, 109, 115,	toy.load.patient, <i>52</i> , <i>66</i> , <i>77</i> , 159
125, 126, 128, 134, 137, 143, 150,	01 107
154, 164, 167	vcgClean, 107
	vcgSmooth, 107, 108
mesh.from.bin, 40, 61, 66, 107	vector.product, 160
mesh.in.new.ref, 109	vol.copy, 160
mesh.repair, 110	vol.create, 161
mesh.spheric.proj, 111	vol.from.bin, 162
	vol.gradient, 163
nesting.bin, 5, 112, 114, 116	vol.in.new.ref, <i>154</i> , 164
nesting.cube, 5, 7, 12, 13, 19, 113, 113, 116	vol.median, 165
nesting.roi, <i>5</i> , <i>113</i> , <i>114</i> , 114	vol.oversampling, 166, <i>170</i>
obj.create, 116	vol.regrid, 167
orientation.create, 117	vol.repair, 168
of Tentacton. create, 117	vol.subsampling, <i>166</i> , 169
pal.rainbow, 118, <i>122</i>	vol.sum, 171
pal.RVV, 119, <i>122</i>	xlsx.from.dcm, 25, 35, 172
plot, 44, 45, 47, 48, 120	xlsx.from.Rdcm, 25, 35, 172 xlsx.from.Rdcm, 25, 35, 173
plot.mesh, 58	X15X.17 OIII. RUCIII, 23, 33, 173
plot.struct, 58	
plot.volume, 58	
Rdcm.inventory, 123	
Rdcm.upgrade, 101, 104-106, 123, 124	
ref.add, 83, 109, 125, 126, 128, 144, 154,	
164, 167	
ref.cutplane.add, 125, 126, 127, 128, 144	
ref.remove, 125, 126, 127, 128	
ref.srctodest.add, <i>125</i> , <i>126</i> , 128	
rt.chi.index, 129, <i>133</i>	