# Package 'GENEAcore'

April 1, 2025
Title Pre-Processing of 'GENEActiv' Data
Version 1.0.0
<b>Date</b> 2025-03-24
Maintainer Jia Ying Chua <jiayingc@activinsights.com></jiayingc@activinsights.com>
Description  Analytics to read in and segment raw 'GENEActiv' accelerometer data into epochs and events.  For more details on the 'GENEActiv' device, see <a href="https://activinsights.com/resources/geneactiv-support-1-2/">https://activinsights.com/resources/geneactiv-support-1-2/</a> .
License GPL (>= 2)
Encoding UTF-8
RoxygenNote 7.3.2
Imports changepoint, signal, methods
<b>Suggests</b> knitr, rmarkdown, testthat (>= 3.0.0), GENEAread (>= 2.0.10), GENEAclassify
Config/testthat/edition 3
VignetteBuilder knitr
NeedsCompilation no
Author Joss Langford [aut], Ian Long [aut], Jia Ying Chua [aut, cre], Activinsights Ltd [cph]
Repository CRAN
<b>Date/Publication</b> 2025-04-01 16:30:10 UTC
Contents
aggregateEpochs aggregateEvents aggregatePeriods apply_AGSA apply_calibration

2 aggregateEpochs

	apply_degrees
	apply_ENMO
	apply_radians
	apply_updown
	binfile_summary
	calc_autocalparams
	createEventMapping
	create_MPI
	create_summary
	detect_nonmovement
	detect_transitions
	geneacore
	get_decimal_places
	get_UniqueBinFileIdentifier
	MPI_summary
	new_cut_times
	sample_binfile
	stepCounter
Index	22

 ${\tt aggregateEpochs}$ 

Aggregate Epochs

### Description

Aggregate Epochs

### Usage

```
aggregateEpochs(
  time_series,
  measure = "AGSA",
  time = "timestamp",
  sample_frequency,
  duration = NA,
  first_epoch_timestamp = NA,
  fun = mean
)
```

### Arguments

time\_series Data frame to be aggregated.

measure Name of the measure columns to be included.

time Name of the time column.

sample\_frequency

Measurement frequency of data.

aggregateEvents 3

```
duration Time duration to aggregate in each epoch.

first_epoch_timestamp

Time to start the first epoch, defaults to first record.

fun Function to apply on aggregation, defaults to mean.
```

#### **Details**

Wrapper function that calls aggregatePeriods for epochs (duration of fixed length).

#### Value

Data frame of aggregated epochs.

### **Examples**

```
timestamp <- c(</pre>
  1619424004, 1619424005, 1619424006, 1619424007,
  1619424008, 1619424009, 1619424010, 1619424011,
  1619424012, 1619424013, 1619424014, 1619424015
)
value <- c(</pre>
  0.729614366, 1.729115871, 0.804973546, 2.510181118,
  2.23764038, 0.613203747, 0.681953275, 0.089566943,
  0.021042388, 2.4780338, 2.437488989, 2.632635727
)
data <- data.frame(timestamp, value)</pre>
aggregated <- aggregateEpochs(data,</pre>
  duration = 5,
  measure = "value",
  sample_frequency = 1,
  first_epoch_timestamp = 1619424005,
  time = "timestamp"
)
```

aggregate Events

Aggregate Events

#### **Description**

Aggregate Events

```
aggregateEvents(
  time_series,
  measure = "AGSA",
  time = "timestamp",
  sample_frequency,
  events = NA,
```

4 aggregateEvents

```
start_time = "start",
end_time = "end",
fun = mean
)
```

#### **Arguments**

time\_series Data frame to be aggregated.

measure Name of the measure columns to be included.

time Name of the time column.

sample\_frequency

Measurement frequency of data.

events Data frame containing the start and end index of each event.

start\_time Name of the column in events containing the start index of the events.

Name of the column in events containing the end index of the events.

fun Function to apply on aggregation, defaults to mean.

#### **Details**

Wrapper function that calls aggregatePeriods for events (duration of variable length).

### Value

Data frame of aggregated events.

```
timestamp <- c(</pre>
  1619424004, 1619424005, 1619424006, 1619424007,
  1619424008, 1619424009, 1619424010, 1619424011,
  1619424012, 1619424013, 1619424014, 1619424015
)
value <- c(</pre>
  0.729614366, 1.729115871, 0.804973546, 2.510181118,
  2.23764038, 0.613203747, 0.681953275, 0.089566943,
  0.021042388, 2.4780338, 2.437488989, 2.632635727
)
data <- data.frame(timestamp, value)</pre>
event_start <- c(1, 5, 10)
event_end <- c(4, 9, 12)
aggregated_events <- aggregateEvents(data,</pre>
  events = data.frame(start = event_start, end = event_end),
  measure = "value",
  time = "timestamp".
  start_time = "start",
  end_time = "end",
  sample_frequency = 1,
  fun = sum
)
```

aggregatePeriods 5

aggregatePeriods Aggregate Periods

### **Description**

Generalised aggregation function generates distinct epochs or events outputs based on the initial parameters provided.

### Usage

```
aggregatePeriods(
  time_series,
  measure = "AGSA",
  time = "timestamp",
  sample_frequency,
  duration = NA,
  first_epoch_timestamp = NA,
  events = NA,
  start_time = "start",
  end_time = "end",
  fun = mean
)
```

### Arguments

time\_series Data frame to be aggregated.

measure Name of the measure columns to be included.

time Name of the time column.

sample\_frequency

Frequency of data.

duration Time duration to aggregate in each epoch.

 ${\tt first\_epoch\_timestamp}$ 

Time to start the first epoch, defaults to first record.

events Data frame containing the start and end index of each event.

start\_time Name of the column in events containing the start index of the events.

end\_time Name of the column in events containing the end index of the events.

fun Function to apply on aggregation, defaults to mean.

#### Value

Data frame of aggregated epochs or events.

6 apply\_calibration

apply\_AGSA

Apply Absolute Gravity-Subtracted Acceleration (AGSA)

### **Description**

Apply Absolute Gravity-Subtracted Acceleration (AGSA)

#### Usage

```
apply_AGSA(x)
```

### Arguments

Χ

Calibrated acceleration data frame.

#### Value

Measure column appended to end of calibrated data frame.

### **Examples**

```
x \leftarrow c(0.14268, 0.21757, -0.529, -0.36383)

y \leftarrow c(0.26385, 0.27295, 0.29220, 0.79510)

z \leftarrow c(0.27722, 0.20296, 0.35092, 0.27459)

calibrated \leftarrow data.frame(x, y, z)

calibrated \leftarrow apply_AGSA(calibrated)
```

apply\_calibration

Apply Calibration

### Description

Apply Calibration

#### Usage

```
apply_calibration(sensor_data, cal_params, measurement_device, use_temp = TRUE)
```

### **Arguments**

```
sensor_data Raw sensor-level data from a bin file in the form (x, y, z, light, button, temp).

cal_params Calibration parameters for acceleration and light from MPI.

measurement_device
```

Name of the measurement device used "GENEActiv 1.1" or "GENEActiv 1.2".

use\_temp Allows auto-calibration to be run with and without temperature compensation.

apply\_degrees 7

### **Details**

Function to apply calibration to sensor-level data from a bin file.

#### Value

Data frame of calibrated sensor data.

### Examples

```
cal_params <- list(</pre>
  scale = c(1.015, 1.017, 1.027),
  offset = c(0.00128, 0.0383, 0.0138),
  temperatureoffset = c(0, 0, 0),
  error = NA,
  lightdenominator = 48,
  lightnumerator = 911
)
rawdata <- data.frame(</pre>
  time = c(rep(1726650857, 5)),
  x = c(0.2421875, 0.24609375, 0.25390625, 0.24609375, 0.23828125),
  y = c(-0.04296875, -0.04687500, -0.03515625, -0.03125000, -0.04296875),
  z = c(-0.9453125, -0.9453125, -0.9531250, -0.9531250, -0.9609375),
  light = c(rep(22, 5)),
  button = c(rep(0, 5)),
  temp = c(rep(21.3, 5)),
  volts = c(rep(4.0896, 5))
)
calibrated <- apply_calibration(rawdata, cal_params, "GENEActiv 1.1")</pre>
```

apply\_degrees

Apply Rotation (degrees)

### Description

```
Apply Rotation (degrees)
```

### Usage

```
apply_degrees(x)
```

### **Arguments**

Х

Calibrated acceleration data frame.

#### Value

Measure column appended to end of calibrated data frame.

8 apply\_radians

#### **Examples**

```
x \leftarrow c(0.14268, 0.21757, -0.529, -0.36383) y \leftarrow c(0.26385, 0.27295, 0.29220, 0.79510) z \leftarrow c(0.27722, 0.20296, 0.35092, 0.27459) calibrated \leftarrow data.frame(x, y, z) calibrated \leftarrow apply_degrees(calibrated)
```

apply\_ENMO

Apply Euclidean Norm Minus One (ENMO)

## Description

Apply Euclidean Norm Minus One (ENMO)

### Usage

```
apply_ENMO(x)
```

### Arguments

x

Calibrated acceleration data frame.

### Value

Measure column appended to end of calibrated data frame.

### **Examples**

```
x \leftarrow c(0.14268, 0.21757, -0.529, -0.36383)

y \leftarrow c(0.26385, 0.27295, 0.29220, 0.79510)

z \leftarrow c(0.27722, 0.20296, 0.35092, 0.27459)

calibrated \leftarrow data.frame(x, y, z)

calibrated \leftarrow apply_ENMO(calibrated)
```

apply\_radians

Apply Rotation (radians)

### Description

Apply Rotation (radians)

```
apply_radians(x)
```

apply\_updown 9

### **Arguments**

x Calibrated acceleration data frame.

#### Value

Measure column appended to end of calibrated data frame.

### Examples

```
x \leftarrow c(0.14268, 0.21757, -0.529, -0.36383)

y \leftarrow c(0.26385, 0.27295, 0.29220, 0.79510)

z \leftarrow c(0.27722, 0.20296, 0.35092, 0.27459)

calibrated \leftarrow data.frame(x, y, z)

calibrated \leftarrow apply_radians(calibrated)
```

apply\_updown

Apply Elevation (updown)

### **Description**

Apply Elevation (updown)

#### Usage

```
apply_updown(x)
```

### Arguments

Χ

Calibrated acceleration data frame.

#### Value

Measure column appended to end of calibrated data frame.

```
x \leftarrow c(0.14268, 0.21757, -0.529, -0.36383)

y \leftarrow c(0.26385, 0.27295, 0.29220, 0.79510)

z \leftarrow c(0.27722, 0.20296, 0.35092, 0.27459)

calibrated \leftarrow data.frame(x, y, z)

calibrated \leftarrow apply_updown(calibrated)
```

10 calc\_autocalparams

binfile\_summary

Bin File Summary

### Description

Bin File Summary

### Usage

```
binfile_summary(input, recursive = TRUE)
```

#### **Arguments**

input

Bin file path.

recursive

TRUE applies the operation to all nested elements.

#### **Details**

Wrapper function that calls create\_summary for bin files only.

#### Value

Data frame of bin file or MPI summary.

calc\_autocalparams

Calculate Auto-calibration Parameters

### **Description**

Function to calculate auto-calibration parameters from known still points from a bin file that create a unitary sphere.

```
calc_autocalparams(
  binfile,
  binfile_path,
  output_folder,
  sphere_points,
  use_temp = TRUE,
  spherecrit = 0.3,
  maxiter = 500,
  tol = 1e-13
)
```

createEventMapping 11

#### Arguments

binfile Text lines read from an open connection to a bin file. Path to the bin file to be processed. binfile\_path output\_folder Path to the folder containing GENEAcore run outputs and Measurement Period Information (MPI) files. List of points that populate a unitary sphere and their associated temperature in sphere\_points the form (x,y,z,temp). Allows auto-calibration to be run with and without temperature compensation. use\_temp The minimum required acceleration value for each axis in both directions for spherecrit auto-calibration to be reliable. maxiter The maximum number of sphere fit iterations attempted during auto-calibration.

The limit of incremental sphere fit improvements before auto-calibration is con-

sidered complete.

#### Value

tol

List of auto-calibration parameters within the measurement period information (MPI).

#### **Examples**

```
binfile_path <- system.file("inst/extdata/10Hz_calibration_file.bin", package = "GENEAcore")
output_folder <- "."
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
MPI <- create_MPI(binfile, binfile_path, output_folder)
nonmovement_list <- detect_nonmovement(binfile, binfile_path, output_folder)
MPI <- calc_autocalparams(binfile, binfile_path, output_folder, nonmovement_list$sphere_points)</pre>
```

createEventMapping

Create Event Mapping

#### **Description**

Create Event Mapping

### Usage

```
createEventMapping(events, start_time, end_time, max_row_number)
```

### Arguments

events Data frame containing the start and end index of each event.

Start\_time Name of the column in events containing the start index of the events.

Name of the column in events containing the end index of the events.

max\_row\_number Number of rows in the source vector the events describe

12 create\_MPI

#### **Details**

Enumerate a vector to identify which event each measurement belongs to.

#### Value

List of mapped events.

### **Examples**

```
events <- data.frame(
   "start" = c(1, 5, 10, 15),
   "end" = c(4, 9, 14, 19)
)
time_series <- rnorm(25)
period_number <- createEventMapping(events, "start", "end", length(time_series))</pre>
```

create\_MPI

Create Measurement Period Information

### **Description**

Create Measurement Period Information

#### Usage

```
create_MPI(binfile, binfile_path, output_folder, out_rds = TRUE)
```

#### **Arguments**

binfile Text lines read from an open connection to a bin file.

binfile\_path Path to the bin file to be processed.

output\_folder Folder to write MPI file in.

out\_rds Allows RDS output to be saved during MPI creation.

#### **Details**

Function to create measurement period information (MPI) from a GENEActiv bin file

#### Value

List of measurement period information.

```
binfile_path <- system.file("inst/extdata/20Hz_file.bin", package = "GENEAcore")
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
MPI <- create_MPI(binfile)</pre>
```

create\_summary 13

create\_summary

Create Summary

#### **Description**

Create Summary

#### Usage

```
create_summary(input, path_type, recursive)
```

### **Arguments**

input Input type of either a bin file path, MPI path or an MPI object.

path\_type The file type within the folder to create summary for.
recursive TRUE applies the operation to all nested elements.

#### **Details**

Function to create a summary of key information of a bin file or MPI path.

#### Value

Data frame of bin file or MPI summary.

detect\_nonmovement

Detect Non-movement

#### **Description**

Detect Non-movement

```
detect_nonmovement(
  binfile,
  binfile_path,
  output_folder,
  still_seconds = 120,
  sd_threshold = 0.011,
  temp_seconds = 240,
  border_seconds = 300,
  long_still_seconds = 120 * 60,
  delta_temp_threshold = -0.7,
  posture_changes_max = 2,
  non_move_duration_max = 12 * 60 * 60
)
```

14 detect\_nonmovement

#### **Arguments**

binfile Text lines read from an open connection to a bin file. binfile\_path Path to the bin file to be processed. output\_folder Path to the folder containing GENEAcore run outputs and Measurement Period Information (MPI) files. still\_seconds The number of seconds included in the rolling standard deviation calculation for stillness to determine the shortest detection duration. sd\_threshold The threshold applied to the rolling standard deviation of combined acceleration to determine stillness. The number of seconds included in the rolling temperature difference calculatemp\_seconds tion or non-wear which also determines the shortest detection duration. border\_seconds The minimum number of seconds of a still event to be classed as a new bout.

border\_seconds The minimum number of seconds of a still event to be classed as a new bout. long\_still\_seconds

The minimum number of seconds of a still bout that is classed as non-wear.

delta\_temp\_threshold

The threshold applied to the rolling temperature difference to determine nonwear.

posture\_changes\_max

The maximum number of adjoining events that make up a single bout.

non\_move\_duration\_max

The maximum number of seconds of a still bout to be classed as non-movement. Still bouts with a duration longer than this number is automatically classed as non-wear.

#### **Details**

Function to detect non-movement bouts, non-wear events and points in a 1Hz downsampled bin file.

#### Value

List of sphere points, non-movement bouts and non-wear events.

```
binfile_path <- system.file("inst/extdata/20Hz_file.bin", package = "GENEAcore")
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
output_folder <- "."
MPI <- create_MPI(binfile, binfile_path, output_folder)
MPI <- detect_nonmovement(binfile, binfile_path, output_folder)</pre>
```

15 detect\_transitions

detect\_transitions

**Detect Transitions** 

### **Description**

**Detect Transitions** 

### Usage

```
detect_transitions(
  binfile,
 binfile_path,
 output_folder,
 minimum_event_duration = 3,
 x_{cpt_penalty} = 20,
 y_cpt_penalty = 30,
  z_{cpt_penalty} = 20,
 CutTime24Hr = "15:00"
)
```

### **Arguments**

binfile

binfile\_path Path to the bin file to be processed. output\_folder Path to the folder containing GENEAcore run outputs and Measurement Period Information (MPI) files. minimum\_event\_duration The minimum interval between changepoint transitions. The manual penalty value applied in the PELT changepoint algorithm for the x x\_cpt\_penalty axis, see cpt.var. The manual penalty value applied in the PELT changepoint algorithm for the y y\_cpt\_penalty axis, see cpt.var.

Text lines read from an open connection to a bin file.

The manual penalty value applied in the PELT changepoint algorithm for the z z\_cpt\_penalty

axis, see cpt.var.

CutTime24Hr Time in 24h to split the days up by.

### **Details**

Function to detect mean and variance changepoints in 1Hz acceleration data from a bin file.

#### Value

List of time, index and day number of each transition within the measurement period information.

16 geneacore

#### **Examples**

```
binfile_path <- system.file("inst/extdata/20Hz_file.bin", package = "GENEAcore")
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
output_folder <- "."
MPI <- create_MPI(binfile, binfile_path, output_folder)
MPI <- detect_transitions(binfile, binfile_path, output_folder)</pre>
```

geneacore

Main GENEAcore Function

### Description

Main GENEAcore Function

#### Usage

```
geneacore(
  data_folder = data_folder,
  CutTime24Hr = "15:00",
  output_epochs = TRUE,
  epoch_duration = 1,
  output_events = TRUE,
  output_steps = FALSE,
  output_csv = FALSE,
  timer = FALSE
)
```

### Arguments

data\_folder Folder that contains raw data bin files to process.

CutTime24Hr Time in 24h to split the days up by.

output\_epochs Create epoch outputs.

epoch\_duration Specify duration of fixed epochs.

output\_events Create event outputs.

output\_steps Include step counts and stepping rate outputs.

output\_csv Allows CSV output to be saved during epoch and event processing.

timer Print elapsed times of each process.

### Value

RDS and CSV of Measurement Period Information, Epoch measures and Event measures.

get\_decimal\_places 17

get\_decimal\_places

Get Decimal Places

### **Description**

Get Decimal Places

#### Usage

```
get_decimal_places(column)
```

### **Arguments**

column

Aggregated data frame column.

#### **Details**

Function to determine the number of decimal places based on column name.

#### Value

Decimal place integer.

### **Examples**

```
epochs_df <- data.frame(
   "x.mean" = c(0.1111, 0.1222, 0.1333, 0.1444),
   "y.mean" = c(0.2111, 0.2222, 0.2333, 0.2444),
   "light.mean" = c(1.25, 1.73, 1.99, 2.02)
)
dp <- get_decimal_places(epochs_df[1])</pre>
```

```
get_UniqueBinFileIdentifier
```

Generate Unique Bin File Identifier

### **Description**

Generate Unique Bin File Identifier

## Usage

```
get_UniqueBinFileIdentifier(binfile)
```

### **Arguments**

binfile

Text lines read from an open connection to a bin file.

MPI\_summary

### **Details**

Function to create a UniqueBinFileIdentifier from a GENEActiv bin file.

#### Value

Single string identifier.

### **Examples**

```
binfile_path <- system.file("inst/extdata/20Hz_file.bin", package = "GENEAcore")
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
UniqueBinFileIdentifier <- get_UniqueBinFileIdentifier(binfile)</pre>
```

MPI\_summary

**MPI Summary** 

### Description

**MPI Summary** 

### Usage

```
MPI_summary(input, recursive = TRUE)
```

### Arguments

input MPI path.

recursive TRUE applies the operation to all nested elements.

### **Details**

Wrapper function that calls create\_summary for MPI only.

#### Value

Data frame of MPI summary.

new\_cut\_times 19

new\_cut\_times

New Cut Times

### **Description**

New Cut Times

#### Usage

```
new_cut_times(df)
```

### Arguments

df

Cut Times data frame.

#### **Details**

Add the timestamps, indexes and day numbers of the cut times and their ends.

#### Value

Data frame with added cut times.

### **Examples**

```
CutTimes_df <- data.frame(
   time = c(1731421000, 1731421100, 1731421362, 1731421480, 1731421525),
   index = c(56, 1, 230, 1, 400), day = c(1, 2, 2, 3, 3)
)
CutTimes_df <- new_cut_times(CutTimes_df)</pre>
```

sample\_binfile

Sample Bin File

### Description

Sample Bin File

```
sample_binfile(
  binfile,
  binfile_path,
  output_folder,
  start_time = NULL,
  end_time = NULL,
  downsample = TRUE,
  output_csv = FALSE
)
```

20 stepCounter

### **Arguments**

binfile	Text lines read from an open connection to a bin file.
binfile_path	Path to the bin file to be processed.
output_folder	Path to the folder containing GENEAcore run outputs and Measurement Period Information (MPI) files.
start_time	Time stamp to start the read from, default start of file.
end_time	Time stamp to end the read from, default end of file.
downsample	Logical to determine whether to downsample the file, default TRUE.
output_csv	Allow outputs of bin file sampling to be saved as CSV.

### **Details**

Function to read in a GENEActiv bin file with option to downsample to 1Hz.

#### Value

List of 1Hz downsampled data or raw sample data.

### **Examples**

```
binfile_path <- system.file("inst/extdata/20Hz_file.bin", package = "GENEAcore")
output_folder <- "."
con <- file(binfile_path, "r")
binfile <- readLines(con, skipNul = TRUE)
close(con)
measurements <- sample_binfile(binfile, binfile_path, output_folder)</pre>
```

stepCounter

Step Counter

### **Description**

Function to calculate the number and variance of the steps in the data.

```
stepCounter(
  StepData,
  samplefreq = 100,
  filterorder = 2,
  boundaries = c(0.5, 5),
  Rp = 3,
  hysteresis = 0.05,
  fun = c("GENEAcount", "mean", "sd")
)
```

stepCounter 21

### **Arguments**

StepData	The data to use for calculating the steps. This should either an AccData object or a vector.
samplefreq	The sampling frequency of the data, in hertz, when calculating the step number (default 100).
filterorder	single integer, order of the Chebyshev bandpass filter, passed to argument n of cheby1.
boundaries	length 2 numeric vector specifying lower and upper bounds of Chebychev filter (default c(0.5, 5) Hz), passed to argument W of butter or cheby1.
Rp	the decibel level that the cheby filter takes, see cheby1.
hysteresis	The hysteresis applied after zero crossing. (default 100mg)
fun	character vector naming functions by which to summarize steps. "count" is an internally implemented summarizing function that returns step count.

### Value

Returns a vector with length fun.

```
d1 <- sin(seq(0.1, 100, 0.1))/2 + rnorm(1000)/10 + 1
Steps4 = stepCounter(d1)
length(Steps4)
mean(Steps4)
sd(Steps4)
plot(Steps4)</pre>
```

## **Index**

```
{\it aggregateEpochs}, {\it 2}
aggregateEvents, 3
aggregatePeriods, 5
apply_AGSA, 6
apply_calibration, 6
apply_degrees, 7
apply_ENMO, 8
{\tt apply\_radians}, {\tt 8}
apply_updown, 9
\verb|binfile_summary|, 10
butter, 21
calc_autocalparams, 10
cheby1, 21
cpt.var, 15
create_MPI, 12
create_summary, 13
{\tt createEventMapping}, 11
detect_nonmovement, 13
detect_transitions, 15
geneacore, 16
get_decimal_places, 17
get_UniqueBinFileIdentifier, 17
MPI_summary, 18
new_cut_times, 19
sample_binfile, 19
stepCounter, 20
```