Package 'emov'

October 13, 2022

Version 0.1.1

Date 2016-04-04
Title Eye Movement Analysis Package for Fixation and Saccade Detection
Author Simon Schwab <schw4b@gmail.com></schw4b@gmail.com>
Maintainer Simon Schwab <schw4b@gmail.com></schw4b@gmail.com>
Depends R (>= 1.8.0)
Description Fixation and saccade detection in eye movement recordings. This package implements a dispersion-based algorithm (I-DT) proposed by Salvucci & Goldberg (2000) which detects fixation duration and position.
License GPL-3
<pre>URL https://github.com/schw4b/emov</pre>
BugReports https://github.com/schw4b/emov/issues
NeedsCompilation no
Repository CRAN
Date/Publication 2016-04-04 18:39:38
R topics documented:
emov.angdia
emov.cart2sphere
emov.filter
emov.idt
emov.read_iviewsamples
fivesec
Index 5

emov.cart2sphere

emov.angdia

Angular size of stimulus.

Description

Angular size of stimulus.

Usage

```
emov.angdia(stimsize, distance)
```

Arguments

stimsize

Size of the stimulus.

distance

Viewing distance from stimulus.

Value

Angular size in degrees.

emov.cart2sphere

Convert Cartesian to Spherical coordinates.

Description

Convert Cartesian to Spherical coordinates.

Usage

```
emov.cart2sphere(x, y, z)
```

Arguments

Χ

х.

y z y. z.

Value

Two angles (radians) and radius

Examples

```
data = emov.cart2sphere(3, 4, 5)
```

emov.filter 3

emov	fi	1	t.	er	

Velocity threshold filter.

Description

Velocity threshold filter.

Usage

```
emov.filter(x, y, threshold)
```

Arguments

x Eye position.y Eye position.threshold Velocity threshold.

Value

Filtered data.

emov.idt

I-DT algorithm.

Description

I-DT algorithm.

Usage

```
emov.idt(t, x, y, dispersion, duration)
```

Arguments

t Vector of timepoints.x horizontal eye positions.y vertical eye positions.

dispersion Maximal dispersion allowed (in units of x and y).

duration Minimal fixation duration allowed (in number of samples)

Value

Fixations: position, start, end.

4 fivesec

References

Salvucci, D. D., & Goldberg, J. H. (2000). Identifying fixations and saccades in eye-tracking protocols. In Proceedings of the 2000 symposium on eye tracking research & applications (pp. 71-78). New York: ACM.

```
emov.read_iviewsamples
```

Read SMI iview sample file.

Description

Read SMI iview sample file.

Usage

```
emov.read_iviewsamples(file, nr_of_headerlines)
```

Arguments

```
file Filename. nr_of_headerlines
```

No. of header lines in datafile.

Value

data file.

fivesec

Eye movement data

Description

Five seconds of eye movement data recorded with an SMI eye tracker 200 Hz

Usage

fivesec

Format

A data.frame that contains time, x and y eye positions.

Source

Simon Schwab

Index

```
* datasets
    fivesec, 4

emov.angdia, 2
emov.cart2sphere, 2
emov.filter, 3
emov.idt, 3
emov.read_iviewsamples, 4

fivesec, 4
```