# Package 'dChipIO'

October 13, 2022

Version 0.1.5

Title Methods for Reading dChip Files
<b>Depends</b> R (>= $2.6.0$ )
<b>Date</b> 2016-01-12
Author Henrik Bengtsson [aut, cre, cph]
Maintainer Henrik Bengtsson <henrikb@braju.com></henrikb@braju.com>
<b>Description</b> Functions for reading DCP and CDF.bin files generated by the dChip software.
License LGPL (>= 2.1)
<pre>URL https://github.com/HenrikBengtsson/dChipIO</pre>
BugReports https://github.com/HenrikBengtsson/dChipIO/issues
LazyLoad TRUE
biocViews Infrastructure, DataImport
NeedsCompilation no
Repository CRAN
<b>Date/Publication</b> 2016-01-13 22:36:24
R topics documented:
dChipIO-package
readCdfBin
readDcp
readDcpRectangle

2 readCdfBin

dChipIO-package

Package dChipIO

# **Description**

Functions for reading DCP and CDF.bin files generated by the dChip software.

The example data used in this package orginates from the Affymetrix Fusion SDK library [3].

# To get started

To get started, see:

- 1. readCdfBin() reads a dChip CDF.bin file.
- 2. readDcp() reads a dChip DCP data file.

#### License

The releases of this package is licensed under LGPL version 2.1 or newer.

#### Author(s)

Henrik Bengtsson.

#### References

- [1] The dChip software, http://www.dchip.org/
- [2] Thread 'DCP File Format', 'dChip Software', Google Groups, December 2008. https://groups.google.com/forum/#!topic/dchip-software/Q7mTJPPpZ5U
- [3] Affymetrix Inc, Fusion Software Developers Kit (SDK), 2008. http://www.affymetrix.com/estore/partners\_programs/programs/developer/fusion/index.affx?terms=no

readCdfBin

Reads a dChip CDF.bin file

# **Description**

Reads a dChip CDF.bin file.

Please note that this method is incomplete as it currently doesn't read all fields. It is only made available so that someelse can continue the development.

# Usage

```
readCdfBin(con, units=NULL, ...)
```

readCdfBinHeader 3

# Arguments

con A connection or a character filename.
units An integer vector specifying the units to be read. If NULL, all units are read.
... Not used.

#### Value

Returns a list structure containing the file header and the unit data.

#### Author(s)

Henrik Bengtsson

#### See Also

To read only the CDF.bin file header, see readCdfBinHeader().

# **Examples**

```
path <- system.file("exData", package="dChipIO")</pre>
chipType <- "Test3"</pre>
filename <- sprintf("%s.CDF.bin", chipType)</pre>
pathname <- file.path(path, filename)</pre>
hdr <- readCdfBinHeader(pathname)</pre>
print(hdr)
data <- readCdfBin(pathname)</pre>
str(data)
# Read a subset of the units
units <- c(10:11, 15:20, 150:105, 2,2,2)
dataT <- readCdfBin(pathname, units=units)</pre>
str(dataT)
# Assert correctness
for (ff in c("unitNames", "numProbes", "CellPos")) {
  stopifnot(length(dataT[[ff]]) == length(units))
  stopifnot(identical(dataT[[ff]], data[[ff]][units]))
}
```

read Cdf Bin Header

Reads the file header of a dChip CDF.bin file

# **Description**

Reads the file header of a dChip CDF.bin file.

readDcp

#### Usage

```
readCdfBinHeader(con, ...)
```

# **Arguments**

con A connection or a character filename.

... Not used.

#### Value

Returns a list structure containing the file header.

#### Author(s)

Henrik Bengtsson

#### See Also

To read the CDF.bin file data, see readCdfBin().

readDcp

Reads a dChip DCP file

#### **Description**

Reads a dChip DCP file.

# Usage

```
readDcp(con, fields=c("rawIntensities", "normalizedIntensities", "calls", "thetas",
   "thetaStds", "excludes"), cells=NULL, units=NULL, .nbrOfUnits=NULL, ...)
```

# **Arguments**

con	A connection or a ch	aracter filename
con	A connection of a ch	aracter mename.

fields A character vector specifying the fields to be read.

cells An integer vector specifying the indices of the cell data to be read.

An integer vector specifying the indices of the unit data to be read.

.nbrOfUnits A integer specifying the number of units available in the file. If NULL, this is

inferred from the file size and the file header. The dChip software itself instead uses the corrsponding value in the CDF.bin file, but that file is specified by the

user leaving room for errors.

... Not used.

readDcpHeader 5

# Value

Returns a list structure containing the file header and the requested data fields.

#### Author(s)

Henrik Bengtsson

#### See Also

To read only the DCP file header, see readDcpHeader().

# **Examples**

```
path <- system.file("exData", package="dChipIO")</pre>
filename <- "Test3-1-121502.dcp"
pathname <- file.path(path, filename)</pre>
hdr <- readDcpHeader(pathname)</pre>
print(hdr)
data <- readDcp(pathname)</pre>
str(data)
# Read a subset of the units
units <- c(10:11, 15:20, 150:105, 2,2,2)
dataT <- readDcp(pathname, units=units)</pre>
str(dataT)
# Assert correctness
stopifnot(length(dataT[[ff]]) == length(units))
 stopifnot(identical(dataT[[ff]], data[[ff]][units]))
}
```

readDcpHeader

Reads the file header of a dChip DCP file

# **Description**

Reads the file header of a dChip DCP file.

# Usage

```
readDcpHeader(con, ...)
```

#### **Arguments**

```
con A connection or a character filename.
... Not used.
```

6 readDcpRectangle

# Value

Returns a list structure containing the file header.

# Author(s)

Henrik Bengtsson

# See Also

To read also the DCP file data, see readDcp().

 ${\it readDcpRectangle} \qquad {\it Reads~a~spatial~subset~of~probe-level~data~from~a~dChip~DCP~file}$ 

# Description

Reads a spatial subset of probe-level data from a dChip DCP file.

# Usage

```
readDcpRectangle(filename, fields=c("rawIntensities", "normalizedIntensities"),
    xrange=c(0, Inf), yrange=c(0, Inf), ..., asMatrix=TRUE)
```

# **Arguments**

filename	The pathname of the DCP file.
fields	The cell fields to be read.
xrange	A numeric vector of length two giving the left and right coordinates of the cells to be returned.
yrange	A numeric vector of length two giving the top and bottom coordinates of the cells to be returned.
asMatrix	If TRUE, the CEL data fields are returned as matrices with element (1,1) corresponding to cell (xrange[1],yrange[1]).
	Additional arguments passed to readDcp().

#### Value

A named list CEL structure similar to what readDcp(). In addition, if asMatrix is TRUE, the CEL data fields are returned as matrices, otherwise not.

# Author(s)

Henrik Bengtsson

readDcpRectangle 7

# See Also

The readDcp() method is used internally. This method was inspired by readCelRectangle() of the **affxparser** package.

# **Examples**

```
path <- system.file("exData", package="dChipI0")

filename <- "Test3-1-121502.dcp"
pathname <- file.path(path, filename)

data <- readDcpRectangle(pathname)

layout(matrix(1:4, nrow=2, byrow=TRUE))
image(data$rawIntensities, main="Raw probe signals")
image(data$normalizedIntensities, main="Normalized probe signals")</pre>
```

# **Index**

```
* IO
    {\sf readCdfBin}, 2
    {\sf readCdfBinHeader}, 3
    readDcp, 4
    readDcpHeader, 5
    readDcpRectangle, 6
* file
    readCdfBin, 2
    readCdfBinHeader, 3
    readDcp, 4
    readDcpHeader, 5
    readDcpRectangle, 6
* package
    dChipIO-package, 2
character, 3-5
connection, 3-5
dChipIO (dChipIO-package), 2
dChipIO-package, 2
integer, 3, 4
list, 3–6
NULL, 3, 4
numeric, 6
readCdfBin, 2, 2, 4
readCdfBinHeader, 3, 3
readDcp, 2, 4, 6, 7
readDcpHeader, 5, 5
readDcpRectangle, 6
TRUE, 6
vector, 3, 4, 6
```