Package 'vfprogression'

October 12, 2022

Type Package
Title Visual Field (VF) Progression Analysis and Plotting Methods
Version 0.7.1
Date 2019-05-21
Author Tobias Elze, Dian Li (documentation), Eun Young Choi (QC)
Maintainer Dian Li lidian@zju.edu.cn>
Description Realization of published methods to analyze visual field (VF) progression. Introduction to the plotting methods (designed by author TE) for VF output visualization. A sample dataset for two eyes, each with 10 follow-ups is included. The VF analysis methods could be found in Musch et al. (1999) <doi:10.1016 s0161-6420(99)90147-1="">, Nouri-Mahdavi et at. (2012) <doi:10.1167 iovs.11-9021="">, Schell et at. (2014) <doi:10.1016 j.ophtha.2014.02.021="">, Aptel et al. (2015) <doi:10.1111 aos.12788="">.</doi:10.1111></doi:10.1016></doi:10.1167></doi:10.1016>
License GPL (>= 2)
Encoding UTF-8
LazyData true
RoxygenNote 6.1.1
Depends R (>= 2.10)
Imports stats, grDevices, graphics
NeedsCompilation no
Repository CRAN
Date/Publication 2019-05-24 15:00:02 UTC
R topics documented:
plotComponentMatrix plotfield.normalized plotTdProbabilities plotTDvalues progression progression.cigts

2 plotComponentMatrix

progression.schel progression.vfi															
vf.cigts															
vf.plr.nouri.2012															
vf.schell2014 .															
vf.vfi															
vfseries	 														

Index **15**

plotComponentMatrix

General plotting function for multiple 24-2 or 30-2 visual field measurements together:

Description

plotComponentMatrix plots the following 24-2 or 30-2 visual field measurement: sensitivity, TD, TD prob, PD, and PD prob:

Usage

```
plotComponentMatrix(componentmatrix, ncomp = ncol(componentmatrix),
  plot.ncols = 5, plot.nrows = NULL,
 plot.annot.topleft.function = toString,
 plot.annot.bottomleft.function = function(i) NULL,
  globaltitle = sprintf("k = %i", ncol(componentmatrix)),
  globalannotright = NULL,
  zmin = -ceiling(max(abs(c(min(componentmatrix),
 max(componentmatrix))))), zmax = -zmin,
  color.pal = colorRampPalette(c("red", "white", "blue"), space =
  "Lab")(256), td.probabilities = FALSE,
  show.colorbar = !td.probabilities, titleheight = 0.2, ...)
```

Arguments

ncomp

```
componentmatrix
```

a matrix or data frame, column represents different eyes and rows are the VF measurements of the same type (sensitivity, TD, TD prob, PD, or PD prob). a numeric variable defines the number of components to be plotted (default: all). plot.ncols a numeric variable defines the number of columns to be plotted (default: 5). plot.nrows a numeric variable defines the number of rows to be plotted (default: NULL (automatically calculated)).

plot.annot.topleft.function

a function(i) that is given to any subplot i to create its top left annotation.

plotComponentMatrix 3

```
plot.annot.bottomleft.function
                  a function(i) that is given to any subplot i to create its bottom left annotation
                  (default: returns NULL).
globaltitle
                  a string for global title (default: k = ncomp; set to NULL to suppress global
globalannotright
                  a string annotation to the right of the global title (default: NULL).
                  minimum value of the color scale (default: auto defined).
zmin
                  maximum value of the color scale (default: auto defined).
zmax
color.pal
                  an object that defines color scale theme (default: colorRampPalette(c("red",
                   "white", "blue"), space = Lab")(256)).
td.probabilities
                  a logic variable indicates whether to plot TD probability symbols instead of TD
                  colors (default: FALSE).
show.colorbar
                  a logic variable indicates whether to show a global colorbar (default: !td.probabilities).
titleheight
                  a numeric variable defines the height of the title relative to height of row one.
                  other variables to be added.
```

Value

heatmap for sensitivity, TD and PD input. Value plot for TD prob and PD prob input.

```
data(vfseries)
componentmatrix = t(vfseries[1:10, grepl('^s[0-9]+', colnames(vfseries))])
globaltitle = paste("Sensitivities, k = ", ncol(componentmatrix), sep = '')
plotComponentMatrix(componentmatrix, globaltitle = globaltitle)
componentmatrix = t(vfseries[1:10, grepl('^td[0-9]+', colnames(vfseries))])
globaltitle = paste("TDs, k = ", ncol(componentmatrix), sep = '')
plotComponentMatrix(componentmatrix, globaltitle = globaltitle)
componentmatrix = t(vfseries[1:10, grepl('^pd[0-9]+', colnames(vfseries))])
globaltitle = paste("PDs, k = ", ncol(componentmatrix), sep = '')
plotComponentMatrix(componentmatrix, globaltitle = globaltitle)
componentmatrix = t(vfseries[1:10, grepl('^tdp[0-9]+', colnames(vfseries))])
globaltitle = paste("TD Probs, k = ", ncol(componentmatrix), sep = '')
plotComponentMatrix(componentmatrix, globaltitle = globaltitle, td.probabilities = TRUE)
componentmatrix = t(vfseries[1:10, grepl('^pdp[0-9]+', colnames(vfseries))])
globaltitle = paste("PD Probs, k = ", ncol(componentmatrix), sep = '')
plotComponentMatrix(componentmatrix, globaltitle = globaltitle, td.probabilities = TRUE)
```

4 plotfield.normalized

plotfield.normalized Single plotting function for one 24-2 or 30-2 visual field measurement:

Description

plotfield.normalized plots the following 24-2 or 30-2 visual field measurement: sensitivity, TD, TD prob, PD, and PD prob:

Usage

```
plotfield.normalized(eigenfields, component = 1,
  zmin = -max(abs(c(min(eigenfields), max(eigenfields)))),
  zmax = max(abs(c(min(eigenfields), max(eigenfields)))),
  color.pal = colorRampPalette(c("red", "white", "blue"), space =
  "Lab")(256), show.colorbar = TRUE, topleftannotation = NULL,
  bottomleftannotation = NULL, labelcex = 2, ...)
```

Arguments

component

color.pal

eigenfields a vector contains Sensitivity/TD/PD measurement. For 24-2 VF eigenfields should have 52 or 54 elements. For 30-2 VF, eigenfields should have 74 or 76 elements.

Number of components to be plotted (default: 1). zmin minimum value of the color scale (default: auto defined).

maximum value of the color scale (default: auto defined). zmax

an object that defines color scale theme (default: colorRampPalette(c("red",

"white", "blue"), space = Lab")(256)).

show.colorbar a logic value to show colorbar (default: TRUE).

topleftannotation

a string annotation shown on the top left side of the plot (default: NULL).

bottomleftannotation

a string annotation shown on the bottom left side of the plot (default: NULL).

labelcex a numeric variable for label size (default: 2).

other variables to be added. . . .

Value

heatmap for sensitivity, TD and PD input

```
data(vfseries)
eigenfields = t(vfseries[1, grepl('^s[0-9]+', colnames(vfseries))])
plotfield.normalized(eigenfields)
title(main = "Sensitivity", line = 3)
```

plotTdProbabilities 5

```
eigenfields = t(vfseries[1, grepl('^td[0-9]+', colnames(vfseries))])
plotfield.normalized(eigenfields)
title(main = "Total Deviation", line = 3)
eigenfields = t(vfseries[1, grepl('^pd[0-9]+', colnames(vfseries))])
plotfield.normalized(eigenfields)
title(main = "Pattern Deviation", line = 3)
```

plotTdProbabilities

Value plotting function for 24-2 or 30-2 visual field measurement:

Description

plotTdProbabilities plots the following 24-2 or 30-2 visual field measurement: TD probs, and PD probs:

Usage

```
plotTdProbabilities(tdprob, cex = 2, rectangle.color = "black", rectangle.width = 0.16, margins = c(2, 1, 2, 2) + 0.1, ...)
```

Arguments

Value

value plot for TD prob and PD prob input.

```
data(vfseries)
tdprob = t(vfseries[1, grepl('^tdp[0-9]+', colnames(vfseries))])
plotTdProbabilities(tdprob)
title(main = "Total Deviation Probability", line = 3)
tdprob = t(vfseries[1, grepl('^pdp[0-9]+', colnames(vfseries))])
plotTdProbabilities(tdprob)
title(main = "Pattern Deviation Probability", line = 3)
```

6 plotTDvalues

$n l \cap t$	TDvalues

Value plotting function for 24-2 or 30-2 visual field measurement:

Description

plotTDvalues plots the following 24-2 or 30-2 visual field measurement: sensitivity, TD, and PD:

Usage

```
plotTDvalues(tds, cex.tds = 1, textcolor = function(x) "black",
    show.lines = T, ...)
```

Arguments

tds	a vector contains sensitivity/TD/PD measurement. For 24-2 VF tds should have 52 or 54 elements. For 30-2 VF, tds should have 74 or 76 elements.
cex.tds	a numeric variable for label size (default: 1).
textcolor	a function defines the label color.
show.lines	a logical variable indicates whether to show the horizontal and vertical lines.
	other variables to be added.

Value

value plot for sensitivity, TD and PD input.

```
data(vfseries)
tds = t(vfseries[1, grepl('^s[0-9]+', colnames(vfseries))])
plotTDvalues(tds)
title(main = "Sensitivity", line = 3)
tds = t(vfseries[1, grepl('^td[0-9]+', colnames(vfseries))])
plotTDvalues(tds)
title(main = "Total Dviation", line = 3)
tds = t(vfseries[1, grepl('^pd[0-9]+', colnames(vfseries))])
plotTDvalues(tds)
title(main = "Pattern Dviation", line = 3)
```

progression 7

progression	general progression function	

Description

progression returns the progression criterion with four methods. plr.nouri.2012, vfi, schell2014, cigts

Usage

```
progression(vfseries, method = c("plr.nouri.2012", "vfi", "schell2014",
    "cigts"))
```

Arguments

vfseries is a data frame. MUST contain the following columns: yearsfollowed', and

'eyeid'. Rows represent the single measurements. Other requirements, such as number of minimum measurements (rows), and necessary VF measurements

could be found in each progression method's documentation

method selected from one or more from: plr.nouri.2012, vfi, schell2014, cigts. Default

it ...

Value

```
"stable", "worsening", or "improving" of measurements in measmatrix
```

See Also

```
https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4495761/
```

```
data(vfseries)
progression(vfseries)
progression(vfseries[vfseries$eyeid == 1,])
progression(vfseries[vfseries$eyeid == 2,])
progression(vfseries, method=c("cigts"))
```

progression.cigts

CIGTS VF progression

Description

progression.cigts returns the progression of visual field test based on 52 or 54 total deviation probabilities (tdp). CIGTS VF progression (Musch et al., 1999).

Usage

```
progression.cigts(measmatrix)
```

Arguments

measmatrix

is a data frame. MUST contain the following columns: 52/54 TD probs (column names MUST be 'tdp1' ~ 'tdp52' or 'tdp1' ~ 'tdp54'), 'yearsfollowed', and 'eyeid'. Rows represent the single measurements. The minimum measurements (rows) is 5.

Value

"stable", "worsening", or "improving" of measurements in measmatrix. Note: If a VF series is temporarily improving and temporarily worsening, it is assumed to be "stable" overall

References

```
http://www.aaojournal.org/article/S0161-6420(99)90147-1/abstract
```

Examples

```
data(vf.cigts)
colnames(vf.cigts)
progression.cigts(vf.cigts)
progression.cigts(vf.cigts[vf.cigts$eyeid == 1,])
progression.cigts(vf.cigts[vf.cigts$eyeid == 2,])
```

```
progression.plr.nouri.2012
```

Nouri-Mahdavi 2012 VF progression

Description

progression.plr.nouri.2012 returns the progression criterion, using Pointwise Linear Regression (PLR) progression detection method according to Nouri-Mahdavi et al. (2012).

progression.schell2014

Usage

```
progression.plr.nouri.2012(measmatrix)
```

Arguments

measmatrix

is a data frame. MUST contain the following columns: 52/54 TD (column names MUST be 'td1' ~ 'td52' or 'td1' ~ 'td54'), 'yearsfollowed', and 'eyeid'. Rows represent the single measurements. The minimum measurements (rows) is 3.

Value

"stable", "worsening", or "improving" of measurements in measmatrix

See Also

```
https://www.ncbi.nlm.nih.gov/pubmed/22427560/
```

Examples

```
data(vf.plr.nouri.2012)
colnames(vf.plr.nouri.2012)
progression.plr.nouri.2012(vf.plr.nouri.2012)
progression.plr.nouri.2012(vf.plr.nouri.2012[vf.plr.nouri.2012$eyeid == 1,])
progression.plr.nouri.2012(vf.plr.nouri.2012[vf.plr.nouri.2012$eyeid == 2,])
```

progression.schell2014

Schell 2014 VF progression

Description

progression.schell2014 returns the progression criterion after Schell et al. 2014, which is essentially like CIGTS but with MD, and only one follow-up is enough to confirm progression.

Usage

```
progression.schell2014(measmatrix)
```

Arguments

measmatrix

is a data frame. MUST contain the following columns: 'md' (mean deviation) and 'eyeid'. Rows represent the single measurements. The minimum measurements (rows) is 4.

Value

"stable", "worsening", or "improving" of measurements in measmatrix. Note: If a VF series is temporarily improving and temporarily worsening, it is assumed to be "stable" overall

10 progression.vfi

See Also

```
https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4495761/
```

Examples

```
data(vf.schell2014)
colnames(vf.schell2014)
progression.schell2014(vf.schell2014)
progression.schell2014(vf.schell2014[vf.schell2014$eyeid == 1,])
progression.schell2014(vf.schell2014[vf.schell2014$eyeid == 2,])
```

progression.vfi

progression according to VFI (significant slope, p <= 0.05)

Description

progression.vfi returns the progression criterion used in Aptel et al. (2015).

Usage

```
progression.vfi(measmatrix)
```

Arguments

measmatrix

is a data frame. MUST contain the following columns: 'vfi' (visual field index), 'yearsfollowed', and 'eyeid'. Rows represent the single measurements. The minimum measurements (rows) is 3.

Value

```
"stable", "worsening", or "improving" of measurements in timepoints
```

See Also

```
https://www.ncbi.nlm.nih.gov/pubmed/26095771/
```

```
data(vf.vfi)
colnames(vf.vfi)
progression.vfi(vf.vfi)
progression.vfi(vf.vfi[vf.vfi$eyeid == 1,])
progression.vfi(vf.vfi[vf.vfi$eyeid == 2,])
```

vf.cigts 11

vf.cigts

Combined Visual Field Series for General Progression Method

Description

Data

Usage

```
data(vf.cigts)
```

Format

A data frame sample for CIGTS progression method, which includes visual field related measurement for two eyes, each with 10 follow-ups. Rows represent the single measurements.

Source

```
eyeid eyeid, labeled as 1,2... for different eyes.
```

yearsfollowed follow-up years. The minimum measurements /rows for one eye is 5.

tdp1-tdp54 52 total deviation probability, or 'tdp' measurements. The minimum measurements, or rows for one eye is 5. ...

Examples

```
data(vf.cigts)
colnames(vf.cigts)
progression.cigts(vf.cigts)
progression.cigts(vf.cigts[vf.cigts$eyeid == 1,])
progression.cigts(vf.cigts[vf.cigts$eyeid == 2,])
```

vf.plr.nouri.2012

Combined Visual Field Series for General Progression Method

Description

Data

Usage

```
data(vf.plr.nouri.2012)
```

Format

A data frame sample for Pointwise Linear Regression (PLR) method according to Nouri-Mahdavi 2012 progression, which includes visual field related measurement for two eyes, each with 10 follow-ups. Rows represent the single measurements.

12 vf.schell2014

Source

```
eyeid eyeid, labeled as 1,2... for different eyes
yearsfollowed follow-up years. The minimum measurements, or rows, for one eye is 3
td1-td54 52 total deviation, or 'td' measurements. The minimum measurements, or rows, for one eye is 3 ...
```

Examples

```
data(vf.plr.nouri.2012)
colnames(vf.plr.nouri.2012)
progression.plr.nouri.2012(vf.plr.nouri.2012)
progression.plr.nouri.2012(vf.plr.nouri.2012[vf.plr.nouri.2012$eyeid == 1,])
progression.plr.nouri.2012(vf.plr.nouri.2012[vf.plr.nouri.2012$eyeid == 2,])
```

vf.schell2014

Combined Visual Field Series for General Progression Method

Description

Data

Usage

```
data(vf.schell2014)
```

Format

A data frame sample for progression method by Schell et al. 2014, which includes visual field related measurement for two eyes, each with 10 follow-ups. Rows represent the single measurements.

Source

```
eyeid eyeid, labeled as 1,2... for different eyes.

md mean deviation measurements. The minimum measurements, or rows, for one eye is 4. ...
```

```
data(vf.schell2014)
colnames(vf.schell2014)
progression.schell2014(vf.schell2014)
progression.schell2014(vf.schell2014[vf.schell2014$eyeid == 1,])
progression.schell2014(vf.schell2014[vf.schell2014$eyeid == 2,])
```

vf.vfi

vf.vfi

Combined Visual Field Series for General Progression Method

Description

Data

Usage

```
data(vf.vfi)
```

Format

A data frame for CIGTS progression example, which includes visual field related measurement for two eyes each with 10 follow-ups.

Source

```
eyeid eyeid, labeled as 1,2... for different eye groups.

yearsfollowed follow-up years. The minimum measurements, or rows, for one eye is 3.

vfi visual field index. The minimum measurements, or rows, for one eye is 3. ...
```

Examples

```
data(vf.vfi)
colnames(vf.vfi)
progression.vfi(vf.vfi)
progression.vfi(vf.vfi[vf.vfi$eyeid == 1,])
progression.vfi(vf.vfi[vf.vfi$eyeid == 2,])
```

vfseries

Combined Visual Field Series for General Progression Method

Description

Data

Usage

```
data(vfseries)
```

Format

A data frame sample including the following visual field related measurement for two eyes, each with 10 follow-ups.

14 vfseries

Source

```
eyeid eyeid, labeled as 1,2... for different eyes.
nvisit number of visits.
yearsfollowed follow-up years.
distprev to be updated.
age in years.
righteye 1 as right eye, 0 as left eye.
malfixrate VF test malfixation rate.
ght glaucoma hemifield test result.
vfi visual field index.
md mean deviation.
mdprob mean deviation probability.
psd pattern standard deviation.
psdprob pattern standard deviation probability.
s1-s54 52 sensitivity measurements.
td1-td54 52 total deviation measurements.
tdp1-tdp54 52 total deviation probability measurements.
pdp1-pdp54 52 pattern deviation probability measurements. ...
```

```
data(vfseries)
progression(vfseries)
progression(vfseries[vfseries$eyeid == 1,])
progression(vfseries[vfseries$eyeid == 2,])
progression(vfseries, method=c("cigts"))
progression.cigts(vfseries)
progression(vfseries, method=c('plr.nouri.2012', 'schell2014', 'vfi'))
```

Index

```
* datasets
    vf.cigts, 11
    vf.plr.nouri.2012, 11
    vf.schell2014, 12
    vf.vfi, 13
    vfseries, 13
plotComponentMatrix, 2
plotfield.normalized, 4
plotTdProbabilities, 5
plotTDvalues, 6
progression, 7
{\tt progression.cigts}, {\tt 8}
progression.plr.nouri.2012, 8
progression.schell2014,9
\verb|progression.vfi|, 10
vf.cigts, 11
vf.plr.nouri.2012, 11
vf.schell2014, 12
vf.vfi, 13
vfseries, 13
```