# Package 'Epoch'

# September 3, 2025

| <u>.</u>  |
|---|
| Title IEEG (Intracranial Electroencephalography) Epoch Data Tools   |
| Version 1.0.3   |
| <b>Description</b> Provides tools for working with iEEG matrix data, including downloading curated iEEG data from OSF (The Open Science Framework <a href="https://osf.io/">https://osf.io/</a> ) (EpochDownloader()), making new objects (Epoch()), processing (crop() and resample()), and visualizing the data (plot()). |
| License GPL (>= 3)  |
| Encoding UTF-8  |
| RoxygenNote 7.3.2   |
| Biarch true   |
| <b>Depends</b> R (>= $4.1$ )  |
| <b>Imports</b> ggplot2, methods, ramify, TableContainer, glue, osfr, jsonlite, rlang, ggtext, gsignal   |
| SystemRequirements GNU make   |
| Suggests knitr, rmarkdown, testthat (>= 3.0.0)  |
| Config/testthat/edition 3   |
| VignetteBuilder knitr   |
| NeedsCompilation no   |
| Author Jiefei Wang [aut, cre]   |
| Maintainer Jiefei Wang <szwjf08@gmail.com></szwjf08@gmail.com>  |
| Repository CRAN   |
| <b>Date/Publication</b> 2025-09-03 21:30:07 UTC   |
|   |
| Contents  |
| .checkIndex         .standardizeIEEG         coltimes         crop         dim,Epoch-method   |

2 .standardizeIEEG

|       | Epoch                     | 5  |
|-------|---------------------------|----|
|       | Epoch-class               |    |
|       | EpochDownloader           | 7  |
|       | EpochDownloader-class     |    |
|       | get_config_data           | 10 |
|       | plot,Epoch,missing-method | 10 |
|       | resample                  | 11 |
|       | show, Epoch-method        | 12 |
| Index |                           | 13 |
|       |                           |    |

.checkIndex

Check and keep valid index only

# Description

Check and keep valid index only

# Usage

.checkIndex(indices, names)

# Arguments

indices Numeric or character index to check

names Character. All names corresponding to the indices

 $. \verb|standardize| IEEG| Standardize| iEEG| row| data for plotting$ 

# Description

Standardize iEEG row data for plotting

# Usage

.standardizeIEEG(data)

# **Arguments**

data Matrix or data frame of iEEG data

# Value

Standardized data matrix

coltimes 3

coltimes

Obtain the time points for the Epoch matrix

# **Description**

Obtain the time points for the Epoch matrix

#### Usage

```
coltimes(x)
## S4 method for signature 'Epoch'
coltimes(x)
```

#### **Arguments**

Χ

An Epoch object

#### Value

A numeric vector of time points, or column indices if time points are not defined

# See Also

Other Epoch methods: crop(), plot, Epoch, missing-method, resample(), show, Epoch-method

# **Examples**

```
# Create an Epoch object
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)
# get the time points of an Epoch object
coltimes(epoch)</pre>
```

crop

Methods for Epoch class

# Description

Truncating iEEG data to a specific time range.

4 dim,Epoch-method

#### Usage

```
crop(x, start, end, ...)
## S4 method for signature 'Epoch'
crop(x, start, end, checkTimeRange = TRUE)
```

#### **Arguments**

x An Epoch object

start Numeric value specifying start of new time range end Numeric value specifying end of new time range

... Not used

checkTimeRange Logical, whether to check the validity of the time range. This includes checking

if the time range is empty, if start is greater than end, and if start or end are

out of bounds. Default is TRUE.

#### Value

clip the time range of the Epoch object

#### See Also

Other Epoch methods: coltimes(), plot, Epoch, missing-method, resample(), show, Epoch-method

#### **Examples**

```
# Create an Epoch object
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)
# crop the epoch
crop(epoch, start = 0.5, end = 1.5)</pre>
```

dim, Epoch-method

Wrapper functions for calling TableContainer methods

#### **Description**

Wrapper functions for calling TableContainer methods

Epoch 5

#### Usage

```
## S4 method for signature 'Epoch'
dim(x)

## S4 method for signature 'Epoch'
dimnames(x)

## S4 method for signature 'Epoch'
x[i, j, ..., drop = TRUE]
```

# Arguments

| X    | An Epoch object   |
|------|---|
| i    | Row indices for subsetting. If only $i$ is provided, it will return the entire $row(s)$ . |
| j    | Column indices for subsetting.  |
|      | Additional arguments.   |
| drop | Not used.   |

#### Value

[: A new Epoch object with the selected data.

# **Examples**

```
# Create an Epoch object
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)

# wrappers
dim(epoch)
dimnames(epoch)
epoch[1]</pre>
```

Epoch

Constructor for Epoch class

# Description

Constructor for Epoch class

Epoch Epoch

# Usage

```
Epoch(
  table,
  electrodes = NULL,
  times = NULL,
  startTime = NULL,
  samplingRate = NULL,
  rowData = NULL,
  colData = NULL,
  metaData = NULL
)
```

# Arguments

| table        | Matrix containing epoch data (rows=electrodes, columns=time points)  |
|--------------|--|
| electrodes   | Optional character vector for electrode names, if not provided, row names of data are used. If row names are also not available, there will be no electrode names. |
| times        | Optional numeric vector of time points.  |
| startTime    | Optional numeric value for start time, if provided, times will be calculated based on this and samplingRate.   |
| samplingRate | Optional numeric value for sampling rate, if provided, times will be calculated based on this and startTime.   |
| rowData      | Optional data frame containing metadata for rows (electrodes).   |
| colData      | Optional data frame containing metadata for columns (time points).   |
| metaData     | Optional list containing metadata for the Epoch object. Element name "SamplingRate" is reserved by the Epoch class.  |

# Value

An Epoch object

# **Examples**

```
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)</pre>
```

Epoch-class 7

Epoch-class

Epoch Class

#### **Description**

S4 class to handle epoch data with electrodes and time points

#### Slots

```
table a matrix containing iEEG data (columns=time points, rows=electrodes) colData a data frame containing metadata for columns (time points) rowData a data frame containing metadata for rows (electrodes) metaData a list containing metadata for the Epoch object
```

**EpochDownloader** 

EpochDownloader constructor

#### **Description**

Manually update the project list from the remote repository. This function will attempt to download the latest configuration from the GitHub repository. If it fails, the current configuration will remain unchanged.

This function returns the list of available projects. If the configuration has not been fetched yet, it will automatically update the project list from the remote repository.

#### Usage

```
EpochDownloader(id = NULL, progress = TRUE, verbose = FALSE, path = NULL)
## S4 method for signature 'EpochDownloader'
names(x)
## S4 method for signature 'EpochDownloader'
x[i]
## S4 method for signature 'EpochDownloader'
x$name
## S4 method for signature 'EpochDownloader'
x[[i]]
## S4 method for signature 'EpochDownloader'
show(object)
```

8 EpochDownloader

```
## $4 method for signature 'EpochDownloader'
length(x)

wiki(x, ...)

## $4 method for signature 'EpochDownloader'
wiki(x)

updateRepos(verbose = FALSE)

EpochRepos(verbose = TRUE)
```

#### **Arguments**

id Either the ID of an OSF project or the name of an iEEG projects (case insensi-

tive). Check the available projects using EpochRepos(). The default points to

the fragility data from the Fragility multi-center retrospective study.

progress Logical indicating whether to show progress during download.

verbose Logical indicating whether to show messages

path The path to the temporary folder where the files will be downloaded.

x An EpochDownloader object.

i Index or name of the files to be accessed.

name The name of the file to be accessed.

object An EpochDownloader object.

... Not used, for future extensibility

#### Value

EpochDownloader: An EpochDownloader object.

names: A character vector of file names.

[: A named list of Epoch objects. The names are the dataset names.

\$: A single Epoch object.

[[: A single Epoch object.

show: Prints a summary of the EpochDownloader object.

length: Returns the number of files in the EpochDownloader object.

wiki: Opens the wiki page in the default browser

updateRepos: No return value, called for side effects.

EpochRepos: A list of project names and their corresponding OSF project IDs.

EpochDownloader-class

#### **Examples**

```
# list all available projects
EpochRepos()

# downloader for the fragility data
dl <- EpochDownloader(id = "fragility")

# list all Epoch objects in the downloader
names(dl)

# download the first Epoch object

dl[1]
 # equivalent to (index by name)
dl[names(dl)[1]]

# download the multiple Epoch objects

dl[c(1, 2)]
 # equivalent to (index by name)
dl[names(dl)[c(1, 2)]]

EpochRepos()</pre>
```

EpochDownloader-class EpochDownloader

#### **Description**

EpochDownloader is a class that allows downloading and accessing files from a OSF project.

#### **Slots**

```
id The ID of the OSF project.
files The files in the OSF project.
dataNames The names of the files in the OSF project.
tmp_folder The temporary folder where the files are downloaded.
progress Logical indicating whether to show progress during download.
```

get\_config\_data

Get configuration data from remote URL

# **Description**

Get configuration data from remote URL

# Usage

```
get_config_data()
```

#### Value

A list of project configurations

```
plot, Epoch, missing-method
```

Plot method for Epoch objects

#### **Description**

Plot method for Epoch objects

#### Usage

```
## S4 method for signature 'Epoch,missing'
plot(
    x,
    y,
    gap = 2,
    groupIndex = NULL,
    timeResolution = 2048,
    maxLabels = 50,
    x.lab.size = 2,
    ...
)
```

# Arguments

x An Epoch object

y Not used (for S4 method compatibility)

gap Numeric value specifying the gap between electrode traces (default: 2)

groupIndex Integer or string. A group of electrodes to show together in a different color. If

NULL(default), all electrodes are shown in the same color.

resample 11

```
timeResolution Maximum number of time points to keep for each electrode (default: 2048)

maxLabels Maximum number of electrode labels to display on the y-axis (default: 50)

x.lab.size Size of the x-axis label text (default: 2)

Additional arguments (not currently used)
```

#### Value

plot: A ggplot object showing iEEG electrode traces

#### See Also

```
Other Epoch methods: coltimes(), crop(), resample(), show, Epoch-method
```

#### **Examples**

```
# Create an Epoch object
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)
# Plot the epoch
plot(epoch)</pre>
```

resample

Generic function for resampling objects

# Description

This function allows you to resample an object to a different sampling frequency.

This function allows you to resample an Epoch object to a different sampling frequency.

#### Usage

```
resample(x, ...)
## S4 method for signature 'Epoch'
resample(x, samplingRate, ...)
```

#### **Arguments**

```
x An Epoch object to be resampled.... Additional arguments passed to gsignal::resamplesamplingRate The new sampling frequency (unit: Hertz).
```

show,Epoch-method

#### Value

An Epoch object with the resampled data.

#### See Also

Other Epoch methods: coltimes(), crop(), plot, Epoch, missing-method, show, Epoch-method

#### **Examples**

```
# Create an Epoch object
epoch_data <- matrix(rnorm(1000), nrow = 10)
rownames(epoch_data) <- paste0("Electrode_", 1:10)
epoch <- Epoch(epoch_data, startTime = 0, samplingRate = 100)
# downsample the epoch to 50 Hz
resample(epoch, samplingRate = 50)
# upsample the epoch to 200 Hz
resample(epoch, samplingRate = 200)</pre>
```

show, Epoch-method

Print the Epoch Object

#### **Description**

Print the Epoch Object

# Usage

```
## S4 method for signature 'Epoch'
show(object)
```

#### **Arguments**

object

Epoch object

# Value

returns an invisible NULL

#### See Also

```
Other Epoch methods: coltimes(), crop(), plot, Epoch, missing-method, resample()
```

# **Index**

```
* Epoch downloader
                                                 names, EpochDownloader-method
    EpochDownloader-class, 9
                                                          (EpochDownloader), 7
* Epoch methods
                                                 plot, Epoch, missing-method, 10
    coltimes, 3
    crop, 3
                                                 resample, 3, 4, 11, 11, 12
    plot, Epoch, missing-method, 10
                                                 resample, Epoch-method (resample), 11
    resample, 11
    show, Epoch-method, 12
                                                 show, Epoch-method, 12
.Epoch (Epoch-class), 7
                                                 show, EpochDownloader-method
.EpochDownloader
                                                          (EpochDownloader), 7
        (EpochDownloader-class), 9
.checkIndex, 2
                                                 updateRepos (EpochDownloader), 7
.standardizeIEEG, 2
[,Epoch-method(dim,Epoch-method),4
                                                 wiki (EpochDownloader), 7
[,EpochDownloader-method
                                                 wiki, EpochDownloader-method
        (EpochDownloader), 7
                                                          (EpochDownloader), 7
[[,EpochDownloader-method
        (EpochDownloader), 7
$,EpochDownloader-method
        (EpochDownloader), 7
coltimes, 3, 4, 11, 12
coltimes, Epoch-method (coltimes), 3
crop, 3, 3, 11, 12
crop, Epoch-method(crop), 3
dim, Epoch-method, 4
dimnames, Epoch-method
        (dim, Epoch-method), 4
Epoch, 5
Epoch-class, 7
EpochDownloader, 7
EpochDownloader-class, 9
EpochRepos (EpochDownloader), 7
{\tt get\_config\_data}, {\tt 10}
length,EpochDownloader-method
        (EpochDownloader), 7
```